

## Product datasheet for **MR228953**

### Sh3gl3 (NM\_001277954) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sh3gl3 (NM\_001277954) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** Sh3gl3  
**Synonyms:** EEN-B2; Sh3d2c; Sh3d2c2; SH3P13  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**ORF Nucleotide Sequence:** >MR228953 representing NM\_001277954  
**Red=Cloning site Blue=ORF Green=Tags(s)**

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCTGAACACTGTGTCGAAGCTCCGAGGGCAGGTGAAGGCCACTGGCTACCCACAGACGGAAGGCTTGC  
TGGGGACTGCATGCTGAAGTATGGCAAGGAGCTCGGAGAAGACTCTGCTTTGGCAACTCGTTGGTAGA  
TGTTGGTGAAGCCCTGAACTCATGGCTGAGGTGAAAGACTCTCTGGATATTAATGTGAAGCAAACCTTT  
ATTGACCCACTGCAGCTACTGCAAGACAAAGATTTAAAGGAGATCGGGCACCCCTGAGAAAGCTGGAAG  
GCCGTCGCTGGATTATGATTATAAAAAGCGGGGTAGGTAAGATCCCCGAGGAAGAAATCAGACAAGC  
AGTAGAGAAGTTGAAGAGTCAAAGGAGTTGGCCGAAAGGAGCATGTTTAATTTTTAGAAAATGATGTA  
GAGCAAGTGAGCCAGCTGGCTGTGTTGTAGAGGGCGCATTAGACTATCACAGGCAGTCCACAGAGATCC  
TCCAGGAGCTGCAGAGCAAGCTGGAGTTGCGAATATCTCTTGATCCAAAGTCCCCAAGCGAGAATTCAT  
GCCAAAGCCTGTGAACATGAGTTCCACCGATGCCAATGGGGTCCGACCCAGCTTTCATCAAAGACACCA  
GGTACTGACACTCCCTCGGACCAGCCCTGCTGTCGTGGTCTCTATGACTTTGAGCCAGAAAATGAAGGAG  
AATTAGGATTTAAGAAGGGGACATCATTACATTAACCAATCAGATAGATGAAAAGTGGTATGAAGGGAT  
GCTTCGTGGGAATCCGGCTTCTCCCCATTAATTACGTGGAAGTCATTGTGCCTTTACCTCCG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR228953 representing NM\_001277954  
Red=Cloning site Green=Tags(s)

MLNTVSKLRGQVKATGYPQTEGLLGDCMLKYGKELGEDSAFGNSLVDVGEALKLMAEVKDSLDINVKQTF  
 IDPLQLLQDKDLKEIGHHLRKLGRRLDYDKRRVVGKIPEEEIRQAVEKFEESEKELAERSMFNFLENDV  
 EQVSQLAVFVEAALDYHRQSTEILQELQSKLELRISLASKVPKREFMPKPVNMSSTDANGVGPSSSSKTP  
 GTDTPSDQPCCRGLYDFEPENEGELGFKEGDIITLTNQIDENWYEGMLRGESGFFPINYVEVIVPLPP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001277954  
**ORF Size:** 834 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001277954.2</a>
<b>RefSeq Size:</b>	1641 bp
<b>RefSeq ORF:</b>	837 bp
<b>Locus ID:</b>	20408
<b>Cytogenetics:</b>	7 47.15 cM
<b>MW:</b>	31.8 kDa
<b>Gene Summary:</b>	Implicated in endocytosis. May recruit other proteins to membranes with high curvature (By similarity).[UniProtKB/Swiss-Prot Function]