

Product datasheet for MR228957

Tm6sf1 (NM_001291282) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tm6sf1 (NM_001291282) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tm6sf1
Synonyms: AI428514; BC023123; C630016D09Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228957 representing NM_001291282
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGACGCACTACTTGAGAGAGGGTGAACCGTATCTGAACACGGCCTATGGGCACATGATCTGCTACTGGG
ATGGCTCTGTTTCATTACCTGATGTACCTAGTGATGGTGGCCGCATAGCCTGGGAGGAAAGCTACAGAAC
CATTGGCCTATACTGGGTCGGATCCATTATAATGAGCATCGTGGTTTTTGTGCCTGGAAACATCGTAGGA
AAGTATGGGACACGGATTTGCCCTGCATTTTTCTTAAGCATACCATATACATGTCTCCGGTCTGGGCTG
GCTTCAGAATCTATAATCAACCATCAGAAACTATAATTATCCCTCAAAGTTCTTCAAGAAGCCCAAGC
AAAGGCCCTGCTGAGAAGACCGTTTGACTTAGTGTGGTTTTGTGTCTTTTCTTGGAACCTGGATTTTGC
CTGTTCAGAGGCTTGATTGCTTTGGATTGCCAGCTGAGCTCTGCCGACTCTACACGCAGTTTCAAGAGC
CCTACCTAAAGGACCCTGCTGCTTATCCTAAGATTGAGATGCTGGCATAACATGTTCTATCCGTCCTTA
CTTTGTGATTGCCCTTTATGGCTTAGTGGTCCCTGGATGCTCATGGATGCCTGATATAACGCTGGTCCAT
GCTGGAGGTCTAGCTCAGGCTCAGTTTTCTCACATTGGTGCCTCCCTCATGCTAGAAGTCTTATGTCT
ACAGAGTCCCCGAAGAAGCGAAAATCTTTTTCTAGCACTAAACATAGCATATGGAGTTCTTCTCAGCT
CTTGGCCTACCGTTGTACCTACAATCCCAGATTCTTCTAAGAAGTAAAGGCAGATGAGAAGTTGGAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001291282.1 , NP_001278211.1
RefSeq Size:	1726 bp
RefSeq ORF:	840 bp
Locus ID:	107769
UniProt ID:	P58749
Cytogenetics:	7 D3
MW:	32.3 kDa
Gene Summary:	May function as sterol isomerase.[UniProtKB/Swiss-Prot Function]