

Product datasheet for MR228968

Nanog (NM_001289830) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nanog (NM_001289830) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nanog
Synonyms:	2410002E02Rik; ecat; ecat4; EN; ENK
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR228968 representing NM_001289830 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCTGCAGTTTTTCATCCCGAGAACTATTCTTGCTTACAAGGGTCTGCTACTGAGATGCTCTGCACAG
AGGCTGCCTCCTCGCCCTTCTCTGAAGACCTGCCTTTCAAGGCAGCCCTGATTCTTCTACCGTCC
CAAACAAAAGCTCTCAAGTCTGAGGCTGACAAGGCCCTGAGGAGGAGGAGAACAAGGTCTTCCAGG
AAGCAGAAGATGCGGACTGTGTTCTCTCAGGCCAGCTGTGTGCACTCAAGGACAGGTTTCAGAAGCAGA
AGTACCTCAGCCTCCAGCAGATGCAAGAATCTCTCCATTCTGAACCTGAGCTATAAGCAGGTTAAGAC
CTGGTTTCAAACCAAAGGATGAAGTGCAAGCGGTGGCAGAAAAACCAAGTGGTTGAAGACTAGCAATGGT
CTGATTCAGAAGGGCTCAGCACCAGTGGAGTATCCAGCATCCATTGCAGCTATCCCCAGGGCTATCTGG
TGAACGCATCTGGAAGCCTTCCATGTGGGGCAGCCAGACTTGGACCAACCCAATTGGAGCAGCCAGAC
CTGGACCAACCCAATTGGAACAACCAGACCTGGACCAACCCAATTGGAGCAGCCAGGCCTGGACCGCT
CAGTCTGGAACGCCAGCCTTGAATGCTGCTCCGCTCCATAACTTCGGGGAGGACTTTCTGCAGCCTT
ACGTACAGTTGCAGAAAATCTCTGCCAGTGATTTGGAGGTGAATTTGGAAGCCACTAGGAAAAGCCA
TGCGCATTTTAGCACCCACAAGCCTTGAATTATTCTGAACTACTCTGTGACTCCACCAAGTGAATA

ACGCGTACGCGGCCGCTCGAGCAGAAAATCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR228968 representing NM_001289830
 Red=Cloning site Green=Tags(s)

MPAVFHpenysclqgsateMLCTEaAsPRPSSedLPLQGSpsdStSPKQKLSSPEADKGPeeeeenKVLAR
 KQKMRTVFSQAQLCALKDRFQKQKYLslQMQELSSILNLsYKQVKTWFQnRMKCKRWQKNQWLKTSNG
 LIQKGSAPVEYPSIHCSYPQGYLVNAsGSLSMWGSQTwNPTWSSQTwNPTWNNQTwNPTWSSQAWTA
 QSWNGQPWNAAPLHNFGEDFLQPYVQLQqNFsASDLeVnLEATRESHAHFSTPQALELFLNYSVTPPPGEI


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001289830
ORF Size: 840 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001289830.1 , NP_001276759.1
RefSeq Size:	2020 bp
RefSeq ORF:	843 bp
Locus ID:	71950
UniProt ID:	Q80Z64
Cytogenetics:	6 F1
MW:	32.3 kDa
Gene Summary:	The protein encoded by this gene is a DNA binding homeobox transcription factor involved in embryonic stem (ES) cell proliferation, renewal, and pluripotency. The encoded protein can block ES cell differentiation and can also autorepress its own expression in differentiating cells. Several transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2015]