

Product datasheet for MR228990

Gjb1 (NM_001302498) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gjb1 (NM_001302498) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gjb1
Synonyms:	Cnx32; connex; connexin-32; Cx32; Cxng
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR228990 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGTGGACAGGTCTATACACCTTGCTCAGTGGCGTGAATCGGCACTCTACAGCCATTGGCCGAGTAT
GGCTGTCTGTATCTTCATCTTCAGAAATCATGGTGTGGTGGTGGCTGCTGAGAGCGTGTGGGGGATGA
GAAGTCTCTTTTCATCTGTAACACCTCCAGCCGGGCTGCAACAGCGTCTGCTATGACCATTTTTTCCCC
ATCTCCACGTGCGCCTATGGTCCCTGCAGCTTATCTTGGTTTCCACCCAGCTCTCCTCGTGGCAATGC
ACGTAGCTCACCAACAGCACATAGAAAAGAAAATGCTACGGCTTGAGGGGATGGGGACCCCTTACCT
GGAAGAGGTAAAGAGACACAAGGTGCACATCTCAGGGACTGTGGTGGACCTATGTCATCAGTGTGGT
TTCCGGTGTGTTCGAGGCTGTCTTCATGTATGTCTTCTATCTGCTCTACCCCGGCTATGCCATGGTGC
GGCTGGTCAAGTGTGAAGCCTTCCCCTGCCCAACACAGTGGACTGCTTCGTGTCCCGCCACCCGAGAA
AACCGTCTTCACTGTCTTATGCTCGCAGCCTCCGGCATCTGCATTATCCTCAACGTGGCGGAGGTGGT
TACCTCATCATCCGGCCTGTGCCCGCGTGTGTCAGCGCCGCTCCAATCCGCCCTCCCGAAGGGCTCGG
GCTTCGGCCACCGCCTCTCACCTGAATAACAAGCAGAATGAGATCAACAAGCTGCTGAGCGAGCAGGATGG
CTCTCTGAAAGACATACTGCGCCGAGCCCTGGCACAGGGGCGGGCTCGCTGAAAAGAGCGACCCGATGC
TCAGCCTGC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTAA



[View online »](#)

Protein Sequence: >MR228990 protein sequence
Red=Cloning site Green=Tags(s)

MNWTGLYTLISGVNRHSTAIGRVLSVIFIFRIMVLVVAESVWGDEKSSFCNTLQPGCNSVCYDHFPP
 ISHVRLWSLQLILVSTPALLVAMHVAHQHIEKKMLRLEGHGDPHLHEEVKRHKVHISGLWWTYVISVV
 FRLLFEAVFMYVFYLLYPGYAMVRLVKCEAFPCPNTVDCFVSRPTEKTVFTVFMLAASGICILNVAEVS
 YLIIRACARRAQRRSNPPSRKGSFGHRLSPEYKQNEINKLLEQDGLKDILRRSPGTGAGLAEKSDRC
 SAC

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001302498

ORF Size: 852 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001302498.1](#), [NP_001289427.1](#)

RefSeq Size: 1546 bp

RefSeq ORF: 852 bp

Locus ID: 14618


UniProt ID: [P28230](#)

Cytogenetics: X 44.06 cM

MW: 32 kDa

Gene Summary: This gene is a member of the gap junction protein (connexin) family. The encoded protein is a component of gap junctions, which are composed of arrays of intercellular channels that provide a route for the diffusion of ions and small molecules between cells. Mutations in a similar gene in human cause X-linked Charcot-Marie-Tooth disease, an inherited peripheral neuropathy. Alternative splicing results in multiple transcript variants encoding the same protein. [provided by RefSeq, Oct 2014]

Product images:



Circular map for MR228990