

Product datasheet for **MR228995**

Tmem200b (NM_001201367) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem200b (NM_001201367) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Tmem200b
Synonyms: EG623230
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228995 representing NM_001201367
Red=Cloning site Blue=ORF Green=Tags(s)

TTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGACGGCCGGGAGCCCCGGAGACGGCGGCGCGCAGGAGCCCCGAGGGCCGCGTTTCGCGCTGGGCC
GCCGCTTGGGCCGGCGGAGGGCGCGCTCCCCGCCAGAGCCGCTGCGGGTGGGGCTCGGCTGCGGCT
GCGTTCGCGCTGGGGGCGTTCGCGGCGCTGGGCGCGCTCGTGGTGTGGTGGGCATGGGCATCGCCGTG
GCCGGCTACTGGCCCCGGCCGCGCCTCCACACCCACGCCCGAGGACCGCCGCGCGCACGGGCCGACG
AGCGCCTGCGGCTGCTGGGACCGGTAATCATGGGCGTCGGCTGTTCGTGTTTCATCTGTGCCAACACATT
GCTCTACGAGAACC GCGACCTGGAGACCCGGCTCTCCGCCGGGGCATGCTGAGGGCTCAGGCGCTGCGG
CCGCCCGACGGCCAGGCTGGGACGCGCTGCTGCCAGTCCC GCGCATGGGTCCCCGGGCGCCGTCGCGG
AACCCGAACTTGGGACCTGGCCCCGAGGAGGGGCCCTCGCCGTCCCGTGGTGGCAGCCTGCGCTC
TGAACCTGCCAACCCTCGATCGGGGCTGCCTGCCCTGCTCCACAGCTACCCGTTGAAGGGCCCGGGGCTG
CCACCACCCTGGGTCTCGGACTCAGACTGGACACGTGATCATCACGGTGCAGCCCTCTGGTTCCTGCA
TCGAGCACTCAAGTCTCTGGACCTGGGCTAGGGGAGCTTCTTCTGGGGCGCCAGCAACTCGGGACTG
TGCCACCGAAGCTGGCCACGGCTGGACCGACTCAGTTTGGGGGCTATGCCAAATTAGGCGGAGATTTG
GGGCCCCGGGTC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR228995 representing NM_001201367
 Red=Cloning site Green=Tags(s)

MTAGSPGDGGARRSPEGRVSRLLGRRLGRRRRPRSPPEPLRVRARLRLRSPSGAFAALGALVVLVGMGIIV
 AGYWPGRASHTHAPRTGRAHGPHERLRLLLGPVIMGVGLFVFICANTLLYENRDLETRRLRRGMLRAQALR
 PPDGPGWDALLPSAHGSPGAVAEPETWDLAPRRGPPSPVSVRSRSEPANPRSGLPALLHSYPLKGPGL
 PPPWGPRTQTHVITVQPSGSCIEHKSLLDLGLGELLGAPATRDCAHRSWPRLDRLSLGGYAKLGGDL
 GARV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001201367

ORF Size: 852 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001201367.1 , NP_001188296.1
RefSeq Size:	1370 bp
RefSeq ORF:	855 bp
Locus ID:	623230
Cytogenetics:	4 D2.3
MW:	30.8 kDa