

Product datasheet for **MR228996**

Kcnip3 (NM_001291005) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kcnip3 (NM_001291005) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Kcnip3
Synonyms: 4933407H12Rik; AI413860; Csen; DREAM; KChIP3; R74849
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR228996 representing NM_001291005
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCCGGATCGCC

ATGAGGCAGCTGCCAGCCGGACCAAGCAGTCTGGCTTGTCTGGGCTGCAAAGCGGGAAGATTAGTGACGG
TCCCTTTCAGCAGCAGAGATGCAGAGGACCAAGGAAGCCGTGAAGGCATCAGATGGCAACCTCCTGGGAG
ATCCTGGGCGCATACCACTGAGCAAGAGGGAAAGCATCAAGTGGCAAAGGCCACGGTTCACCCGCCAGGC
CCTGATGCGTTGCTGCTTAATCAAGTGGATCCTGTCCAGTGTGCCCCACAAGGCTCAGCAGTGACAGTG
AACTGGAGTTATCCACGGTGCGCCATCAGCCAGAGGGCTTGGACCAGCTACAAGCTCAGACCAAGTTCAC
CAAGAAGGAGCTGCAGTCCCTTTACCGAGGCTTCAAGAATGAGTGTCCCACAGGCTGGTGGATGAAGAC
ACCTTCAAACCTATTTATCCAGTTCTTCCCTCAGGGAGATGCCACCACCTATGCACACTTCTCTTCA
ATGCCTTCGATGCTGATGGGAACGGGCCATCCACTTTGAGGACTTTGTGGTTGGGCTCTCCATCCTGCT
TCGAGGGACGGTCCATGAGAAGCTCAAGTGGGCCTTCAATCTCTATGACATTAACAAGGATGGTTGCATC
ACCAAGGAGGAGATGCTGGCCATCATGAAGTCCATCTACGACATGATGGGCCGCCACCTACCCATCC
TGCGGGAGGATGCACCCCTGGAGCATGTGGAGAGGTTCTTTCAGAAAATGGACAGGAACCCAGGATGGAGT
GGTGACCATTGATGAATTTCTGGAGACTTGTGAGAAGGATGAGAACATCATGAACTCCATGCAGCTGTTT
GAGAACGTCATC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR228996 representing NM_001291005
 Red=Cloning site Green=Tags(s)

MRQLPAGPSSLACSGCKAGRLVTVPFSSRDAEDQGSREGIRWQPPGRSWAHTTEQEGKHQVAKATVHPPG
 PDALLLNQVDPVQCCPTRLSSDSELELSTVRHQPEGLDQLQAQTKFTKKELQSLYRGFKNECPTGLVDED
 TFKLIYSQFFPQGDATTYAHFLFNAFDADGNGAIHFEDFVVGLSILLRGTVHEKLNKWFNLVDINKDGC
 TKEEMLAIMKSIYDMMGRHTYPIILREDAPLEHVERFFQKMDRNQDGVVTIDEFLETQKQDENIMNSMQLF
 ENVI


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001291005

ORF Size: 852 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001291005.1 , NP_001277934.1
RefSeq Size:	2771 bp
RefSeq ORF:	855 bp
Locus ID:	56461
UniProt ID:	Q9QXT8
Cytogenetics:	2 F1
MW:	32.6 kDa
Gene Summary:	Calcium-dependent transcriptional repressor that binds to the DRE element of genes including PDYN and FOS. Affinity for DNA is reduced upon binding to calcium and enhanced by binding to magnesium. Seems to be involved in nociception.[UniProtKB/Swiss-Prot Function]