

## Product datasheet for MR229048

### Ccnd3 (NM\_001081635) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Ccnd3 (NM_001081635) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Ccnd3
Synonyms:	9230106B05Rik; AA682053; AL024085; AW146355; C78795
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR229048 representing NM_001081635 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGGAGCTGCTGTGTTGCGAGGGCACCCGGCACGCGCCCCGGGCCGGCCGACCCGCGGCTGCTTGGGG  
ACCAGCGTGTCTGCAGAGTTTACTCCGCCTGGAGGAGCGCTACGTGCCGCGAGCCTCTACTTCCAGTG  
CGTGCAAAGGAGATCAAGCCGCACATGCGGAAGATGCTGGCATACTGGATGCTGGAGGTGTGTGAGGAG  
CAGCGCTGCGAGGAGGATGTCTCCCTCTGGCTATGAACCTGGATCGCTACCTGTCTCGCTCCCA  
CCCGAAAGGCGCAATTGCAGTTCTAGGTACTGTCTGCCTGTTGCTGGCCTCCAAGCTGCGCGAAACCAC  
GCCCTGACTATTGAGAAGCTTTGCATCTATACGGACCAGGCTGTGGCTCCATGGCAGTTGCGGGAGTGG  
GAGGTGCTGGTCTGGGGAAGCTCAAGTGGGACCTGGCTGCTGTGATTGCGCACGACTTCTGGCCTTGA  
TTCTGCACCGCTGTCTCTGCCAGTGACCGGCAGGCTTTGGTCAAAAAGCATGCCAGACCTTTTTGGC  
CCTCTGTGCTACAGATTACACCTTTGCGATGTATCCTCCATCCATGATCGCCACAGGCAGCATTGGGGCA  
GCCGTGCTAGGCCTAGGCGCCTGCTCTATGTCTGCGGATGAGCTCACAGAGTTGCTGGCCGGGATCACAG  
GCACTGAAGTGGACTACCTGCGAGCCTGCCAGGAACAGATCGAAGCTGCCCTCAGGGAGAGCCTCAGGGA  
AGCTGCTCAGACAGCCCCAGCCAGTGCCCAAAGCCCCCGGGGCTCTAGTAGCCAGGGGCCAGTCAG  
ACCAGCACTCCACAGATGTCACAGCCATTCACCTG

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR229048 representing NM\_001081635  
Red=Cloning site Green=Tags(s)

MELLCCGTRHAPRAGDPDLLGDQRLVQSLRLLEERYVPRASYFQCQVQKEIKPHMRKMLAYWMLEVCEE  
 QRCEEDVFPLAMNYLDRYLSCVPTRKAQLQLLGTVCLLLASKLRETTPLTIEKLCIYTDQAVAPWQLREW  
 EVLVLGKLGKLDAAVIAHDFLALILHRLSLPSDRQALVKKHAQTFALCATDYTFAMYPSPMIATGSIGA  
 AVLGLGACSMSADEL TELLAGITGTEVDYL RACQEQIEAALRESLREAAQTAPSPVPKAPRGSSSQGPSQ  
 TSTPTDVTAIHL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001081635

**ORF Size:** 876 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001081635.1](#), [NP\\_001075104.1](#)

**RefSeq Size:** 2086 bp

**RefSeq ORF:** 879 bp

**Locus ID:** 12445


**UniProt ID:** [P30282](#)

**Cytogenetics:** 17 23.37 cM

**MW:** 32.5 kDa

**Gene Summary:** Regulatory component of the cyclin D3-CDK4 (DC) complex that phosphorylates and inhibits members of the retinoblastoma (RB) protein family including RB1 and regulates the cell-cycle during G(1)/S transition. Phosphorylation of RB1 allows dissociation of the transcription factor E2F from the RB/E2F complex and the subsequent transcription of E2F target genes which are responsible for the progression through the G(1) phase. Hypophosphorylates RB1 in early G(1) phase. Cyclin D-CDK4 complexes are major integrators of various mitogenic and antimitogenic signals. Also substrate for SMAD3, phosphorylating SMAD3 in a cell-cycle-dependent manner and repressing its transcriptional activity. Component of the ternary complex, cyclin D3/CDK4/CDKN1B, required for nuclear translocation and activity of the cyclin D-CDK4 complex.[UniProtKB/Swiss-Prot Function]

### Product images:



Circular map for MR229048