

## Product datasheet for **MR229058**

### Jund (NM\_001286944) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Jund (NM_001286944) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Jund
Synonyms:	Jund1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR229058 representing NM_001286944 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGCTGAAGAAAGACGCGCTGACGCTCAGCCTGGCGGAGCAGGGAGCGGGGATTGAAACCAGGGTCGG  
CCACTGCACCTTCTGCGCTGCGCCCGACGGCGCCCGACGGGCTGCTGGTTTCGCCGGATCTTGGGCT  
GCTCAAACCTCGCGTCGCGGAGCTGGAGAGGCTGATCATCCAGTCCAACGGGCTGGTGACCACTACCCCG  
ACCAGTACGCAGTTCCTCTACCCGAAGGTGGCAGCTAGCGAGGAGCAGGAGTTCGCCGAAGGCTTCGTCA  
AGGCGCTGGAGGACCTGCACAAGCAAAGCCAGCTGGGTGCGGCCACCGCGCCACCTCAGGGGCTCCCGC  
GCCTCCCGCGCCCGCCGACCTGGCCGCCACCCCGGGGCCACGGAGACCCCGGTCTACGCCAACCTGAGC  
AGTTTCGCGGGTGGCGCCGGGCCCTGGGGCGCGGCCACCGTGGCTTTCGCCGCGGAGCCAGTGCCCT  
TCCCGCCGCCCCGGGCGCGCTGGGGCCGCCACCTCCGCATCCACCGCGCTGGCCGCGCTCAAGGA  
CGAGCCGACACCGTGCCGGACGTGCCGAGCTTCGGCGACAGCCCTCCGCTGTCGCCATCGACATGGAC  
ACGCAAGAACGCATCAAGGCGGAGCGCAAGAGGCTGCGCAACCGCATCGCCGCTCCAAATGCCGAAGC  
GCAAGCTGGAGCGTATCTCGCGCTGGAGGAGAAAGTCAAGACCCTCAAAGCCAGAACACCGAGCTGGC  
GTCCACCGCCAGCCTGTCGCGAGCAGGTGGCGCAGCTCAAACAGAAAGTCTCAGCCAGTCAACAGC  
GGTGCCAGCTGCTGCCACGACCAGGTCCCGCGCTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR229058 representing NM\_001286944  
 Red=Cloning site Green=Tags(s)

MLKKDAL T L S L A E Q G A A G L K P G S A T A P S A L R P D G A P D G L L A S P D L G L L K L A S P E L E R L I I Q S N G L V T T T P  
 T S T Q F L Y P K V A A S E E Q E F A E G F V K A L E D L H K Q S Q L G A A T A A T S G A P A P P A D L A A T P G A T E T P V Y A N L S  
 S F A G G A G P P G G A A T V A F A A E P V F P P P P G A L G P P P P H P P R L A A L K D E P Q T V P D V P S F G D S P P L S P I D M D  
 T Q E R I K A E R K R L R N R I A A S K C R K R K L E R I S R L E E K V K T L K S Q N T E L A S T A S L L R E Q V A Q L K Q K V L S H V N S  
 G C Q L L P Q H Q V P A Y

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

Sgfl-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:** NM\_001286944

**ORF Size:** 879 bp

<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001286944.1</a> , <a href="#">NP_001273873.1</a>
<b>RefSeq Size:</b>	2894 bp
<b>RefSeq ORF:</b>	882 bp
<b>Locus ID:</b>	16478
<b>Cytogenetics:</b>	8 34.15 cM
<b>MW:</b>	31 kDa
<b>Gene Summary:</b>	The protein encoded by this intronless gene is a member of the JUN family, and a functional component of the AP1 transcription factor complex. This protein has been proposed to protect cells from p53-dependent senescence and apoptosis. Alternative translation initiation site usage results in the production of different isoforms (PMID:12105216). [provided by RefSeq, Nov 2013]