

## Product datasheet for MR229121

### PrIr (NM\_001253782) Mouse Tagged ORF Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** PrIr (NM\_001253782) Mouse Tagged ORF Clone  
**Tag:** Myc-DDK  
**Symbol:** PrIr  
**Synonyms:** A1987712; Pr-1; Pr-3; PrIr-rs1  
**Mammalian Cell Selection:** Neomycin  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**ORF Nucleotide Sequence:** >MR229121 representing NM\_001253782  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGTCATCTGCACTTGCTTACATGCTGCTTGTCTCAGCATCAGCCTCCTGAATGGACAGTCACCACCTG  
 GAAAACCTGAAATCCACAAATGTCGTTCCCTGACAAGGAAACATTCACCTGCTGGTGAATCCTGGGTC  
 AGATGGAGGACTCCCACCAATTATTCATTGACATACAGCAAAGAAGGAGAGAAAAACCTATGAATGT  
 CCAGACTACAAAACCAAGTGGCCCAATTCCTGTTCTTTAGCAAGCAGTACACTTCCATATGGAAAAAT  
 ACATCATCACAGTAAATGCCACGAACGAAATGGGAAGCAGTACCTCGGATCCACTTTATGTGGATGTGAC  
 TTACATTGTTGAACCAGAGCCTCCTCGGAACCTGACTTTAGAAGTGAAACAATAAAGACAAAAAACA  
 TATCTGTGGGTAAAATGTTGCCACCTACCATAACTGATGTA AAAACTGGTTGGTTTACAATGGAATATG  
 AAATTCGATTAAGTCTGAAGAAGCAGATGAGTGGGAGATCCACTTCACAGGTCATCAAACACAATTTAA  
 GGTTTTTGACTTATATCCAGGACAAAAGTATCTTGTCCAGACTCGCTGCAAGCCAGACCATGGATACTGG  
 AGTAGATGGGGCCAGGAGAAATCTATTGAAATACCAATGACTTCACCTTGAAAGACACAACCTGTGTGGA  
 TCATTGTGGCCGTTCTCTGCTGTCATCTGTTGATTATGGTCTGGGCAGTGGCTTTGAAGGGTTATAG  
 CATGATGACCTGCATCTTCCACCAGTTCGGGGCCAAAAATAAAAGGATTTGATACTCATCTGCTAGAG  
 CTCTGGTGTCTATTCTCCAGTTGACTTCCCTGGTGAAGATTCCAACAACCTGAGTTTTTATGTGATTTA

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR229121 representing NM\_001253782  
Red=Cloning site Green=Tags(s)

MSSALAYMLLVLSISLLNGQSPPGKPEIHKCRSPDKETFTCWNP GSDGGLPTNYSLTSYKEGEKNTYEC  
 PDYKTS GPN SCFFSKQYTSIWKIYIITVNATNEMGSSTSDPLYVDVTYI VEPEPPRNLTLEVKQLKDKKT  
 YLWVKWLPPTITDVKTGWFTMEYEIRL KSEEADWEIHF TGHQTQFKVFDLYPGQKYL VQTRCKPDHGYW  
 SRWGQEK SIEIPNDFTLKDTTVWII VAVLSAVICLIMVVAVALKGYSMMTCIFPPVPGPKIKGFDTHLLE  
 LWCSILQLTSLVKIPTTEFLCDL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:** SgfI-MluI

**Cloning Scheme:**

Cloning sites used for ORF Shuttling:



\* The last codon before the Stop codon of the ORF

**ACCN:** NM\_001253782

**ORF Size:** 909 bp


**OTI Disclaimer:** Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at [custsupport@origene.com](mailto:custsupport@origene.com) or by calling 301.340.3188 option 3 for pricing and delivery.

The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_001253782.2</u> , <u>NP_001240711.1</u>
<b>RefSeq Size:</b>	1079 bp
<b>RefSeq ORF:</b>	912 bp
<b>Locus ID:</b>	19116
<b>UniProt ID:</b>	<u>Q08501</u>
<b>Cytogenetics:</b>	15 5.23 cM
<b>MW:</b>	35.2 kDa
<b>Gene Summary:</b>	This is a receptor for the anterior pituitary hormone prolactin.[UniProtKB/Swiss-Prot Function]

**Product images:**



Circular map for MR229121