

Product datasheet for MR229149

Tex264 (NM_001286498) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tex264 (NM_001286498) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tex264
Synonyms:	TEG-264
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR229149 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCCGGATCTCCTACTACTGGCCTGATTGGGGCCCTGACGCTGCTGTTGCTGCTGACGCTGCTGGCCT
TTGCTGGTTATTCAGGACTGCTGACTGGGGTGACAGTGAGCGCTGGATCACCCCAATCCGCAACATAAC
TGTGGCCTACAAGTTCACGTGGGGTCTATGGTGACACTGGGCACCTTTTCACAGAGAGCTGCAGCATC
TCTCCCAAGCTCCGTTCCATCGCTGTCTACTATGACAACCCCATACGGTGCCTCCTGAGAAGTGCCGCT
GTGCAGTCGGCAGCATCCTGAGTGAGGGGAGGAGTCGCCTTACCTGAGCTCATCCACCTTATCAGAA
ATTTGGCTTCAAGATATTCTCCTCCCAGCACCTAGCCATGTGGTCATAGCTACCTCCCTTACACCACC
CCCATATCCATCTGGCTGGCTGCCCGCCGAGTCCATCCTGCCTTGATACCTACATCAAGGAGCGGAAGC
TGTGTGCTCACCTCGCCTGGAGATCTACCATCAAGACAAGATCCATTTTCATGTGCCCACTGGCAAGGCA
AGGAGATTTCTACGTGCCAGAGGTGAAGGAGACAGAGCGGAAATGCCGGGAGCTTGGGAGGCCACTGAC
ACCCAGACGGATGGCACAGGAGCTGATACAAGTATGCAAGTTCTGTGAGCCTGGATGTTCCGCCCTGGCA
GCCGGGAGACTTCAGCCACCACACTTTCTCCTGGGGCAGGCAACCGTGGCTGGGACGACGGTGACAACCG
CAGCGAGCACAGCTACAGTGAATCGGGTCCAGTGGCTCATCCTTTGAGGAGCTGGACCTGGAGGGCGAG
GGACCTTGGGAGAACCCTGACTGAACCCTGAAGCCAAGCTTCTGGGGCCCCCTCGGAGCTCAGCACCC
CTGAGAGGGGTGAGGAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229149 protein sequence
Red=Cloning site Green=Tags(s)

MPDLLLLGLIGALTLLLLLLTLLAFAGYSGLLTGVTVSAGSPPIRNITVAYKFHVGSYGDTHLFTESCSI
 SPKLRSIAVYYDNPHTVPPEKCRCAVGSILSEGEESPSPELIHL YQKFGFKIFSFAPSHVVIATFPYTT
 PISIWLAARRVHPALDTYIKERKLC AHPRL E IYHQDKIHFMCLARQGD F YVPEVKETERKRELAETD
 TQTDGTGADTSDASSVSLDVRPGSRETSATTLSPGAGNRGWDDGDNRSEHSYSSES G ASGS SFEELDLEGE
 GPLGEPRLNPEAKLLGPPRELSTPERGEE

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001286498

ORF Size: 927 bp


OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001286498.1 , NP_001273427.1
RefSeq Size:	1651 bp
RefSeq ORF:	930 bp
Locus ID:	21767
UniProt ID:	E9Q137
Cytogenetics:	9 F1
MW:	33.6 kDa
Gene Summary:	<p>Major reticulophagy (also called ER-phagy) receptor that acts independently of other candidate reticulophagy receptors to remodel subdomains of the endoplasmic reticulum into autophagosomes upon nutrient stress, which then fuse with lysosomes for endoplasmic reticulum turnover. The ATG8-containing isolation membrane (IM) cradles a tubular segment of TEX264-positive ER near a three-way junction, allowing the formation of a synapse of 2 juxtaposed membranes with trans interaction between the TEX264 and ATG8 proteins. Expansion of the IM would extend the capture of ER, possibly through a 'zipper-like' process involving continued trans TEX264-ATG8 interactions, until poorly understood mechanisms lead to the fission of relevant membranes and, ultimately, autophagosomal membrane closure. Also involved in the repair of covalent DNA-protein cross-links (DPCs) during DNA synthesis: acts by bridging VCP/p97 to covalent DNA-protein cross-links (DPCs) and initiating resolution of DPCs by SPRTN.[UniProtKB/Swiss-Prot Function]</p>

Product images:



Circular map for MR229149