

Product datasheet for **MR229151**

Art5 (NM_001291354) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Art5 (NM_001291354) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Art5
Synonyms: ARTC5; Yac-2
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR229151 representing NM_001291354
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGATTCTGGAGGATCTGCTGATGGTCCTCAGCTGCCTCAGCTTGACGCCCCTCTGGAAGGTTTCGAGCTG
TTCCTATCCTGCCCTGAGCCTGGTTCCAGATACCTTTGATGATGCCTATGTGGGCTGCTCTGAGGAGAT
GGAGGAGAAAGCAGGCCTGCTGCTAAAGGAAGAGATGGCAGCCATGCCCTGCTGCGGGAATCCTGGGAA
GCAGCACAAGAGGCTGGGCACACCGGCGTCACAAGCTCACTCTACCTCCTGGCTTCAAAGCCCAGCAGC
GAGTAGCGATTATGGTGTACACCAACTCATCCAACACCTTGTACTGGGAGCTGAACCAGGCCGTACGGAC
AGGAGGTGGCTCCCGGAGCTCTACATGAGGCACTTCCCCTTCAAGGCCCTGCATTTCTACTTGACCCGG
GCTCTGCAGCTGCTGCGGGGCTCTGGAGGGTGCAGTAGGGGACCTGGGGAGGTGGTGTCCGAGGTGTGG
GCAGTCTTCACTTTGAACCCAAGAGGCTGGGGGACTCTGTCCGATTAGGACAGTTTACCTCCAGCTCTGT
GGATGAGAGAGTGGCTCGCAGGTTTGGGAATGCCACCTTCTTCAATCTAAGGACTTGTTTGGGGCTCCT
ATCCAGGCGTTGTCTGTCTTCCCTGAGGAGCGTGAGGTGCTGATACCCACACGAAGTCTTCTTGGTCA
CTGGGTTCTCCAGGATGGAGCCAGAGCATAGTACTCTCTGGAGTTATGATCAGACCTGCAGCCACTT
TAACTGCGCCTATTTGGGTGGGAGAAGAGACGTGGCTGTGTCTCAAGAGCAGTCGGGCAGCCAGAG
GCACCCCTCACAGAGGCTCTGGCTCTGCAATCAGGAAAGACACTGCTCTTGGACCAAGGAAGCTGCAGC
TCTCCAGAGCTGGACCC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229151 representing NM_001291354
Red=Cloning site Green=Tags(s)

MILEDLLMVL SCLSLHALWKVRAVPILPLSLVPDTFDDAYVGCSEEMEEKAGLLLKEEMARHALLRESWE
 AAQEAWAHRRHKLTLPPGFKAQHGVAIMVYTNSNTLYWELNQAVRTGGGSREL YMRHFPFKALHFYLTR
 ALQLLRGSGGCSRGPGEVFRGVGSLHFEPKRLGDSVRLGQFTSSVDERVARRFGNATFFNLRTCFGAP
 IQALSVFPEEREVLIPPHEVFLVTGFSQDGAQSIIVTLWSYDQTC SHFNCAYLGGEKRRGCVSSRAVGQPE
 APSTEALALQSGKTL LLDPRKQLSRAGP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001291354

ORF Size: 927 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291354.1](#), [NP_001278283.1](#)

RefSeq Size: 1662 bp

RefSeq ORF: 930 bp


Locus ID: 11875

UniProt ID: [P70352](#)

Cytogenetics: 7 54.71 cM

MW: 34.3 kDa

Product images:



Circular map for MR229151