

Product datasheet for MR229328

P2ry14 (NM_001287119) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: P2ry14 (NM_001287119) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: P2ry14
Synonyms: A330108O13Rik; Gpr105; P2Y14
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
ORF Nucleotide Sequence: >MR229328 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGAACAACCTCCACCACCACAGACCCTCCAAACCAGCCCTGCTCCTGGAACACCCTGATCACAAGCAGA
 TCATTCCTGTTGTACGGTATGGTCTTCATCACGGGGCTCCTCCTCAATGGGATATCAGGATGGATATT
 CTTTTATGTGCCAGCTCCAAGAGTTTTATCATCTATCTCAAGAACATAGTGGTGGCTGACTTTCTCATG
 GGCTGACTTTCCCTTTCAAAGTCCTTGGTGACTCAGGCCTCGGCCCTGGCAGGTGAATGTGTTGTGT
 GCAGGGTCTCTGCCGTCATCTTCTATGTTAATATGTACGTCAGCATCGTGTCTTTGGGCTCATCAGCTT
 TGACAGGTAATAAAAATTGTGAAGCCCTTCTGACGCTATTGTGACGTCGGTGAATAAGCAAGCTG
 CTTTCTGTGCTCGTGTGGATGCTCATGCTTCTCCTTGCTGTCCAAACATCATCCTGACAAACCAGGGTG
 TCAAGGAGGTCACGAAGATACAGTGCATGGAGCTCAAAAACGAGCTGGGGCGGAAGTGGCACAAGGCGTC
 TAACTATATCTTCGTGAGTATCTTCTGGGTCGTGTTTCTTCTGCTAATCGTCTTACACGGCCATCAGC
 AGGAAGATCTTCAAGTCTCACCTCAAGTCCAGGAAGAATTCCACCTCCGTCAAGAGGAAGTCCAGCCGCA
 ATATCTTCAGCATCGTCTGCTTTTTGTGCTGCTTTGTGCCTTACCACATTGCCAGAATCCCCTACAC
 AAAGAGTCAGACGAAGGTCACCTACAGTGCCTGGACGAAGGAGACCCTGCTCTATGCCAAGAATCACT
 CTGCTACTCTGGCTGCCAATGTGTCTGGACCCATTATTTATTTCTTCTTATGCCAGCCATTTAGAG
 AAGTCTTAAATAAGAAGTTACACATGTCACTCAAAGTCCAGAATGACCTAGAGGTTTCCAAAACAAAAG
 GGAAAATGCGATTTCATGAAAGCACAGATACTTTG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229328 protein sequence
Red=Cloning site Green=Tags(s)

MNNSTTTDPPNQPCSWNTLITKQIIPVLYGMVFITGLLLNGISGWIFFYVPSKSFIIYLKNIIVADFLM
 GLTFPFKVLGDSGLGPWQVNVFVCRVSAVIFYVNMVVSIVFFGLISFDRIYKIVKPLLTISIVQSVNYSKL
 LSVLVWMLMLLLAVPNIIILTNQGVKEVTKIQCMELKNELGRKWHKASNYIFVSIWVVFLLLIIVFYTAIT
 RKIFKSHLKSARKNSTSVKRKSSRNIFSVLVFVVCVFPYHIARIPYTKSQTEGHYSCRTKETLLYAKEFT
 LLLSAANVCLDPIIYFFLCQPFREVLNKKLHMSLKVQNDLEVSCTKRENAIHESDTL

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001287119

ORF Size: 1017 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001287119.1](#), [NP_001274048.1](#)

RefSeq Size: 1617 bp

RefSeq ORF: 1017 bp

Locus ID: 140795


UniProt ID: [Q9ESG6](#)

Cytogenetics: 3 28.96 cM

MW: 38.9 kDa

Gene Summary: Receptor for UDP-glucose coupled to G-proteins.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR229328