

Product datasheet for **MR229410A1V**

Mouse Pdlim2 (NM_001253736) AAV Particle

Product data:

Product Type: AAV Particles
Product Name: Mouse Pdlim2 (NM_001253736) AAV Particle
Tag: Myc-DDK
Symbol: Pdlim2
Synonyms: 4732462F18Rik; Slim
Mammalian Cell Selection: None
Vector: pAAV-AC-Myc-DDK (PS100089)
ORF Nucleotide Sequence: >MR229410 representing NM_001253736
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCCGGATCGCC

ATGGCGTTGACTGTGGATGTGGCAGGACCAGCACCTTGGGGCTTCCGAATTAGCGGGGGCAGAGATTTCC
 ACACACCCATCATTGTGACCAAGGTCACAGAGCGGGGCAAGGCTGAAGCAGCTGATCTCCGGCCTGGCGA
 CATCATTGTGGCCATCAATGGACAGAGTGCAGAGAACATGCTACACGCGGAGGCCAAAGCAAGATCCGA
 CAGAGCGCCTCACCCCTAAGACTGCAGCTGGACCGTCCCAAACAGCCTCTCCTGGGCAGACCAATGGGG
 AGGGCTCCTTGAAGTGCTGGCAACCAGATTCCAGGGTCCCTGAGGACACACCGTGACAGCCAGTCTTC
 CCAGAGGTCTGCCTGCTTCAGCCAGTCTCTCTCAGCCCCAGGCCTTGCAGCCCCCTTCTCCACCCACCC
 CCTACCAGCCCAGTTGCCCTTTCTAAAGAGGATATGATTGGCTGTAGTTTCCAGAGTCTGACACACTCTC
 CAGGCCTTGCTGCTGCTCACCATTGACCTACCCTGGCCACCCACCAGCCAACAGGCCGGCCACAGCAG
 CCCAAGCGACTCCGCAGTGAGGGTGCTGCTCCATTCCCCAGGACGGCCCTCCAGCCCTAGGTTTCAGCAGT
 TTGGATCTGGAGGAAGACTCAGAGGTGTTCAAGATGCTGCAGGAGAACCGCCAGGGACGGGCCGCCCAA
 GGCAGTCCAGCTTTTTCGACTTTACAGGAAGCCTTGGAGGCTGAGGAGAGAGGTGGCACACCTGCCTT
 TGTGCCAGCTCGCTGAGCTCCAGGCTTCTTGGCCACCTCCAGGGCCTTGGCCACTCCACCAAGCTC
 CACACCTGTGAGAAATGCAGCGTCAACATCTCGAACCAGCGGTCCGCATCCAGGAGGGGAGGTACCGAC
 ACCCTGGCTGCTACACTTGCAGACTGTGGGCTGAACCTGAAGATGCGCGGCCACTTCTGGGTGGGCAA
 TGAGTTGTACTGCGAGAAGCATGCCCGCAGCGCTACTCTATGCCTGGAAGTCTCAACTCTCGAGCC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence:	>MR229410 representing NM_001253736 Red=Cloning site Green=Tags(s) MALTVDVAGPAPWGFRISSGGRDFHTPIIVTKVTERGKAEAADLRPGDIIVAINGQSAENMLHAEQSKIR QSASPLRLQLDRSQTASPGQTNGEGSLEVLATRFQGLRTHRDSQSSQRSACFSPVLSPRPCSPFSTPP PTSPVALSKEDMIGCSFQSLTHSPGLAAHHLTYPGHPTSQQAGHSSPSDSAVRVLLHSPGRPSSPRFSS LDLEEDSEVFKMLQENRQGRAAPRQSSSFRLLEALEAEERGGTPAFVPSLSSQASLPTSRALATPPKL HTCEKCSYINISQAVRIQEGRYRHPGCYTCADCGLNLKMRGHFWVGNELYCEKHARQRYSMPTLNSRA TRTRPLEQKLI SEEDLAANDILDYKDDDDKV
Species:	Mouse
Serotype:	AAV-2
ACCN:	NM_001253736
ORF Size:	1047 bp
Buffer:	PBS with 0.001% Pluronic F68
Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001253736.1 , NP_001240665.1
RefSeq Size:	1587 bp
RefSeq ORF:	1050 bp
Locus ID:	213019
UniProt ID:	Q8R1G6
Cytogenetics:	14 D2
MW:	37.7 kDa