

Product datasheet for **MR229432**

Pbx3 (NM_001290576) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pbx3 (NM_001290576) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Pbx3
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR229432 representing NM_001290576
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGGACGATCAATCCAGGATGCTGCAGACTCTGGCCGGGGTGAACCTGGCCGGCCACTCGGTCCAGGGG
GTATGGCTCTGCCGCTCCCCCTCACGGCCACGAAGGGGCGGACGGCGACGGCAGGAAGCAGGACATCGG
CGACATCCTCCACCAGATCATGACCATCACCGACCAGAGCTTGGACGAGGCGCAAGCAAAGAAACATGCT
CTGAACCTGTCACAGAATGAAACCGCGCTCTTCAGCGTCTTGTGTGAGATCAAAGAGAAAACAGGTCTCA
GCATCAGAGGAGCCAGGAGGAAGACCTCCCGACCCAGCTAATGAGACTGGACAATATGCTTTTGGC
AGAAGGGGTTTCAGGTCTGAGAAAGGTGGGGATCGGCGGCAGCAGCTGCAGCCGCGGCAGCCTCTGGA
GGTTCTTCAGATAACTCTATTGAACACTCAGATTACAGAGCCAAATTGACCCAGATCAGACAAATCTATC
ACACAGAACTGGAGAAATATGAACAGGCATGTAATGAATTCACAACACATGTGATGAACCTTCTCCGAGA
ACAGAGTAGAACGCGTCCCATTCCCGAAAGAGATCGAAAGGATGGTGGGCATCATCCATCGGAAATTT
AGTTCATTTCAGATGCAGCTCAAGCAAAGCACGTGTAAGCAGTTATGATTTTGAGATCAAGGTTCTTG
ATGCCAGACGGAAAAGGCGTAATTCAGTAAACAGGCTACGGAATCTTGAATGAATATTTTACTCACA
CCTCAGCAACCCCTACCCAGTGAAGAAGCCAAAGAGGAGCTGGCCAAGAAATGCAGCATTACAGTATCA
CAGGTATCCAACCTGGTTTGGCAACAAACGAATCAGGTACAAGAAGAACATCGGCAAAATTTCAAGAGGAGG
CCAATCTCTATGCTGCAAAAACAGCCGTGACAGCCGCACATGCAGTGGCTGCAGCCGTGCAGAACAACCA
GACCAACTCACCCACCACGCCAAACTCTGGTGGATACCCCTCCGTATGTTATCAATCAGACGGGAGGCTA
CAG

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA



[View online »](#)


Protein Sequence: >MR229432 representing NM_001290576
 Red=Cloning site Green=Tags(s)

MDDQSRMLQTLAGVNLAGHSVQGGMALPPPHGHEGADGDGRKQDIGDILHQIMTITDQSLDEAQAKKHA
 LNCHRMKPALF SVLCEIKEKTGLSIRGAQEEDPPDPQLMRLDNMLLAEGVSGPEKGGGSAASAAAAAASG
 GSSDNSIEHSDYRAKL TQIRQIYHTELEKYEQACNEFTTHVMNLLREQSRTRPISPKEIERMVGIIHRKF
 SSIQMLKQSTCEAVMILRSRFLDARRKRRNF SKQATEILNEYFYSHLSNPYPSEEAKEELAKKCSITVS
 QVSNWFGNKRIRYKKNIGKFQEEANLYAAKTAVTAAHVAHVAAVQNNQNTNSPTTPNSGGYPPSCYQSDGRL
 Q


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-RsrII

Cloning Scheme:



Plasmid Map:



ACCN: NM_001290576

| | |
|-------------------------------|---|
| ORF Size: | 1053 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001290576.1 , NP_001277505.1 |
| RefSeq Size: | 2757 bp |
| RefSeq ORF: | 1056 bp |
| Locus ID: | 18516 |
| UniProt ID: | O35317 |
| Cytogenetics: | 2 22.84 cM |
| MW: | 39.2 kDa |
| Gene Summary: | Transcriptional activator that binds the sequence 5'-ATCAATCAA-3'. [UniProtKB/Swiss-Prot Function] |