

Product datasheet for MR229507

Gsg1 (NM_010352) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gsg1 (NM_010352) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gsg1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR229507 representing NM_010352 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGTTCAGAAGGGTCTCTGACCAGCGGACCTTCATTTCTGCCATCCTCAACATGTTATCACTCG
GCCTCTCCACAGCATCCTTGCTCAGCAGCGAGTGGTTTGTGGGCACACAGAAGGTGCCAAGCCCTGTG
CGGGCAAAGTCTGGCAGCCAAGTGCTTTGACATGCCGATGCTTGGATGGGGGCATCGCCAACACATCA
GCCCAGGAGGTGGTACAATACACCTGGGAGACTGGGGATGACCGTTCTCCTTCCCTTGCCTCCGCAGTG
GCATGTGGCTATCCTGTGAGAAACCATGGAAGAGCCAGGGGAGAAGTGCAGACGTTTCATTGAACTCAC
ACCACCAGCCCAGAGAGGTGAGAAAGGACTACTGGAATTTGCCACGTTGCAAGGCTCATGTACCCCACT
CTCCGATTTGGAGGGGAGTGGTTAATGGAGAAGGCTTCTCTCCTCCACCTCCCTTGGGGGCCGTGGCAA
AGGTCTTTGGCTGTCACTGGGAGCCAGACTGCCTATATCGGACTTCAACTCATCAGCTTCTCCTGCT
ACTGACGGATCTGTTGCTCACCACGAACCCTGGCTGTGGGCTCAAGCTAAGCGGTTTGCAGCCGTCTCC
TTGGTCTGTGAGGACTTCTGGGATGGTGGCTCATATGCTATATCACAAAGTCTCCAGGCAACTGCCA
ACTTAGGTCGGGAGGACTGGAGACCACACTCTTGAATTACGGCTGGGCCTTCTACACAGCGTGGGTTTC
CTTCACCTGCTGCATGGGTCAGCGGTACCACCTTCAACATGTACACGAGGATGGTGTGGAGTTCAAG
TGACGGACAGCAAGAGCTTTAACCAACCCAGCTGCCTGGCGCAGCACCACCCTGTTTCTTCCCTC
CTCCGCTGACGTGCAACCCACGCAGGGGAACCTTTGTCCAGCTGCCATCAGTACCCAGCCATCCCAT
CCGCTCTGTCTCTGAAGCTATTGACCTCTACTCGCGCTACAGGACAAAGAATTTCAACAAGGGATCAGC
CAGGAGCTAAAGGAAGTGGTCGAGCCATCTGTAGAAGAGCAGCGT

ACGCGTACGCGGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229507 representing NM_010352
 Red=Cloning site Green=Tags(s)

MEFQKSSDQRTFISAILNMLSLGLSTASLLSSEWVFGTQKVPKPLCGQSLAAKCFDMPMSLDGGIANTS
 AQEVVQYTWETGDDRFSLAFRSGMWLSCEETMEEPGEKCRRFIELTPPAQRGEKGLLEFATLQGSCHPT
 LRFGGEWLMEKASLLHLPWGPVAKVFWLSLGAQTAYIGLQLISFLLLLTDLTTPNPGCGLKLSAFAAVS
 LVL SGLLGMVAHMLYSQVFQATANLGPEDWRPHSWNYGWAFYTAWVSFTCCMASAVTTFNMYTRMVLEFK
 CRHSKSFNTNPSCLAQHRCFLPPFLTCTTHAGEPLSSCHQYPSHPIRSVSEAIDLYSALQDKEFQQGIS
 QELKEVVEPSVEEQR


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_010352

ORF Size:	1095 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_010352.2 , NP_034482.2
RefSeq Size:	1448 bp
RefSeq ORF:	1098 bp
Locus ID:	14840
UniProt ID:	Q8R1W2
Cytogenetics:	6 66.07 cM
MW:	41.1 kDa
Gene Summary:	May cause the redistribution of PAPOLB from the cytosol to the endoplasmic reticulum. [UniProtKB/Swiss-Prot Function]