

Product datasheet for **MR229581**

Camk1d (NM_001290376) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Camk1d (NM_001290376) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Camk1d
Synonyms: A630059D12Rik; CaMKIdelta; CKLiK; E030025C11Rik; mCKLiK
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR229581 representing NM_001290376
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGGAGTTTGTGCTTGGTCTGCTTAAATTCAGATGGAGCTGGATAAAGGGTCCAGAAACAGTG
GGCCTTTCTGAAGTTGTTTTAGCCGAGGAGAAAGCTACTGGGAAGCTCTTCGAGTGAAGTGCATCCC
GAAGAAGCGCTGAAGGGCAAGGAGAGCAGCATCGAGAACGAGATTGCCGTGCTTAGAAAGATTAAGCAT
GAAAACATTGTTGCCTTGAAGATATTTATGAAAGCCAAATCACCTCTACCTGGTCATGCAACTTGTGT
CTGGTGGAGAAGCTCTCGATCGGATAGTGGAGAAGGGTTTTACACAGAGAAAGATGCCAGCACTCTCAT
CCGCCAGGTCCTGGATGCCGTATACTATCTCCACAGAATGGGCATTGTCCACAGGGACCTCAAGCCGGAG
AATCTCTTATACTACAGTCAAGACGAGGAGTCCAAAATAATGATCAGTGAAGTGGCTTGTGCAAAATGG
AGGGCAAAGGAGATGTGATGTCCACGGCTGCGGGACCCAGGCTATGTTGCTCCGGAAGTTCTCGCCCA
GAAACCGTACAGCAAAGCTGTGGACTGCTGGTCCATCGGGGTGATCGCCTATATCTTGTCTGTGGTTAC
CCTCCTTTTTATGATGAAAATGACTCGAAGCTGTTTGAACAGATCCTCAAGGCAGAATATGAGTTTGATT
CCCCCTACTGGGATGACATCTCGACTCTGCCAAAGACTTCATTGCGAATCTGATGGAGAAGACCCAAA
TAAAAGATACACTGTGAGCAGGAGCTCGACACCCATGGATTGCTGGTGACACAGCCCTTAGCAAAAAC
ATTCACGAATCTGTGAGTCCCAGATCCGGAAGAATTTGCAAGAGCAAATGGAGACAAGCGTTTAAACG
CCACGGCAGTCGTGAGACATATGCGGAGGCTCCAGCTTGGCAGCAGCCTGGACAGTTCAAATGCAAGTGT
CTCAAGCAACCTCAGTTTGGCCAGCCAAAAGATTGTCTGGCACCTTCCACGCTCTGTAGTTTCTTTCT
TCTTCGTCGGGGTTCGAGGAGTTCGGAGCTGAGAGGAGACCCAGGCCACCCTGTGACAACAGGGCACA
CTGGAAGCAAG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229581 representing NM_001290376
 Red=Cloning site Green=Tags(s)

MAEFVSWSCLNFRWSWIKGSRNSGAFSEVVLAEKATGKLFVAVKCIKKKALKGKESSIENEIAVLRKIKH
 ENIVALEDIYESPNHL YLVMQLVSGGELFDRIVEKGFYTEKDASTLIRQVLDVAVYVLRHMGIVHRDLKPE
 NLLYYSQDEESKIMISDFGLSKMEGKGDVMSTACGTPGYVAPEVLAQKPYSKAVDCWSIGVIAYILLCGY
 PPFYDENDSKLFEQILKAEYEFDSPYWDDISDSAKDFIRNLMKDPNKRYTCEQAARHPWIAGDTALSKN
 IHESVSAQIRKNFAKSKWRQAFNATAVVRHMRRLQLGSSLDSSNASVSSNLSLASQKDC LAPSTLCSFLS
 SSSGVAGVGAERRPRPTTVTTGHTGSK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001290376

| | |
|-------------------------------|--|
| ORF Size: | 1131 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001290376.1 , NP_001277305.1 |
| RefSeq Size: | 6742 bp |
| RefSeq ORF: | 1134 bp |
| Locus ID: | 227541 |
| UniProt ID: | Q8BW96 |
| Cytogenetics: | 2 A1 |
| MW: | 42.6 kDa |
| Gene Summary: | Calcium/calmodulin-dependent protein kinase that operates in the calcium-triggered CaMKK-CaMK1 signaling cascade and, upon calcium influx, activates CREB-dependent gene transcription, regulates calcium-mediated granulocyte function and respiratory burst and promotes basal dendritic growth of hippocampal neurons. In neutrophil cells, required for cytokine-induced proliferative responses and activation of the respiratory burst. Activates the transcription factor CREB1 in hippocampal neuron nuclei. May play a role in apoptosis of erythroleukemia cells. In vitro, phosphorylates transcription factor CREM isoform Beta (By similarity). Isoform 1 but not isoform 2 activates CREB1.[UniProtKB/Swiss-Prot Function] |