

Product datasheet for MR229660

Atf2 (NM_001284373) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Atf2 (NM_001284373) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Atf2
Synonyms:	Atf-2; CRE-BP; Creb2; D18875; D130078H02Rik; mXBP; Tg(Gzma-Klra1)7Wum
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR229660 representing NM_001284373 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGTGATGACAAACCTTTCTATGCACTGCCCTGGGTGTGGCCAGCGTTTTACCAACGAGGATCATT
TGGCTGTCCATAAACATAAACATGAGATGACACTGAAATTTGGTCCAGCACGTAATGACAGTGTCATTGT
GGCTGATCAGACTCCAACGCCAACAAAGATTCTAAAAAAGTGTGAAGAAGTGGGTTTGTTCATGAGTTG
GCAAGTCCATTTGAAAAATGAATTCAGAAAGGCTTCCGAAGATGACATTAATAAATGCCTCTAGATTTGT
CCCTCTTGCAACACCCATCATAAGAAGCAAAATGAGGAGCCTTCTGTTGTAGAAAACAACCTACCAGGA
CAGCCCTTTACCTCACCCGAGTCGACTACCAGTGATGAAAAGCTTGTTCGGCCAGTCACCATGGTGCT
AGTGTTCAGGAATCCCAGGCCCTTCTCTCCTCAACAGTCCAGTCAGAAGCAAAAATGAGATTAAG
CTGCTTTGACCCAGCAACACCCCTCCAGTTACCAATGGTGATACTGTAAGGCCATGGCAGTGGATTGGT
TAGGACTCAGTCAGAAGAGTCTCGCCACAGTCTTGCAGCAGCCAGCCACCTCCACTACAGAACTCCG
GCTTCTCCAGTCCACAACTCCTCAGACCCAAAATAAAGTGGCCGTCGAAGAAGAGCAGCTAATGAAG
ATCCTGATGAGAAAAGGAGGAAGTTTCTAGAAGCAAAATAGAGCAGCAGCTTCAAGATGCCGACAAAAAG
GAAAGTGTGGTTTCAGTCTTAGAGAAGAAAGCAGAAGACTTGAAGTCACTAAATGGCCAGCTGCAGAGC
GAAATCACCTGCTGAGAAATGAAGTGGCCAGCTGAAAACAGCTTCTTCTGGCTCATAAAGATTGCCCTG
TAACTGCCATGCAGAAGAAGTCTGGCTATCATACTGCTGATAAAGATGACAGTTCAGAAGACCTTTCTGT
GCCAAGCAGTCCACATACAGAAGCGATCCAGCACAGCTCTGTGAGCAGTCCCAATGGAGTCAAGTCAACA
TCAAAAGCAGAAGCTGTAGCCACTTCAGTCTCACCCAGATGGCCGACCAGAGCACGGAGCTGCACCTT
CACAGATTGTCATGGCTCTCCCTCCAGGCACAGCCCTCAGGAAGT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229660 representing NM_001284373
 Red=Cloning site Green=Tags(s)

MSDDKPFLLCTAPGCGQRFNEDHLAVHKHKHEMTLKFGPARND SVIVADQTPPTPTRLKNCEEVGLFNEL
 ASPFENEFFKASEDDIKKMPDLDSLPLATPIIRSKIEEPSVVETTHQDSPLPHPESTTSDEKLVRPVTMVP
 SVPGIPGPSSQPVQSEAKMRLKAAL TQQHPPVTNGD TVKGHGSGLVRTQSEESRPQSLQQPATSTTETP
 ASPAHTTPQTQNTSGRRRRAANEDPDEKRRKFLERNRAAASRCRQKRKVVVQSLEKKAEDLSSLNGQLQS
 EVTLLRNEVAQLKQLLLAHKDCPVTAMQKKSgyHTADKDDSSEDLSVPSSPHTEAIQHSSVST SNGVSST
 SKAEAVATSVLTQMDQSTEPALSQIVMAPPSSQAQPSGS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001284373

ORF Size:	1167 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001284373.1 , NP_001271302.1
RefSeq Size:	3943 bp
RefSeq ORF:	1170 bp
Locus ID:	11909
Cytogenetics:	2 C3
MW:	42.8 kDa
Gene Summary:	Transcriptional activator which regulates the transcription of various genes, including those involved in anti-apoptosis, cell growth, and DNA damage response. Dependent on its binding partner, binds to CRE (cAMP response element) consensus sequences (5'-TGACGTCA-3') or to AP-1 (activator protein 1) consensus sequences (5'-TGACTCA-3'). In the nucleus, contributes to global transcription and the DNA damage response, in addition to specific transcriptional activities that are related to cell development, proliferation and death. In the cytoplasm, interacts with and perturbs HK1- and VDAC1-containing complexes at the mitochondrial outer membrane, thereby impairing mitochondrial membrane potential, inducing mitochondrial leakage and promoting cell death. The phosphorylated form (mediated by ATM) plays a role in the DNA damage response and is involved in the ionizing radiation (IR)-induced S phase checkpoint control and in the recruitment of the MRN complex into the IR-induced foci (IRIF). Exhibits histone acetyltransferase (HAT) activity which specifically acetylates histones H2B and H4 in vitro. In concert with CUL3 and RBX1, promotes the degradation of KAT5 thereby attenuating its ability to acetylate and activate ATM. Can elicit oncogenic or tumor suppressor activities depending on the tissue or cell type (By similarity).[UniProtKB/Swiss-Prot Function]