

Product datasheet for MR229700

Oprm1 (NM_001302793) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Oprm1 (NM_001302793) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Oprm1
Synonyms:	m; M-OR-1; MO; MOP-R; mor; MOR-1; MOR-1O; mu; muOR; O; Oprm
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR229700 representing NM_001302793 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGACAGCAGCGCCGGCCAGGAACATCAGCGACTGCTCTGACCCCTTAGCTCCTGCAAGTTGCTCCC
CAGCACCTGGCTCCTGGCTCACTTGTCCACGTTGATGGCAACCAGTCCGACCCATGCGGTCTAAACCG
CACGGGGCTTGGGGGAGCCACAGCCTGTGCCCTCAGACCGGCAGCCCTCCATGGTCACAGCCATCACC
ATCATGGCCCTCTATTCTATCGTGTGTGTAGTGGCCTCTTTGGAACTCCTGGTCATGTATGTGATTG
TAAGATATACAAAATGAAGACTGCCACCAACATCTACATTTTCAACCTTGCTCTGGCAGATGCCTTAGC
CACTAGCACGCTGCCCTTTCAGAGTGTTAACTACCTGATGGGAACGTGGCCCTTTGGAAACATCCTCTGC
AAGATCGTGATCTCAATAGACTACTACAACATGTTCAACAGTATCTTACCCTCTGCACCATGAGTGTAG
ACCGCTACATTGCCGTCTGCCACCGGTCAAGGCCCTGGATTTCCGTACCCCGGAAATGCCAAAATTGT
CAATGTCTGCAACTGGATCCTCTCTTCTGCCATTGGTCTGCCGTAATGTTTCATGGCAACCACAAAATAC
AGGCAGGGGTCCATAGATTGCACCCTCAGTTCCTCTCATCCACATGGTACTGGGAGAACCTGCTCAAAA
TCTGTGTCTTTCATCTTCCGCTTTCATCATGCCGTCCTCATCATCACTGTGTGTTATGGCATGATCCTT
ACGACTCAAGAGTGTCCGATGCTGTCCGGCTCCAAGAAAAGGACAGGAACCTGCGCAGGATCACCCGG
ATGGTGTGTTGGTCTGTGCTGATTTATTGTCTGCTGGACCCCATCCACATCTATGTCATCATCAAAG
CACTGATCACGATTCCAGAAACCACTTTCCAGACTGTTTCTGGCACTTCTGCATTGCCCTGGGTTACAC
AAACAGCTGCCTGAACCCAGTCTTTATGCGTTCCTGGATGAAAACCTCAAACGATGTTTATAGAGAGTTC
TGCATCCCAACTCCTCCACAATCGAACAGAAAACCTCTGCTCGAATCCGTCAAACACTAGGGAACACC
CCTCCACGGCTAATACAGTGGATCGAATAACCACCAGCTAGAAAATCTGGAAGCAGAACTGCTCCATT
GCCC

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229700 representing NM_001302793
Red=Cloning site Green=Tags(s)


MDSSAGPGNISDCSDPLAPASCSAPAGSWLNLSHVDGNQSDPCGPNRTGLGGSHSLCPQTGSPSMVTAIT
 IMALYSIVCVVGLFGNLFVMYVIVRYTKMKTATNIYIFNLALADALATSTLPFQSVNYLMGTWPFGNILC
 KIVISIDYNNMFTSIFTLCTMSVDRIYAVCHPVKALDFRTPRNAKIVNVCNWILSSAIGLPVMFMATTKY
 RQGSIDCTLTFSHPTWYWENLLKICVFIFAFIMPVLIITVCYGLMILRLKSVRMLSGSKEKDRNLRRITR
 MVLVVAVFIVCWTPIHIVYVIKALITIPETTFQTVSWHFCIALGYTNSCLNPVLYAFLDENFKRCRFRE
 CIPTSSTIEQQNSARIRQNTREHPSTANTVDRTNHQLENLEAETAPLP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001302793

ORF Size:	1194 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001302793.1 , NP_001289722.1
RefSeq Size:	2230 bp
RefSeq ORF:	1197 bp
Locus ID:	18390
UniProt ID:	P42866
Cytogenetics:	10 1.85 cM
MW:	44.9 kDa
Gene Summary:	This gene encodes the mu opioid receptor which is where drugs such as morphine and other opioids have pharmacological effects. Several alternatively spliced transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Nov 2014]