

Product datasheet for MR229706

Clasp2 (NM_001286599) Mouse Tagged ORF Clone

Product data:

| | |
|--------------------------|--|
| Product Type: | Expression Plasmids |
| Product Name: | Clasp2 (NM_001286599) Mouse Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Clasp2 |
| Synonyms: | 1500004F14Rik; 8030404L10Rik; C77448; CLASP2beta; mKIAA0627 |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Cell Selection: | Neomycin |
| ORF Nucleotide Sequence: | >MR229706 representing NM_001286599 Red=Cloning site Blue=ORF Green=Tags(s) |

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGGAGCCCCGCGCGCCGAGTACTTCTGCGCCCAAGTGTGCAGAAGGACGTGAGCGGCCGGCTGCAGG
CGGGCGAGGAGCTGCTGCTCTGCCTCGGCACGCCCGGAGCCATCCCGGACCTGGAGGACGATCCGAGCCC
CCTGGCTAAGACCGTGGACGCGCTCACCCGCTGGGTGGGCTCGAGCAACTACCGGTATCGTTGTTAGGA
TTGGAGATTTTAAAGTGCCTTTGTGGATAGATTGTCAACACGATTCAAATCTTATGTAACAATGGTTACCA
CAGCTTTAATAGACAGAATGGGAGATGTCAAAGACAAAGTTTCGAGAAGAAGCTCAGAACCTGACACTGAA
GCTGATGGATGAGGTGGCGCCGCCTATGTACATCTGGGAGCAGCTAGCTTCTGGCTTCAAGCACAAGAAC
TTTCGATCTCGAGAAGGCGTGTGCCTGTGTCTTATTGAGACCTTGAACATCTTTGGGACTCAACCACTGG
TCATCAGCAAGTTGGTGCCGATTTATGTGTCTTATTGGAGACTCTAACAGTCAGGTGAGAAATGCTGC
ACTATCAGCTGTAGTGGAGATTTATAGACACGTAGGAGAGAAGCTGAGGATTGATTTGTGTAAGAGAGAC
ATCCCCCTGCTAGATTAGAAAATGGTACTTGCCAAATTTGATGAAGTTCAGAATCTGGCCGTATGATTT
TGAGCGTCTGCAAAGATAAAAGCTTTGATGATGAGGAATCAGTGGATGGAATAGGCCGTCTGAGCTGC
TTCAGCCTTCAAGTTCCCTGCACCTAAACACCTGGGAATCCTGTGAGCAGTGCAAGAAAGCCTGGCTCA
GCAGGTGGCCCTAAGTTGGAGTCTTCTAAAGAAGGAGGGGCTGGAGCAGTTGATGAAGATGACTTTA
TAAAGCTTTTACAGATGTTCTTCTGTTGAGTCTATTCTAGTCGAGAAGTGAAGAGACGTTAAATAA
GATCAGGGAAATTTGTGATGACAAACATGACTGGGACCAGCGTGCCAATGCGTGCATGCTCTGCAGC
CTGGTAGCCAGTGAAGTTGAGAGAGCTCTGGCCTGGCCTCTGTGGAGTCACCTAGCTACTTATCAGCACC
ATTGTCATTGCACCCTGAGTCATGAGGACTTGCCAGAGAAGAGATTGACTTCTCCTGTGAGCTCGGCATT
TGTTCAACAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online >](#)

Protein Sequence: >MR229706 representing NM_001286599
 Red=Cloning site Green=Tags(s)

MEPRGAEYFCAQVLQKDVSGRLQAGEELLLCLGTPGAIPDLEDDPSRLAKTVDALTRWVGSSNYRVSLLG
 LEILSAFVDRLSTRFKSYVTMVTALIDRMGDVKDKVREEAQLTLKLMDEVAPPMYIWEQLASGFKHKN
 FRSREGVCLCL IETLNIFGTQPLVISKLVPHLCVLF GDSNSQVRNAALSAVVEIYRHVGEKLRIDLCKRD
 IPPARLEMVLAKFDEVQNSGGMILSVCKDKSFDDEESVDGNRPSSAASAFKVPAPKTPGNPVSSARKPGS
 AGGPKVGGPSKEGGAGAVDEDDFIKAFDTVPSVQIYSSRELETLNKIREILSDDKHDWDQRANACMSCS
 LVASEVERALAWPLWSHLATYQHHCHCTL SHEDLPEKRLTSPVSSAFVQH

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001286599

| | |
|-------------------------------|--|
| ORF Size: | 1200 bp |
| OTI Disclaimer: | The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info |
| OTI Annotation: | This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene. |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001286599.1 , NP_001273528.1 |
| RefSeq Size: | 2279 bp |
| RefSeq ORF: | 1203 bp |
| Locus ID: | 76499 |
| UniProt ID: | Q8BRT1 |
| Cytogenetics: | 9 F3 |
| MW: | 44.5 kDa |
| Gene Summary: | Microtubule plus-end tracking protein that promotes the stabilization of dynamic microtubules. Involved in the nucleation of noncentrosomal microtubules originating from the trans-Golgi network (TGN). Required for the polarization of the cytoplasmic microtubule arrays in migrating cells towards the leading edge of the cell. May act at the cell cortex to enhance the frequency of rescue of depolymerizing microtubules by attaching their plus-ends to cortical platforms composed of ERC1 and PHLDB2. This cortical microtubule stabilizing activity is regulated at least in part by phosphatidylinositol 3-kinase signaling. Also performs a similar stabilizing function at the kinetochore which is essential for the bipolar alignment of chromosomes on the mitotic spindle. Acts as a mediator of ERBB2-dependent stabilization of microtubules at the cell cortex.[UniProtKB/Swiss-Prot Function] |