

Product datasheet for MR229712

Kcnj15 (NM_001271687) Mouse Tagged ORF Clone

Product data:

Product Type: Expression Plasmids
Product Name: Kcnj15 (NM_001271687) Mouse Tagged ORF Clone
Tag: Myc-DDK
Symbol: Kcnj15
Synonyms: 4930414N08Rik; AI182284; AI267127; IRKK; Kir4.2
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
ORF Nucleotide Sequence: >MR229712 representing NM_001271687
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGTAGCCAGGTGGGAGAAGGGGAGCGAGGACGCGCCACTTACCCTGCAGAAGATTCTGACCTACAGA
 GTGGTCCAAGATGGATGCCATTCACCTTGGCATGTCCAGTGCCCACTGGTGAAGCATACCAACGGGGT
 TGGACTCAAGGCCACAGACCCCGAGTCATGTCAAAGAGTGGGCACAGTAATGTGAGAATCGATAAGGTA
 GACGGAATCTATTTACTCTACCTCCAGGACTTGTGGACAACCGTCATCGACATGAAGTGGCGATACAAGC
 TCAACCTATTTGCTGCCACCTTTGTGATGACCTGGTTTCTGTTTGGAGTGGTCTACTATGCCATAGCCTT
 TATTCATGGTGACTTACAACCTTGGGGAATCTAATTCCAACCACACACCCTGCATTATGAAAGTGGACTCT
 CTCACAGGAGCATTCTCTTTTCTTGGAACTTCAGACAACCAATTGGCTACGGGGTCCGTTCCATCACAG
 AGGAGTGTCCCATGCTATCTTCCTCTTAGTCGCCCAACTGGTCATCACCACATTGATTGAGATCTTCAT
 TACGGGGACCTTTCTGGCTAAAATTGCAAGACCCAAAAGCGAGCCGAGACCATTAAGTTCAGCCACTGT
 GCTGTATCAGCAAGCAGAATGGAAAGCTATGCCTGGTCATCCAGGTGGCCAACATGAGGAAGAGTCTCC
 TGATTCAGTGCCAGCTCTCTGAAAACCTCTGCAGACACAGTCACCAAAGAGGGAGAACGCATTCTCCT
 CAACCAGGCCACTGTCAAATCCAGTGGACTCCTCTTCCGAGAGTCCCTTCTCATCTGCCCCATGCC
 TTCTACCACGTGTTGGATGAGACAAGCCCCCTGCGGGACCTCACACCCAAAACCTAAAGGAGAAGGAGT
 TTGAGCTGGTGGTACTTCTCAACGCCACGGTGGAGTCTACCAGCGCCGTCTGCCAGAGCCGAACGTCTTA
 CATCCCGGAGGAGATCTACTGGGCTTTGAGTTTGTGCCTGTGGTTTCTCTCTCCAAAATGGAAAGTAT
 GTGGCTGATTTCAAGTCAATTTGAGCAGATCAGGAAGAGCCGGATTGTACCTTCTACTGTCCGATTCTG
 AGAAGCAGAAGCTTGAAGAACAGTACAGGCAAGAGGACCAGAGGGAGCGGGAGCTGAGGAGCCTCCTGCT
 ACAGCAGAGCAATGTC

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229712 representing NM_001271687
 Red=Cloning site Green=Tags(s)


MVARWEKGSSEDAPLTLQKIPDLQSGPKMDAHLGMSSAPLVKHTNGVGLKAHRPRVMSKSGHSNVRIDKV
 DGIYLLYLQDLWTTVIDMKWRYKLTFLAATFVMTWFLFGVYVYAI AF IHGDLQLGESNSNHTPCIMKVDS
 LTGAFLFSLESQTTIGYGVRSITEECPHAIFLLVAQLVITTLIEIFITGTFLAKIARPKKRAETIKFSHC
 AVISKQNGKLC LVIQVANMRKSLLIQCQLSGKLLQTHVTKEGERILLNQATVKFHVDSSESFPFLILPMT
 FYHVLDETSPLRDLTPQNLKEKEFELVLLLNATVESTSAVCQSRTSYIPEEIYWGFEFVPPVVSLSKNGKY
 VADFSQFEQIRKSPDCTFYCADSEKQKLEEQYRQEDQRERELRLLLLQQSNV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001271687

ORF Size:	1206 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001271687.1 , NP_001258616.1
RefSeq Size:	5210 bp
RefSeq ORF:	1209 bp
Locus ID:	16516
Cytogenetics:	16 55.86 cM
MW:	46 kDa
Gene Summary:	Inward rectifier potassium channels are characterized by a greater tendency to allow potassium to flow into the cell rather than out of it. Their voltage dependence is regulated by the concentration of extracellular potassium; as external potassium is raised, the voltage range of the channel opening shifts to more positive voltages. The inward rectification is mainly due to the blockage of outward current by internal magnesium.[UniProtKB/Swiss-Prot Function]