

Product datasheet for MR229760

Tead2 (NM_001285500) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tead2 (NM_001285500) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tead2
Synonyms:	Etdf; ETF; TEAD-2; TEF-4; TEF4
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
ORF Nucleotide Sequence:	>MR229760 ORF sequence Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGGGGATCCCCGGACTGGGGCCCCCTGGATGATGGCGGGGGCTGGACAGGTAGCGAGGAAGGCAGCG
AAGAGGGCACCGGAGGCAGCGAGGGCGTTGGAGGTGATGGCAGCCCTGATGCAGAAGGCCGGAACGAACT
CATTGCCCGTTACATCAAGCTGAGGACAGGGAAGACGAGAACCGCAAAGCAGGTCTCCAGCCATATTCAG
GTTTTGGCTCGAAGGAAATCGAGAGAAATTCAGTCCAAGCTGAAGGACCAAGTCTCCAAGGACAAGGCCCT
TCCAGACGATGGCCACCATGTCTTCGGCACAGCTCATCTCGCCCTTCCCTCCAGGCCAAGCTGGGCC
TTCTGGCCCTCAGGCCACTGAGCTTTTCCAGTTCTGGTCAGGGAGCTCTGGGCCACCATGGAATGTTCCA
GACGTGAAGCCCTTCTCACAGGCACCGTTCTCCGTGTCAGTACGCCCCAGCCTCTGACCTACCAGGGT
ACGAGCCGCCCCAGCCCTCTCACCCCTGCCCCACCCCGTCCGTCTCCCCAGCCTGGCAGGCTCGGGC
CCTGGGCACTGCCCGCTGCAGCTGATAGAGTTCTCAGCGTTTGTGGAACCGCCAGACGCAGTTGACTCG
TTCCAGAGGCATCTGTTGTCCACATCAGTCAGCAGTGTCCAGCCCTGGAGCACCACCCCTAGAGAGTG
TGGACGTGCGGCAGATCTACGACAAATCCCTGAGAAGAAGGGCGGCCTCCGCGAGCTGTATGACCGAGG
GCCACCACATGCCTTCTCCTCGTCAAGTTCTGGCGGACCTGAAGTGGGCCCCAGTCCGAGGAGGCA
GGGAGCAGCGGAGGTGGCGGTGGCTTCTATGGAGTGAGCAGCAGTATGAGAGCCGGGAGCTCATGACAC
TCACCTGCTCCTCAAGGTCTGCTCCTTTGGCAAGCAAGTGGTAGAGAAGGTGGAGACGGAACGGGCCCA
GCTGGAGGACGGGCGCTTTGTGTACCGTCTGCTGCGCTCTCCATGTGTGAGTACCTGGTTAATTTCTG
CACAAGCTCCGTGAGCTGCCTGAACGCTACATGATGAACAGTGTCTGGAGAAGTTCACCATCCTCCAGG
TTGTGACAAACAGGGACACTCAGGAACTGCTGCTGTGTAAGTCTGCTACGTCTTTGAAGTCTCCACCAGTGA
ACGAGGAGCCAGTACCACATCTACCGCTGGTCAGGGAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Protein Sequence: >MR229760 protein sequence
Red=Cloning site Green=Tags(s)

MGDPRTGAPLDDGGGWTGSEEGSEEGTGGSEGVGGDGPDAEGRNELIARYIKLRTGKTRTRKQVSSHIO
 VLARRKSREIQSKLKDQVSKDKAFQTMATMSSAQLISAPSLQAKLGPSPQATELFFQWSSGSPWNVP
 DVKPFSSQAPFSVSLTPPASDLPGYEPAL SPLPPPAPSPPAWQARALGTARLQLIEFSAFVEPPDAVDS
 FQRHLFVHISQQCPSPGAPPLESVDRQIYDKFPEKKGGLREL YDRGPPHAFLLVKFWADLNWGPSAEEA
 GSSGGGGGFYGVSSQYESREMLTLTCCSSKVCSEFGKQVVEK VETERAQLEDGRFVYRLLRSPMCEYLVNFL
 HKLRQLPERYMMNSVLENFTILQVVTNRDTQELLLCTAYVFEVSTSERGAQYHIYRLVRD

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

ACCN: NM_001285500

ORF Size: 1233 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001285500.2](#)

RefSeq Size: 2010 bp

RefSeq ORF: 1233 bp


Locus ID: 21677

Cytogenetics: 7 29.19 cM

MW: 45.1 kDa

Gene Summary: Transcription factor which plays a key role in the Hippo signaling pathway, a pathway involved in organ size control and tumor suppression by restricting proliferation and promoting apoptosis. The core of this pathway is composed of a kinase cascade wherein MST1/MST2, in complex with its regulatory protein SAV1, phosphorylates and activates LATS1/2 in complex with its regulatory protein MOB1, which in turn phosphorylates and inactivates YAP1 oncoprotein and WWTR1/TAZ. Acts by mediating gene expression of YAP1 and WWTR1/TAZ, thereby regulating cell proliferation, migration and epithelial mesenchymal transition (EMT) induction (By similarity). Binds to the SPH and GT-IIC 'enhancers' (5'-GTGGAATGT-3'). May be involved in the gene regulation of neural development. Binds to the M-CAT motif.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR229760