

## Product datasheet for **MR229931**

### Atf2 (NM\_001284376) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Atf2 (NM_001284376) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Atf2
Synonyms:	Atf-2; CRE-BP; Creb2; D18875; D130078H02Rik; mXBP; Tg(Gzma-Klra1)7Wum
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
ORF Nucleotide Sequence:	>MR229931 representing NM_001284376 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGCAGTCCCCTGGGTGTGGCCAGATCAGACTCCAACGCCAACAAGATTCCTAAAAACTGTGAAGAAG  
TGGTTTTGTTCAATGAGTTGGCAAGTCCATTTGAAAATGAATTCAGAAGGCTTCCGAAGATGACATTAA  
AAAAATGCCTCTAGATTTGTCCCCTCTTGCAACACCCATCATAAGAAGCAAAATTGAGGAGCCTTCTGTT  
GTAGAAACAACCTACCAGGACAGCCCTTACCTCACCCGAGTCGACTACCAAGTATGAAAAGGAAGTAC  
CATTGGCACAACCTGCACAGCCACATCAGCTATCGTTTCGTCACGATCATTACAGGTTCCCAATGTGCT  
GCTCACAAGTTCTGACTCAAGTGAATTAATTCAACAAGCAGTACCTTCACCAACCTCAAGTACTGTAATC  
ACCCAGGCACCATCCTCTAACAGGCCAATTGTTCTGTACCAGGCCCATTCCTCTTCTTACATCTTC  
CTAATGGACAACCATGCCCGTTGCTATTCCTGCATCAATTACAAGTTCTAATGTGCATGTTCCAGCTGC  
AGTCCCCTTGCAGTTGTTCCAGGAATCCCAGGCCCTTCTCCTCAACCAGTCCAGTCAGAAGCAAAA  
ATGAGATTAAGAGCTGCTTTGACCCAGCAACACCCTCCAGTTACCAATGGTGATACTGTAAGGCCATG  
GCAGTGGATTGGTTAGGACTCAGTCAGAAGAGTCTCGCCACAGTCTTGCAGCAGCCACCCACCTCCAC  
TACAGAAACTCCGGCTTCTCCAGCTCACACAACCTCCTCAGACCCAAAATACAAGTGGCCGTCGAAGAAGA  
GCAGCTAATGAAGATCCTGATGAGAAAAGGAGGAAGTTTCTAGAACGAAATAGAGCAGCAGCTTCAAGAT  
GCCGACAAAAAGGAAAGTGTGGTTTCAGTCTTAGAGAAGAAAGCAGAAGACTTGAAGTCACTAAATGG  
CCAGCTGCAGAGCGAAGTCAACCTGCTGAGAAAATGAAGTGGCCAGCTGAAACAGCTTCTTCTGGCTCAT  
AAAGATTGCCCTGTAAGTCCATGCAGAAGAAGTCTGGCTATCATACTGCTGATAAAGATGACAGTTCAG  
AAGACCTTTCTGTGCCAAGCAGTCCACATACAGAAGCGATCCAGCACAGCTCTGTCAGCACATCCAATGG  
AGTCAGTTCAACATCAAAAGCAGAAGCTGTAGCCACTTCACTCCTCACCCAGATGGCGGACCAGAGCAGC  
GAGCCTGCACTTTCACAGATTGCATGGCTCCTCCCTCCAGGCACAGCCCTCAGGAAGT

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

**Protein Sequence:** >MR229931 representing NM\_001284376  
Red=Cloning site Green=Tags(s)


MHCPWVWPDPQTPTRFLKNCEEVGLFNELASPFENEFFKASEDDIKMPLDLSPLATPIIRSKIEEPSV  
 VETTHQDSPLPHPESTTSDEKEVPLAQT AQPTSAIVRPASLQVPNVLLTSSDSSVVIQQAVPSPTSSTVI  
 TQAPSSNRPIVPVPGPFLLLLHLPNGQTMPVAIPASITSSNVHVAAPLALVPGIPGPSSPQPVOSEAK  
 MRLKAAL TQQHPVVTNGD TVKGHGSGLVRTQSEESRQSLQQPATSTTETPASPAHTTPQTQNTSGRRRR  
 AANEDPDEKRRKFLERNRAAASRCRQKRKVVWQSLKKAEDLSSLNGQLQSEVTLRNEVAQLKQLLLAH  
 KDCPVTAMQKKSgyHTADKDDSSEDLSVPSSPHTeAIQHSSVSTsNGVSSTSKAEAVATSVLTQMADQST  
 EPALSQIVMAPPsQAQPSGS

TRRPLEQKLISEEDLAANDILDYKDDDDKV


**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**



**ACCN:**

NM\_001284376

<b>ORF Size:</b>	1320 bp
<b>OTI Disclaimer:</b>	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. <a href="#">More info</a>
<b>OTI Annotation:</b>	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_001284376.1</a> , <a href="#">NP_001271305.1</a>
<b>RefSeq Size:</b>	4182 bp
<b>RefSeq ORF:</b>	1323 bp
<b>Locus ID:</b>	11909
<b>Cytogenetics:</b>	2 C3
<b>MW:</b>	47.7 kDa
<b>Gene Summary:</b>	Transcriptional activator which regulates the transcription of various genes, including those involved in anti-apoptosis, cell growth, and DNA damage response. Dependent on its binding partner, binds to CRE (cAMP response element) consensus sequences (5'-TGACGTCA-3') or to AP-1 (activator protein 1) consensus sequences (5'-TGACTCA-3'). In the nucleus, contributes to global transcription and the DNA damage response, in addition to specific transcriptional activities that are related to cell development, proliferation and death. In the cytoplasm, interacts with and perturbs HK1- and VDAC1-containing complexes at the mitochondrial outer membrane, thereby impairing mitochondrial membrane potential, inducing mitochondrial leakage and promoting cell death. The phosphorylated form (mediated by ATM) plays a role in the DNA damage response and is involved in the ionizing radiation (IR)-induced S phase checkpoint control and in the recruitment of the MRN complex into the IR-induced foci (IRIF). Exhibits histone acetyltransferase (HAT) activity which specifically acetylates histones H2B and H4 in vitro. In concert with CUL3 and RBX1, promotes the degradation of KAT5 thereby attenuating its ability to acetylate and activate ATM. Can elicit oncogenic or tumor suppressor activities depending on the tissue or cell type (By similarity).[UniProtKB/Swiss-Prot Function]