

Product datasheet for **MR230198**

Gpcpd1 (NM_001291051) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gpcpd1 (NM_001291051) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Gpcpd1
Synonyms:	2310004G06Rik; 2310032D16Rik; AU015213; Gde5; mKIAA1434; Prei4
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230198 representing NM_001291051
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATG**TCCAACAGCCTGGAGATATCCTTAATAAGTGACAATGAGTTCAAGTGCAGGCACTCACAGCCAGAAT**
GTGGGTATGGCTTACAGCCCGATCGTTGGACAGAGTACAGCATACAGACAATGGAACCAGATAATCTGGA
GCTCATCTTTGACTTTTTTGAGGAAGATCTCAGT**GAGCATGTAGTTCAGGGT**GATGTTCTTCTG**GACAC**
GTGGGCACAGCATGCCTCCTGTCTTCTACCATTGCTGAGAGTGGAAGAAGCGCTGGAATCCTTACTCTTC
 CCATCATGAGCAGAAATCCAGAAAACTATAGGCAAAGT**CAGAGTTGATTTTATCATCATCAAGCCATT**
 ACCTGGATATAGTTGTTCTATGCAGTCTTCATTTTCCAAGTATTGGAAACCAAGAATACCATTGGACGTT
 GGACATCGTGGTGCAGGAACTCAACAACGACTGCCAAGCTAGCTAAAGTACAGGAAAATACTATCGCTT
 CTTTAAGAAATGCTGCCAGT**CATGGCGCAGCATT**TGTAGAAATTTGATGTCCACCTTCAAAGGACTTTGT
 GCCCGTGGTGTATCATGACCTCACCTGCTGTCTGACCATGAAGAGGAAAATGAAGCTGATCCAGTTGAA
 TTGTTTGAATCCAGTAAAAGAATTAACATTTGACCAACTCCAGTTATTGAAGCTTTCTCATGTGACTG
 CATAAAAACCAAGACCGGAAACAATCTTTGTATGAGGAGGAAAATTTCTTTTCTGAAAATCAGCCATT
 TCCTTCTCTTAAGATGGTTTTAGAAATCATTGCCAGAAAATGATAGGATTTAATATAGAAAATAAAATGGATT
 TGCCAACACAGGGATGGAGTATGGGATGGCAACTTATCAACATATTTTGTATGAATGTGTTTTTGGATA
 TAATTTTAAAACTGTTTTAGAAAATCTGGGAAGAGGAGAATAGTGTTTTCTTCTTTTGTATGCAGATAT
 TTGTACAATGGTTCGGCAGAAGCAGAACAAATATCCCATATTATTTT**GACCCAAGGAAAGTCTGATATT**
 TACCCGAAC**T**CATGGACCTCAGATCTCGGACAACACCCATTGCAATGAGTTTTGCACAGTTT**GAAAATA**
 TTTTGGGATAAATGCCATACTGAAGACCTCCTTAGAAAACCCATCCTATGTCCAAGGCAAAAAGCTAA
 GGGATTGGT**CATATTCTGCTGGGGT**GATGATACCAACGATCCTGAAAACAGAAAGGAAACTGAAGGAATTT
 GGAGTAAATGGTCTAATATATGATAGGATATATGATTGGATGCCTGAACAGCCAAATATATTCCAAGTGG
 AGCAGTTGGAGCGCTGAAGCAAGAATTGCCAGAGCTTAAGAAGTGTGTGTCCCACTGTAGCCACTT
 CATCCCTTCTTCTTTCTGTGTGGAGCCTGATATCCACGTGGATGCCAACGGCATTGATAGTGTGGAGAAC
 GCT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230198 representing NM_001291051
 Red=Cloning site Green=Tags(s)

MSNSLEISLISDNEFKCRHSQPECGYGLQPDRWTEYSIQTMEDPNLELIFDFFEEEDLSEHVVGQDVLPGH
 VGTACLLSSTIAESGRSAGILTLPIMSRNSRKTIGKVRVDFIIKPLPGYSCSMQSSFSKYWKPRIPLDV
 GHRGAGNSTTTAKLAKVQENTIASLRNAASHGAAFVEFDVHLSKDFVPVVYHDLTCCLTMKRKYEADPVE
 LFEIPVKELTFDQLQLLKLSHVTALKTKDRKQSLYEEENFFSENQPFPSLKMVLESLPENVGFNIEIKWI
 CQHRDGVWDGNLSTYFDMNVFLDIILKTVLENSGKRRIVFSFDADICTMVRQKQNKYPILFLTQGKSDI
 YPELMDLRSRTTPIAMSFQFENILGINAHTEDLLRNPSYVQEAkakGLVIFCWGDDTNDPENRRKLKEF
 GVNGLIYDRIYDWMPQPNIFQVEQLERLQELPELKNCLCPTVSHFIPSSFCVPEPIHVDANGIDSVEN
 A

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI


Cloning Scheme:

Cloning sites used for ORF Shuttling:



* The last codon before the Stop codon of the ORF

Plasmid Map:



ACCN: NM_001291051

ORF Size: 1473 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291051.1](#), [NP_001277980.1](#)

RefSeq Size: 3622 bp

RefSeq ORF: 1476 bp

Locus ID: 74182

UniProt ID: [Q8C0L9](#)

Cytogenetics: 2

MW: 56.4 kDa

Gene Summary: May be involved in the negative regulation of skeletal muscle differentiation, independently of its glycerophosphocholine phosphodiesterase activity.[UniProtKB/Swiss-Prot Function]