

Product datasheet for **MR230210**

Mef2a (NM_001291192) Mouse Tagged ORF Clone

Product data:

| | |
|---------------------------|---|
| Product Type: | Expression Plasmids |
| Product Name: | Mef2a (NM_001291192) Mouse Tagged ORF Clone |
| Tag: | Myc-DDK |
| Symbol: | Mef2a |
| Synonyms: | A430079H05Rik |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |



[View online »](#)

ORF Nucleotide Sequence:

>MR230210 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGGGCGAAAGAAGATACAAATCACACGCATAATGGATGAGAGGAACCGACAGTTACTTTTACAAGC
 GAAAGTTTGGATTGATGAAGAAAAGCCTATGAACTCAGTGTGCTCTGTGACTGTGAAATAGCACTCATCAT
 TTTCAACAGCTCTAATAAGCTGTTTCAGTACGCTAGCACTGACATGGACAAAGTCTTCTCAAATACACT
 GAGTATAACGAGCCTCATGAAAGCAGGACGAACTCGGATATCGTTGAGGCTCTGAACAAGAAGGAACACA
 GAGGGTGCAGACGCCAGACCCTGATACTTCATATGTGCTAACTCCACATACAGAAGAAAAATATAAAAA
 AATTAATGAGGAATTTGATAATATGATGCGGAATCATAAAATCGCACCTGGCTTGCCACCTCAGAACTTC
 TCAATGTCTGTCACAGTCCCTGTGACCAGCCCTAATGCTTTGTCGTACACTAACCCAGGGAGTTCACCTCG
 TGTACCCTCTTTGGCAGCCAGCTCAACATTAGCAGATTGAGCATGCTGTCTCCACCCCGCCAGCCAGCT
 ACATAGAAATGTGTCCCAGGAGCTCCGCAGAGACCACCAAGTACCGGCAGTGCAAGTGGGATGTTGAGC
 ACTACAGACCTCACGGTGCACAAATGGAGCTGGAATAGTCTGTGGGGAATGGATTTGTAAACTCAAGGG
 CCTCTCCAAATTTGATTGGAAATACTGGTGCAAAATAGTTTAGGCAAAGTCATGCCTACAAAGTCTCCCC
 GCCACCAGGTGGTGGCAGTCTTGAATGAACAGTCGGAACCAGATCTACGAGTTGTCAATCCCCATCA
 AGCAAGGGCATGATGCCCCACTTAATGCCAAAGGATAAGCAGTTCTCAAGCCACTCAACCTCTTGCCA
 CGCCTGTAGTGTCTGTGACAACCTCCGAGCTTGCCTCCACAGGGACTAGTGTACTCAGCAATGCCGACAGC
 CTACAACACTGATTACTCACTGACCAGCGCTGACCTGTCTGCCCTGCAAGGCTTTACTTCCCCTGGAATG
 CTGTCTCTGGGGCAGGCTTCAGCCTGGCAGCACCCTAGGACAAGCAGCCCTCAGCTCTCTGGTGG
 TGGAGGGCAGTTATCTCAGGGTTCAAATTTGTCCATTAAACACCAACCAGAACATCAACATTAAGTCAGA
 GCCAATTTACCTCCCCGTGATAGAATGACCCCATCAGGCTTCCAACAGCAGCAGCAGCCACAACAG
 CAGCCACCGCCACAACCTCCACAGCCCAACCCCGGCAGGAAATGGGGCGCTCCCCTGTGGACAGTCTGA
 GCAGCTCCAGCAGCTCCTATGATGGCAGTACCAGGGAGGATCCACGGGGTACTTCCATTCTCAATTGT
 GCTTGGCCGACCTCCAAATACTGAGGACAGAGAAAGCCCTTCTGTAAAACGCATGAGGATGGACACATGG
 GTGACC

ACGGTACGGCGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230210 protein sequence
 Red=Cloning site Green=Tags(s)

MGRKKIQITRIMDERNRQVFTTKRKFGLMKKAYELSVLCDCEIALIIFNSSNKL FQYASTDMDKVL LKYT
 EYNPHESTRNSDIVEALNKKEHRGCDSPDPDTSYVLTPHTEEKYKKINEEFDNMMRNHKAIPGLPPQNF
 SMSVTVPVTPSPNALSYTNPGSSLVSPSLAASSTLADSSMLSPPPATLHRNVSPGAPQRPPSTGSASGMLS
 TTDLTVPNGAGNSPVNGFVNSRASPNIIGNTGANS LGKVMPTKSPPPPGGSLGMNSRKPDLRVVIPP
 SKGMMPPLNAQRISSSQATQPLATPVVSVTTPSLPPQGLVYSAMPTAYNTDYSLSADLSALQGFTSPGM
 LSLGQASAWQHHLGQAALSSLVAGGQLSQGSNLSINTNQINIKSEPI SPPDRMTPSGFQQQQQPQQ
 QPPPQPQPQPRQEMGRSPVDSLSSSSSYDGSREDPRGDFHSPIVLGRPPNTEDRESVSKRMRMDTW
 VT

TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-Mlul

Cloning Scheme:


ACCN: NM_001291192

ORF Size: 1479 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291192.1](#), [NP_001278121.1](#)

RefSeq Size: 5554 bp

RefSeq ORF: 1479 bp

Locus ID: 17258


UniProt ID: [Q60929](#)

Cytogenetics: 7 36.72 cM

MW: 53 kDa

Gene Summary: Transcriptional activator which binds specifically to the MEF2 element, 5'-YTA[AT](4)TAR-3', found in numerous muscle-specific genes. Also involved in the activation of numerous growth factor- and stress-induced genes. Mediates cellular functions not only in skeletal and cardiac muscle development, but also in neuronal differentiation and survival. Plays diverse roles in the control of cell growth, survival and apoptosis via p38 MAPK signaling in muscle-specific and/or growth factor-related transcription. In cerebellar granule neurons, phosphorylated and sumoylated MEF2A represses transcription of NUR77 promoting synaptic differentiation. Associates with chromatin to the ZNF16 promoter (By similarity).[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR230210