

Product datasheet for **MR230218A1V**

Mouse Wipfl (NM_001289723) AAV Particle

Product data:

Product Type:	AAV Particles
Tag:	Myc-DDK
Symbol:	Wipfl
Synonyms:	A1115543; D2Ertd120e; Waspip; WIP
Mammalian Cell	None
Selection:	
Vector:	pAAV-AC-Myc-DDK (PS100089)



ORF Nucleotide Sequence: >MR230218 representing NM_001289723
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

[illegible]

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230218 representing NM_001289723
 Red=Cloning site Green=Tags(s)

MPVPPPPAPPPPTFALANTEKPTLNKTEQAGRNALLSDISKGKKLKKTVTNDRSAPILDKPKGAGASAG
GYGGGGGGGGGGGGGGGGSGGNFGGGGPPGLGGLFQAGMPKL RSTANRDNDSGGRPPILPPGGRATSAK
PFSPSPGPRFPAPSPGHRSGPPEPPRNRMPPPRPDVGSKPDSLPPVPNTPRPVSSLHNRGSPAGLGA
PRPPFPNGRAAFGAGSARQNPSGSSSPFRPPPLPTPSRALDDKPPPPPPVGNRPSMHREAVPPPPSQ
TSKPPVPSTPRPGLGSQAPPPPPPSRGPPLPPASNDEIPRLPQRNLSLTSSAPPLSPSGRSGPLPPP
PSRPPPPVRDPPGRSGPLPPPPPINRNGSTARALPATPQLPSRSGMDSRSGRPPPLPPDRPGAGAPP
PPSTSVRNGFQDSSCEDEWESRFYFHIPIDLPPPEYVPTTKTYPSKLARNESRSGSNRRERGAPPLP
IPR

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Species: Mouse

Serotype: AAV-2

ACCN: NM_001289723

ORF Size: 1479 bp

Buffer: PBS with 0.001% Pluronic F68

Stability:	AAV is stable for 1 year when stored at -80°C (long-term storage) or 2-3 weeks when stored at -20°C (short-term storage). Thaw the vial of AAV on ice prior to use and keep it on ice during the experiment. Thawed AAV can be stored at 4°C for 1-2 weeks. Whenever possible, particles should be aliquoted into single use portions to avoid repeated freeze/thaw cycles. Please aliquot at least 10ul per tube and use low protein binding tubes to avoid loss of virus.
RefSeq:	NM_001289723.1 , NP_001276652.1
RefSeq Size:	4451 bp
RefSeq ORF:	1482 bp
Locus ID:	215280
UniProt ID:	Q8K1I7
Cytogenetics:	2 43.68 cM
MW:	50.1 kDa