

Product datasheet for **MR230312**

Cstf2 (NM_001290399) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cstf2 (NM_001290399) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cstf2
Synonyms:	C630034J23Rik; Cstf64; mir-3112
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR230312 ORF sequence
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGCGGGTTTGCCAGTGAGAGACCCTGCGGTGGATCGTTCCTGCGCTCCGCTTCGTTGGGAAACATCC
 CTTATGAAGCTACTGAAGAACAACCTGAAGGACATTTCTCTGAGGTTGGGCCTGTTGTTAGTTTCAGATT
 GGTATACGATAGGGAGACAGAAAACAAAGGGCTATGGCTTCTGCGAGTACCAAGACCAAGAGACAGCA
 CTTAGTGCCATGCGGAACCTGAATGGGCGTGAATTCAGCGGGAGAGCACTTCGAGTGGATAATGCTGCCA
 GTGAAAAGAACAAGGAGGAACCTGAAGAGCCTTGGCACTGGTGCACCTGTCATAGAGTCACCTTATGGAGA
 AAGCATCAGCCCTGAGGATGCCCTGAATCCATTAGCAAAGCAGTTGCCAGTCTTCCACCAGAACAGATG
 TTTGAGTTGATGAAACAATGAAGCTCTGTGTCCAGAATAGTCCCAGGAAGCACGAAACATGTTGCTTC
 AGAATCCACAGCTGGCTTATGCTTTGTTGCAAGCACAGGTAGTGATGAGAATTGGGATCCTGAGATTGC
 GCTGAAAATTCTGCATCGCCAGACAAATATCCCAACGCTGATTTACAGCAACCCTCAGCCAGTTCATGTT
 GCTGGGCCTGGCTCAGGACCCAATGTTTCAATGAACCAGCAGAATCCGCGAGGCCCTCAGGCTCAATCTT
 TGGGTGGAATGCATGTCATGGCGCACCTCCTATGATGCAAGCTTCTATGCCAGGTGGAGTTCACGCCCC
 AGTGCAGATGGCAGCTGCTGTAGGAGGACCTGGCCCTGGTTTCTTAGCCCTGCAGGAGTTATGCAAGCC
 CAAGTTGGAATGCAAGGAGCTGGACCACTGCCATGGAACGAGGACAAGTGCCAATGCAAGACCCAGAG
 CAGCTATGCAGCGGGGAGCCTTGCTACCAACGTCCCACTCCTCGTGGTCTTTTAGGTGATGCTCCCAA
 TGACCCGCGGGGAGGCACCTTAACTGACTGTAACGGAGACGTAGAACCTAGAGCTTACCTGGGACCACCA
 CCACCACCTCATCAGGGCCACCCATGCACCATGTTCTGGTTCATGAAGGTCGTGGACCCCTCCACATG
 ACATGAGAGGAGGACCCCTAGCTGAGCCGAGACCTTAAATGGCAGAACCAAGAGGACCCATGCTAGATCA
 GAGGGGCCACCTCTGGATGCCAGAGGTGAAGAGATCCCGAGGATTAGATGCACGGGGGATGGAGGCC
 AGAGCCATGGAAGCAAGAGGACTAGATGCCAGAGGACTGGAAGCCCGTCCATGGAGGCCGTGCCATGG
 AGGCCCGTCCATGGAAGCCCGTCCATGGAAGCCCGTCCATGGAGGCTCGTCCATGGAGGCCAGAGG
 CATGGACACCAGAGGCCAGTACCTGGACCTAGAGGACCTATGCCTAGTGAATCCAAGTCCCAATCCA
 ATGAACATGGGGCAGTTGTTCCCGAGGATCCAGACAGGTAATGCTTGTAGCATATACA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230312 protein sequence
 Red=Cloning site Green=Tags(s)

MAGLPVRDPAVDRSLRSVFGNIPYEATEEQKIDIFSEVGPVVSFRLVYDRETKPKGYGFCEYQDQETA
 LSAMRNLNGREFSGRALRVDNAASEKNKEELKSLGTGAPVIESPYGESISPEDAPESISKAVASLPPEQM
 FELMKQMKLCVQNSPQEARNMLLQNPQLAYALLQAQVVMRIVDPEIALKILHRQTNIPTLISGNPQPVHV
 AGPGSGPNVSMNQNPQAPQAQSLGGMHVNGAPPMQASMPGGVPAPVQMAAAVGGPGPSLAPAGVMQA
 QVGMQGAGVPMERGQVPMQDPRAAMQRGALPTNVPTPRGLLDAPNDPRGGTLMVTGDVEPRAYLGPP
 PPPHQGPPMHVPGHEGRGPPPHDMRGGPLAEPRLMAEPRGPMLDQRGPPLDARGGRDPRGLDARGMEA
 RAMEARGLDARGLEARAMEARAMEARAMEARAMEARAMEARAMEARAMEARGMDTRGPVPGPRGPMPSGIQGNP
 MNMGAVVPQGSRVMLVAYT

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:

ACCN: NM_001290399

ORF Size: 1533 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)
OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290399.1](#), [NP_001277328.1](#)
RefSeq Size: 2772 bp

RefSeq ORF: 1533 bp


Locus ID: 108062

UniProt ID: [Q8BIQ5](#)
Cytogenetics: X E3

MW: 54.1 kDa

Gene Summary: One of the multiple factors required for polyadenylation and 3'-end cleavage of mammalian pre-mRNAs. This subunit is directly involved in the binding to pre-mRNAs (By similarity). [UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR230312