

Product datasheet for **MR230458**

Nf2 (NM_001252253) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nf2 (NM_001252253) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nf2
Synonyms:	merlin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230458 representing NM_001252253
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGAAATGGAAGGGAAGGACCTGTTTGTATTTGGTGTGCCGACTGGGGCTTCGGAAACCTGGTTCT
 TTGGACTGCAGTATACAATCAAGGACACGGTGGCTGGCTCAAAATGGACAAGAAGGTGTTGGATCATGA
 TGTTTCGAAGGAAGAACCAGTTACCTTTCACTTCTGGCCAAATTTTATCCTGAAATGCTGAGGAGGAG
 CTAGTTCAAGAGATCACGCAACACTTATTTTTCTTACAGGTGAAGAAGCAGATTTTGGATGAAAAGGTCT
 ACTGCCCTCCCGAGGCGTCCGTGCTCTTGGCGTCATATGCTGTCCAGGCCAAGTATGGCGACTATGACCC
 CTCTGTGCACAAGCGGGATTTTAGCCCAAGAGGAATTGCTCCCGAAAAGGGTGATAAATCTCTATCAG
 ATGACTCCGAAATGTGGGAGGAGAGAATTACGGCTTGGTATGCGGAGCACCGGGGCAGAGCCAGGGATG
 AAGCTGAAATGGAGTATTTGAAGATAGCTCAGGACCTGGAGATGTATGGTGTGAACACTTTACAATCCG
 GAATAAAAAGGGCACAGAGTTGCTGCTTGGAGTGGATGCTCTTGGGCTTCATATCTATGACCTGAGAAC
 AGGCTGACCCCAAGATCTCCTTCCCATGGAATGAAATCCGAAACATCTCTACAGCGACAAGGAGTTTA
 CTATTAACCCTGGATAAGAAAATTGATGCTTCAAATTTAACTCCTCAAAGCTTCGTGTTAATAAGCT
 GATTTCTCAGCTATGTATTGGGAACCATGACCTATTTATGAGGCGACGAAAAGCTGACTCTTTAGAAGTT
 CAGCAGATGAAAGCCAGGCCAGGGAAGAGAAGGCTAGAAAGCAGATGAAAGGCCAGCGGCTGGCTCGAG
 AGAAGCAGATGCGGGAGGAGGCCAGCGTACAAGAGATGAGTTAGAGAGGAGGCTCCTGCAGATGAAAGA
 AGAAGCAACGATGGCCATGAAGCTCTGATGCGCTCTGAGGAGACAGCTGATCTGTTGGCTGAAAAGGCT
 CAGATCACAGAGGAGGAGGCCAAGCTTTTGGCACAAAAGGCTGCAGAGGCTGAGCAAGAGATGCAGCGAA
 TCAAGGCCACGGCCATTCGGACAGAGGAGGAGAAGCGCCTGATGGAGCAGAAGGTCTGGAGGCTGAAGT
 GCTGGCATTGAAGATGGCTGAGGAGTCAAGAGGAGGAGGCCAAGGAGGCTGATCAGTTAAAGCAAGACTTG
 CAAGAAGCCAGAGAAGCAGAGCGAAGAGCCAAGCAGAAGCTCTTAGAAATCGCCACCAAGCCACCTATC
 CACCCATGAACCAATTCCACCACCTGCCTCCTGACATACCGAGCTTCGACATTATTGCTGACAGCTT
 GTCATTTGACTTCAAGGATACGGACATGAAGCGACTTTCCATGGAGATAGAGAAAAGAAAAGTGGAGTAC
 ATGGAGAAGCAAGCACCTGCAGGAGCAGCTCAACGAGCTCAAGACGGAGATCGAGGCTTGAACCTCA
 AAGAGCGGAGACGGCTTGGACGCTCTACACAGCGAGAGCTCAGACAGAGGCGGCCCCAGCAGCAAGCA
 TAATACCATTA AAAAGGTACCTGAAATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230458 representing NM_001252253
 Red=Cloning site Green=Tags(s)

MKWKGKDLFDLVCRTLGLRETWFFGLQYTIKDTVAVLKMDDKVLHDVSKKEEPVTFHFLAKFYPENAEIE
 LVQEITQHLFFLQVKKQILDEKVCPEASVLLASYAVQAKYGDYDPSVHKRGFLAQEELLPKRVINLYQ
 MTPPEMWEERITAWYAEHRGRARDEAEEMEYLKIAQDLEMYGVNYFTIRNKKGTLLLLGVDALGLHIYDPEN
 RLTPKISFPWNEIRNISYSDFEFTIKPLDKKIDVFKFNSSKLRVNLILQLCIGNHDLFMRRRKADSLLEV
 QQMKAQAREEKARKQMERQRLAREKQMRREEAERTRDELEERRLLQMKEEATMANEALMRSEETADLLAEKA
 QITEEEAKLLAQAAAEQEMQRIKATAIRTEEEKRLMEQKVLAEVLALKMAEESERRAKEADQLKQDL
 QEAREAEERRAKQLLEIATKPTYPPMNPIPPPLPPDIPSPFDIIADSLSFDFKDTDMKRLSMEIEKEKVEY
 MEKSKHLQEQLNELKTEIEALKKERETALDVLHSESSDRGGPSSKHNTIKKVPPEM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

ACCN:	NM_001252253
ORF Size:	1638 bp
OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001252253.1 , NP_001239182.1
RefSeq Size:	1846 bp
RefSeq ORF:	1641 bp
Locus ID:	18016
UniProt ID:	P46662
Cytogenetics:	11 3.0 cM
MW:	64.6 kDa
Gene Summary:	Probable regulator of the Hippo/SWH (Sav/Wts/Hpo) signaling pathway, a signaling pathway that plays a pivotal role in tumor suppression by restricting proliferation and promoting apoptosis. Along with WWC1 can synergistically induce the phosphorylation of LATS1 and LATS2 and can probably function in the regulation of the Hippo/SWH (Sav/Wts/Hpo) signaling pathway. May act as a membrane stabilizing protein. May inhibit PI3 kinase by binding to AGAP2 and impairing its stimulating activity. Suppresses cell proliferation and tumorigenesis by inhibiting the CUL4A-RBX1-DDB1-VprBP/DCAF1 E3 ubiquitin-protein ligase complex (By similarity). Plays a role in lens development and is required for complete fiber cell terminal differentiation, maintenance of cell polarity and separation of the lens vesicle from the corneal epithelium.[UniProtKB/Swiss-Prot Function]