

Product datasheet for **MR230590**

Nf2 (NM_001252252) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nf2 (NM_001252252) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nf2
Synonyms:	merlin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230590 representing NM_001252252
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGCCGGAGCCATCGTTCTCGCATGAGCTTCACTCAAGAGGAAGCAGCCCAAGACATTACCGG
 TCGGATCGTCACCATGGACGCCGAGATGGAGTTCAACTGCGAGATGAAATGGAAGGGAAGGACCTGTT
 TGATTTGGTGTGCCGGCACTGGGGCTTCGGGAAACCTGGTTCTTTGGACTGCAGTATACAATCAAGGAC
 ACGGTGGCTGGCTCAAATGGACAAGAAGGTGTTGGATCATGATGTTTCGAAGGAAGAACCAGTTACCT
 TTCCTTCTGGCCAAATTTTATCCTGAAAATGCTGAGGAGGAGCTAGTTCAAGAGATCACGCAACACTT
 ATTTTTCTTACAGGTGAAGAAGCAGATTTGGATGAAAAGGTCTACTGCCCTCCCGAGGCGTCCGTGCTC
 TTGGCGTCATATGCTGTCCAGGCCAAGTATGGCGACTATGACCCCTCTGTGCACAAGCGGGGATTTTATG
 CCCAAGAGGAATTGCTCCGAAAAGGGTGATAAATCTCTATCAGATGACTCCGGAATGTGGGAGGAGAG
 AATTACGGCTTGGTATGCGGAGCACCGGGCAGAGCCAGGGATGAAGCTGAAATGGAGTATTTGAAGATA
 GCTCAGGACCTGGAGATGTATGGTGTGAACTACTTTACAATCCGGAATAAAAAGGACACAGAGTTGCTGC
 TTGGAGTGGATGCTCTTGGGCTTCATATCTATGACCCGTGAGAACAGGCTGACCCCAAGATCTCCTTCCC
 ATGGAATGAAATCCGAAACATCTCTACAGCGACAAGGAGTTTACTATTAACCCTGGATAAGAAAATT
 GATGCTTCAAATTTAACTCCTCAAAGCTTCGTGTTAATAAGCTGATTCTTCAGCTATGTATTGGGAACC
 ATGACCTATTTATGAGGCGACGGAAGCTGACTCTTTAGAAGTTCAGCAGATGAAAGCCAGGCCAGGGA
 AGAGAAGGCTAGAAAGCAGATGGAAGGCGAGCGGCTGGCTCGAGAGAAGCAGATGCGGGAGGAGGCCGAG
 CGTACAAGAGATGAGTTAGAGAGGAGGCTCCTGCAGATGAAAGAAGAAGCAACGATGGCCAATGAAGCTC
 TGATGCGCTCTGAGGAGACAGTATCTGTTGGCTGAAAAGGCTCAGATCACAGAGGAGGACCAAGCT
 TTTGGCACAAAAGGCTGCAGAGGCTGAGCAAGAGATGCAGCGAATCAAGGCCACGGCCATTCGGACAGAG
 GAGGAGAAGCGCCTGATGGAGCAGAAGGTGCTGGAGGCTGAAGTCTGGCATTGAAGATGGCTGAGGAGT
 CAGAGAGGAGGCCAAGGAGGCTGATCAGTTAAAGCAAGACTTGAAGAAGCCAGAGAAGCAGAGCGAAG
 AGCCAAGCAGAAGCTCTTAGAAATCGCCACCAAGCCACCTATCCACCCATGAACCAATTCACCACCA
 CTGCTCTGACATACCGAGCTTCGACATTATTGCTGACAGCTTGTCTTTGACTTCAAGGATACGGACA
 TGAAGCGACTTTCCATGGAGATAGAGAAAGAAAAGTGGAGTACATGGAGAAGAGCAAGCACCTGCAGGA
 GCAGCTCAACGAGCTCAAGACGGAGATCGAGGCCCTTAAAAGTCAAAGAGCGGGAGACGGCCTTGGACGTC
 CTACACAGCGAGAGCTCAGACAGAGGCGGCCCCAGCAGCAAGCATAATACCATTAAGGTACCTGAAA
 TG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230590 representing NM_001252252
 Red=Cloning site Green=Tags(s)


MAGAIASRMSFSSLKRKQPKFTFVRIVTMDAEMEFNCEMKWKGDLDLVCRTLGLRETWFFGLQYTIKD
 TVAWLKMDKKVLDHDVSKEEPVTFHFLAKFYPENAEELVQEITQHLFFLQVKKQILDKQVYCPPEASVL
 LASYAVQAKYGDYDPSVHKGFLAQEELLPKRVINLYQMTPEMWEERITAWYAEHRGRARDEAEMEYDKI
 AQDLEMYGVNYFTIRNKKGTELLLGVDALGLHIYDPENRLTPKISFPWNEIRNISYSDKEFTIKPLDKKI
 DVFKFNSSKLRVNLILQLCIGNHDLFMRRRKADSLEVQMQMAQAREEKARKQMERQRLAREKQMRREEAE
 RTRDELERRLLQMKEEATMANEALMRSEETADLLAEKAQITEEEAKLLAQKAAEAQEMQRIRKATAIRTE
 EEKRLMEQKVLAEVLAKMAEESERRAKEADQLKQDLQEAAREARRAKQKLEIATKPTYPPMNPPIPPP
 LPPDIPSFDIADSLSFDFKDTDMKRLSMEIEKEKVEYMEKSKHLQEQNLKTEIEALKLKERETALDV
 LHSESSDRGGPSSKHNTIKKVP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001252252

ORF Size: 1752 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001252252.1 , NP_001239181.1
RefSeq Size:	4748 bp
RefSeq ORF:	1755 bp
Locus ID:	18016
UniProt ID:	P46662
Cytogenetics:	11 3.0 cM
MW:	68.9 kDa
Gene Summary:	Probable regulator of the Hippo/SWH (Sav/Wts/Hpo) signaling pathway, a signaling pathway that plays a pivotal role in tumor suppression by restricting proliferation and promoting apoptosis. Along with WWC1 can synergistically induce the phosphorylation of LATS1 and LATS2 and can probably function in the regulation of the Hippo/SWH (Sav/Wts/Hpo) signaling pathway. May act as a membrane stabilizing protein. May inhibit PI3 kinase by binding to AGAP2 and impairing its stimulating activity. Suppresses cell proliferation and tumorigenesis by inhibiting the CUL4A-RBX1-DDB1-VprBP/DCAF1 E3 ubiquitin-protein ligase complex (By similarity). Plays a role in lens development and is required for complete fiber cell terminal differentiation, maintenance of cell polarity and separation of the lens vesicle from the corneal epithelium.[UniProtKB/Swiss-Prot Function]