

Product datasheet for **MR230658**

Inpp5e (NM_001290437) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Inpp5e (NM_001290437) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Inpp5e
Synonyms:	72kDa; 1200002L24Rik; mKIAA0123
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230658 representing NM_001290437
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCCATCCAAGTCAGCTTGCCCTGCGTCACACTGAGGCACCTGGGCAGCTAGAAGGCAGAATGCTCCAGG
 GACAGCCTCCGAACACGGAAAAGAAGCTCATCCCAACACCAGGCTTCTTACCAGCCTCTGACTCTCAGGG
 TTCAGAGACAAAACCCATGCCTCCTTTACGATTCCAGCAAAAACAAGCAACCAGAACCCACAACTAAG
 GCAAACTCATCACCCACAGCCACCCATTAGACCCAAGTTGGAGCGAACCTGTCCCTCGATGACAAGG
 GCTGGAGAAGGAGGCGTTCCGGGGTAGTCAGGAGGATCTGACTGTCCAGAATGGGGCCAGTCCCTGCAG
 GGGCTCCTTGCAAGACTCAGTAGCCAGTCCCCTGCCTACAGCCGTCCCCTGCCCTGCCTCAGCACATCC
 TTGCAGGAGATACCTAAGTCCCAGGGCCACAGGCAGTGAGGGAGGGAGCCCATCCTTGTGGAGTGACT
 GTCTTTCTGGAATGATCAGCACCTCCTTGACCTCCTGCACAGAGATGCTGCCTCAGGTGGGCCCCCTC
 CAGGTTGGCAAGTTTGCATGCCTCACACACACCACCAGCTATGGACCTCAGCATAGCCTCCAGCTCCCTG
 AGAACAGCAAACAAGGTTGACCCCGAGCACACGGACTACAAGCTCCGCATGCAGACCAGACTGGTACAGG
 CCCACAGCAACCTGGGCCCTAGCAGACCCCGGAGCCCTTGGCTGGTGACGATCACTCCATTCACTCAGC
 CAGATCTTTCAGCCTTCTGGCCCCATCCGCACCAAGGACATCAGAAGCAGGAGCTATCTGGAGGGAAGT
 CTTCTGGCCAGTGGGGCCCTGCTAGGAGCAGAAGAGCTGGCCAGGTACTTCCAGACCGAAACATGGCTC
 TCTTCTGGTACCTGGAACATGCAGGGCCAGAAGGAGCTCCAGCGAGCCTGGATGAGTTTCTGCTCCC
 CACCGAGGCTGACTACACTCAGGACCTGTATGTCATTGGAATTCAGGAGGGCTGCTCTGACAGCGGGAG
 TGGGAGACACCGCTGCAGGAGACTGGGCCCTCAGTATGTAAGTCTGCTCATCAGCAGACATCCACAGTAA
 CAGTACATGTCCTGTTTATCCGTAGGACCTCATCTGGTCTGCTCAGAGGTCGATCACTCCACAGTAA
 TACACGCATCGTGTCTCAGATCAAGACCAAGGGGCCCTGGGCGTCAGCTTACCTTTTTTGGCACCTCC
 TTCTTCTCATCATCTCACTTCACTCTGGAGATGGGAAGGTAGCAGAGCGGCTACTGGACTACAGCA
 GAACCATCAAGCCCTAGCCCTGCCCGGAATGTGCCAGACAAAATCCCTACCGCTCTAGTGCAGGGGA
 TGTCACTACCCGGTTTGTAGAGTCTTCTGGTTTGGGACTTCAACTCCGCCTGAGTGGTGACGAGTG
 GCTGTGGAGGCCCTTCTGAAGCAGAACTGAGGTGGATGTGCTGGCTCTCTCCAACACGACCAGCTCA
 CCCGGGAGATGAAGAAAGGTCATCTTCAAGGGCTTTGAGGAGGAGAGATTCACTTTCTCCATCCTA
 CAAGTTTGACATTGGGAAGGACACTACGACAGCACCTCAAGCAAAGGACACCTCCTACACAGACCGA
 GTCCTATACAAAAGCCGTACAAGGGTACATCTGTCCATGAAGTATTCCTTGTCTGGGATCAAGA
 CTTCAGACCACCTCCTGTGTATGGCTTGTCCAGGTGAAAGTGAGGCCAGGACGAGACAAGAAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230658 representing NM_001290437
 Red=Cloning site Green=Tags(s)


MPSKSACLRHTAPGQLEGRMLQGQPPNTEKKLIPTPGFLPASDSQGSSEINPMPPIFIPAKTISNQNPKTK
 ANLIPTQPPIRPKLERTLSLDDKWRRRFRGSQEDLTVQNGASPCRGSLQDSVAQSPAYSRLPCLSTS
 LQEIPKSRRATGSEGGPSLWSDCLSGMISTSLDLLHRDAASGGPPSRLASLHASHTPPAMDLSIASSL
 RTANKVDPEHTDYKLRMQTRLVRAHSNLGSPRPSPLAGDDHSIHSARFSLLAPIRTKDIRSRSYLEGS
 LLASGALLGAEELARYFPDRNMFVATWNMQGQKELPASLDEFLLPTEADYTQDLYVIGIQEGCSDRRE
 WETRLQETLGPQYVLLSSAAHGVLMSLFIIRDLIWFCSEVEYSTVTTRIVSQIKTKGALGVSFTFFGTS
 FLFITSHFTSGDGKVAERLLDYSRTIQALALPRNVPDTPYRSSAGDVTRFDEVFVWFDFNFRLSGGRV
 AVEAFLKQKPEVDVLLALQHDQLTREMKGSIIFRGFEEAEIHFLPSYKFDIGKDYDSTSKQRTPSYTRD
 VLYKSRHKGDICPMKYSSCPGIKTSDHRPVYGLFQVKVRPGRDKN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001290437

ORF Size: 1815 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001290437.1 , NP_001277366.1
RefSeq Size:	3855 bp
RefSeq ORF:	1818 bp
Locus ID:	64436
UniProt ID:	Q9JII1
Cytogenetics:	2 A3
MW:	67.7 kDa
Gene Summary:	Phosphatidylinositol (PtdIns) phosphatase that specifically hydrolyzes the 5-phosphate of phosphatidylinositol-3,4,5-trisphosphate (PtdIns(3,4,5)P3), phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5)P2) and phosphatidylinositol 3,5-bisphosphate (PtdIns(3,5)P2). Specific for lipid substrates, inactive towards water soluble inositol phosphates. Specific for lipid substrates, inactive towards water soluble inositol phosphates (By similarity) (PubMed:10806194). Plays an essential role in the primary cilium by controlling ciliary growth and phosphoinositide 3-kinase (PI3K) signaling and stability (PubMed:19668215). [UniProtKB/Swiss-Prot Function]