

Product datasheet for **MR230660**

Cacnb2 (NM_001252533) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cacnb2 (NM_001252533) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cacnb2
Synonyms:	AW060387; CAB2; Cavbeta2; Cchb2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230660 representing NM_001252533
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGCTTGACAGGCAGTTGGTGTCTTCTCAGACTCAATCCAGTATTCCTGGGGTTTCGGCAGACTCCTACA
 CCAGCCGCCATCCGATTCAGATGTGTCTTTGGAAGAGGACCGGAGGCAGTACGTAGAGAAGCTGAGCG
 GCAGGCACAGGCACAGTTGAAAAAGCAAAGACAAAACCTGTTGCATTTGCAGTTTCGACCAATGTCAGA
 TACAGCGCAGCCAGGAGGATGACGTTCCGGTGCCTGGCATGGCCATCTCCTTCGAGGCCAAAGATTTTC
 TGCATGTTAAAGAAAAATTAATAATGACTGGTGGATAGGACGGCTGGTTAAAGAAGGCTGTGAAATCGG
 ATTTATTCGAAGCCAGTCAAACCTAGAGAACATGAGGCTACAGCACGAACAGAGACCAAGCAAGGGAAA
 TTCTATCCAGTAAATCAGGAGGAATTCATCATCCAGTTTGGGTGACATAGTCCCAGTTCAGAAAAAT
 CGACACCGCGTCATCTGCTATAGACATAGATGCTACTGGCTTAGATGCAGAAGAGAATGATATTCAGC
 AAACCACCGCTCCCCTAAGCCCAGTCAAACAGTGAACATCACCCCACTCCAAAGAGAAAAGAATGCC
 TTCTTTAAGAAGACAGAGCACACTCCTCCTTATGATGTGGTACCTTCCATGCGACCAAGTGGTGTGGTGG
 GCCCGTCTCTGAAGGGTATGAGGTCACAGATATGATGCAAAAAGCACTGTTTGATTTTCTAAAACACAG
 ATTTGAAGGGCGGATATCCATCACAAGAGTCACTGCTGACATCTCCCTTGCCAAACGCTCAGTATTAAC
 AATCCCAGTAAGCATGCAATCATAGAAAGATCCAAACAAAGTCCAGCTTAGCGGAAGTTTCAGAGTGAAA
 TTGAAAGGATTTTGAACCTGCAAGAACACTGCAATTGGTAGTCCTTGACGCGGATACAATTAATCACCC
 AGCTCAACTCAGTAAGACCTCTTGGCCCTATCATAGTGTACGTAAGATTTCTTCTCCAAGGTTTTA
 CAAAGGTTGATAAAATCTCGAGGGAAATCTCAAGCAAAACACCTCAATGTCCAGATGGTAGCAGTGATA
 AACTGGCCCAATGCTCCTCCGAGGAATCATTGATGTGATCTTGGATGAGAACCAGCTGGAGGACGTTG
 TGAGCATCTTGGCGACTATCTGGAGGCATCTGGAAGGCCACCCACCCTCCAGCGGTAACCTCCCAAC
 CCTCTGCTTAGCCGACTTTAGCCTCCTCAACTTTACCTCTGAGCCCCACCCTCGCTCTAATTCACAGG
 GTTCTCAAGGTGATCAAAGGCCTGATCGCTCTGCTCCCCGGTCTGCTTCCCAAGCTGAAGAAGAACCTTG
 CCTAGAACCCGTCAAAAATCCCAACACCGTCTCCTCAGCCACACACCAAAACCACCGCAGCGGGACA
 GGTGAGGCTCTCTAGGCAAGAGACGTTTACTCTGAAACCAAGAGAGCCGAGACTCTGCCTACGTGG
 AGCCAAAGGAAGATTATTCATGAACATGTGGACCGGTATGTCCACACCGCGAGCATAACCACAGAGA
 GGAGACCCACAGCAATGGCCACAGACACAGGGAGTCTCGCCACCGCTCTAGGACATGGGTGAGAC
 CAGGACCAATGAGTGCATCAAACAACGAAGCCGGCATAAATCTAAGGATCGCTACTGTGACAAGGAGG
 GGAAGTAAATCAAGAGAAGGAATGAGGCTGGCGAGTGAACAGGGATGTATACATCCGCCAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:

>MR230660 representing NM_001252533
 Red=Cloning site Green=Tags(s)


MLDRQLVSSQTQSSIPGGSADSYSRPSDSVLSLEEDREAVRREAERQAQAQLEKAKTKPVAFVRTNVR
 YSAAQEDDVPVPGMAISFEAKDFLHVKEKFNNDWWIGRLVKEGCEIGFIPSPVKLENMRLQHEQRAKQK
 FYSSKSGNSSSLGDIVPSSRKSTPPSSAIDIDATGLDAEENDIPANHRSPKPSANSVTSPPHSKEKRM
 FFKKTEHTPPYDVVPSMRPVVLVGPVSLKGYEVTMMQKALFDLKHREGRISITRVADISLAKRSVLN
 NPSKHAIERSNTRSSLAEVQSEIERIFELARTLQLVVLADTINHPAQLSKTSLAPIIVYVKISSPKVL
 QRLIKSRGKSQAKHLNVQMVAADKLAQCPPQESFDVILDENQLEDACEHLADYLEAYWKATHPPSGNLPN
 PLLSRTLASSTLPLSPTLASNSQGSQDQRPDRSAPRSASQAEPEPCLEPVKKSQHRSSATHQNHRS
 GRGLSRQETFDSETQESRDSAYVEPKEDYSHEHVDVRYVPHREHNHREETHSSNGHRHRESRHRSDMGRD
 QDHNECIKQSRHKS KDRYCDKEGEVISKRRNEAGEWNRDYYIRQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001252533

ORF Size: 1815 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001252533.1 , NP_001239462.1
RefSeq Size:	3814 bp
RefSeq ORF:	1818 bp
Locus ID:	12296
Cytogenetics:	2 10.7 cM
MW:	68.5 kDa
Gene Summary:	The beta subunit of voltage-dependent calcium channels contributes to the function of the calcium channel by increasing peak calcium current, shifting the voltage dependencies of activation and inactivation, modulating G protein inhibition and controlling the alpha-1 subunit membrane targeting.[UniProtKB/Swiss-Prot Function]