

Product datasheet for **MR230669**

Tubgcp2 (NM_001286011) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tubgcp2 (NM_001286011) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tubgcp2
Synonyms:	1700022B05Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230669 representing NM_001286011
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAGAACCCTAGTAAAGGAATACCTGATTCTGGTCACACAGCTGGAACAGCTGCACAGGCAGGGCCTCC
 TCTCACTGCAGAAGCTCTGGTTCTACATCCAGCCAGCCATGCGCACCATCGACATCCTTGCCTCCCTTGC
 CACCTCAGTGGATAAAGGCGAGTGTGTGGCGGGTCGACACTGAGCCTTCTTCATGATAGGAGCTTCAAC
 TACACCGGGGACAGCCAGGCGCAGGAGCTGTGCCTCTATCTGACCAAGGCAGCCAGTGCACCTTACTTTCG
 AGATCCTGGAGAAGTGGATCTATAGGGGCATCATTGACCCATACAGTGAAGTTCATGGTGGAGGAACA
 CGAGCTGCGGAAGGAGAAGATCCAGGAGGACTACAACGACAAGTACTGGGACCAGCGCTACACTGTGCTG
 CCGCAGCAGATCCCATCCTTCTGCAGAAGGTAGCAGGCAAGATCCTCAGCACAGGAAAACTACAAATG
 TGGTCAGAGAGTGTGGTCATGATGTACCTGTCCCGTAGCCAAAGAAATCATCTATACCCTTAAAGAGCG
 GGCATATGTGGAGCAGATAGAGAAAGCGTTCAATTATGCCAGCAAAGTGTGCTGGACTTCTCATGGAG
 GAGAAGGAGCTGGTGGCCACCTGAGGTCCATCAAGCGCTACTTCTGATGGACCAGGGGACTTCTTCG
 TGCAATTCATGGACCTGACTGAGGAGGAGCTCAGGAAGCCAGTGGAGGACATCACTCACACGCCTGGA
 GGCCCTGTTGGAGCTGGCCCTGCGAATGAGCACTGCCAACACTGACCCTTCAAAGACGACCTCAAGATT
 GAACTGATGCCTCAGACCTAATCACTCAGCTTTTGGGGTTCTGGCCATCGAGACCAAGCAGGAAAAGG
 CAATGACCCATGCTGACCCACGGAGCTCACACTGAGTGGTCTGGAGGCCTTCTCCTTGGACTACATGGT
 CAAATGGCCCTGTCAGTATTATTAACAGGAAGGCACTACCCGCTACCAGATGCTCTCCGACACATG
 TTCTACTGCAAGCAGTGGAGCGGCAGCTCTGCAGTGTCTGGATCAGCAACAAGCAGGCAAGCAGCATG
 CACTGCACCTCCGCAAGTGGTTTGTGGTCTTCACTCTGCGGCAGGCAATGCTCAACTTGTCCAGAA
 TATTCAGTACTACATGATGTTTGAAGTATGGAACCGACCTGGCACATCCTTGAGAAAACTGAAATCT
 GCCTCCAACATTGACGATGTCCTTGGCCACCACACAAGCTTCTTGACAAGTGCCTGAAGGACTGCATGC
 TGACCAACCTGAGCTGTGAAGTTTTCTCCAAGCTCATGTCCGTATGCGTCATGTTTACCAACTGTAT
 GCAGAAATTTACACAGAGTATGAACTGGACAGTGAAGTGGAGTGGAGTGGAGTGGAGTGGAGTGGAGT
 CAGGGGCCACCCACAGAGACTGAACGGATTGAGGATCGCCCCGGAAGGAGCTTACCAGGAAGCATCTGT
 CTGAGCATGTGGATGCTCCTCAGCTGGCCTCTGGTTTCGAGGCCACCATCAACAACCTTGGACAAAACTT
 CTCTGCCATCTTCTTGATTTGCTGGCCCGACTAAGTGTCTACAGCACCAGTGAAGCATGGCATG
 GCCAGTGTCTCCAGACTAGATTTCAATGGTTTCTATGCTGAACGCCTGGAACGTCTGTCTGCAGAAA
 GGAGTCAGAAGGCTGCCCCCAAGTGCCTGTCCCACGGGGACCCTCAGCACCAGCACCAGGGTGGCCAT
 TCTGCGCAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR230669 representing NM_001286011
 Red=Cloning site Green=Tags(s)

MRTLVEYLILVTQLEQLHRQGLLSLQKLFYIIPAMRTIDILASLATSVDKGEVGGSTLSLLHDSRFN
 YTGDSQAQELCLYLTKAASAPYFEILEKWIYRGIHDPYSEFMVEEHLELRKEIQEDYNDKYWDQRYTVL
 PQQIPSFQKLVAGKILSTGKYLNVVRECGHDVTCVPAKEIIYTLKERAYVEQIEKAFNYASKVLLDFLME
 EKELVAHLRSIKRYFLMDQGDFFVHFMDLTEEELRKPVEDIILTRLEALLELALRMSTANTDPFKDDLKI
 ELMPHDLITQLLRVLAIETKQEKAMTHADPTLTLGLEAFSFDYMKWPLSLIINRKALTRYQMLFRHM
 FYCKHVERQLCSVWISNKTAKQHALHSKWFAGAFTLRQMLNFVQNIQYMMFVEMPTWHILEKNLKS
 ASNIDDLVGHHTSFLDNCLKDCMLTNPELLKVF SKLMSVCMFTNCMQKFTQSMKLDSELGRLTLDQGSV
 QGPPTETERIEDRPRKELTRKHLSEHVDAPQLASGFEATINNFDKNFSAHLLDLLARLSVYSTSDCEHGM
 ASVISRLDFNGFYAERLERLSAERSQKAAPQVPVPRGPSAPAPRVAIPAQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001286011

ORF Size: 1830 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001286011.1 , NP_001272940.1
RefSeq Size:	2925 bp
RefSeq ORF:	1833 bp
Locus ID:	74237
Cytogenetics:	7 F4
MW:	71 kDa
Gene Summary:	Gamma-tubulin complex is necessary for microtubule nucleation at the centrosome. [UniProtKB/Swiss-Prot Function]