

Product datasheet for **MR230705**

Tiam2 (NM_001286757) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Tiam2 (NM_001286757) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tiam2
Synonyms:	3000002F19Rik; mKIAA2016; STEF
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230705 representing NM_001286757
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAAGCACCAACAGAGGCCATGACCCACCTCCAGGCCTCTGGCTCGTCACCTCTCAGATGCAGATC
 GCCTCCGAAAGTCATCCAGGAGCTTGTAGACACAGAGAAGTCTTACGTGAAGGATCTGAGCTGCCTCTT
 TGAATACTTGGAGCCACTTCAGAATGAGACCTTTCTACCCAAGATGAGATGGAGTCACTTTTTGGG
 AGCCTGCCAGAGATGCTGGAGTTTCAAAGGTGTTCTGGAGACGTTGGAGGATGCGATCTCCGCTTCT
 CGGACTTTAGTGTCTGGAAACCCCTCACAGTTTCGAAAATTGCTGTTCTCCCTGGAGGTTCTTTCT
 TACTATGCGGATCACTTAAGCTATACAGTGGGTTCTGTGCCAACACATTAAGTACAGAGGGTTCTA
 GAGCGAGCTAAAACGGACAAGGCCTTCAAGGCTTTCTGGATGCCGAAATCCCACCAAGCAGCACTCT
 CCACGCTGGAGTCTATCTCATCAAGCCTGTTCCAGAGAGTGTCAAGTATCTCTGCTTCTCAAGGAGT
 AGTGTCACTGACTGACCATGAGAGTGAAGAACAATCACCTGACAGAAGCACTAAAGGCCATGGAAAA
 GTCGCCAGTCACTCAATGAGATGCAGAAGATCTACGAGGATTACGGGATGGTGTGGTGGACAGCTGGTGG
 CAGAGCAGAGTGGCACAGAGAAGGAAGTACAGAGCTGTCCATGGGGAACTTCTGATGACTCTACAGT
 TTCTGTTGAATCCGTTCTGTCTTAGGAAAAGCCAGGAAGGACATTGAGCTCACAGTATTTGTTTTT
 AAGAGAGCTGTACTACTGGTTTATAAAGAAAAGTCAAGCTGAAAAGAAAAGTGCCTCGAATCCCGGC
 CTGCTCACAACCTGTGACTGGATCCATTTAAATCCGCTGGTTGATTTCCATATCTGCGCTTCAAGT
 TAGACTGGGGAACACGGCAGGGACTGAAAATAATCCACGTGGGAGCTGATTCATACCAAGTCGAAAT
 GAAGGACGGCCAGAAACCATCTTCAACTGTGCTGCAGTGACAGCGAGAACAACCAAGCATTGTTAAGG
 TGATTCGTTCTATTCTGAGAGAGAACTCCGGCGCCACATAAAGTGTGAGCTGCCACTGGAGAAGACGTG
 TAAGGACCGGCTAGTACCTCTTAAGAACCGAGTTCCTGTTTCAGCCAAATTAGCCTCGTCCAGGTCGTTG
 AAGGGCTCAGAACATCTCCAGCAGCAGTGGCCAGCGAGCCAGCAAGGGCAACTACTGGACTCAG
 ATGAGTGCAGCCTGAGCAGTGGCACCCAGAGTAGCGGCTGCCCGTAGCCGAGAGCAGGCGAGACTCTAA
 GAGCACCAGCTGGAGAAAGACGCTCAGGAGGGCTGGCGGAGTTTCCAGATGGTCTTATCAAAGAAAGC
 GACATTCTGAGTATGAAGATGAGGACTCCACCACCCTCTGAAACAGGGTAGCCCTACTAAGGACATTG
 AGATTCAGTTCCAGAGACTGAAAATCTCTGAGGAATCCGACGTGCACCCAGTTGGCAGCAGCCTCTCAC
 AGAGTCAGGTGAACAGCCCAAGCTGGTCAGGGGCCATTTTGGCCATTAACGGAAAGCAAACAGCACC
 AAGAGGGGCAGAGAACTTTGCTCAAGGCGCAGACTCGTACCAGTCCCTGGACAGCCACCCAGAACTG
 CCAGCATTGATCTAACTTGGTCTGGAGAGAGAATTGAGTCCAGAGCTTAACCTCAGTCGTCATGA
 GGAGGGTTTTATGAAACACAGAGCCATGGCAAATCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR230705 representing NM_001286757
 Red=Cloning site Green=Tags(s)

MEAPTESHDPPRPLARHLSDADRLRKVIQELVDTEKSYVKDLSCLFELYEPLQNETFLTQDEMESLFG
 SLPEMLEFQKVFLETLEDAISASDFSVLETSPQFRKLLFSLGGSFLYYADHFKLYSGFCANHIKVRVL
 ERAKTDKAFKAFLDARNPTKQHSSTLESYLKPVQRVLKYPLLLKELVSLTDHESEEHYHLTEALKAMEK
 VASHINEMQKIYEDYGMVFDQLVAEQSGTEKEVTELSMGELLMHSTVSWLNPFLSLGKARKDIELTVFVF
 KRAVILVYKENCKLKKLPSNSRPAHNSADLDPFKFRWLIPISALQVRLGNTAGTENNSTWELIHTKSEI
 EGRPETIFQLCCSDSENKTSIVKVIKRSILRENFRRIKCELEPKTCKDRLVPLKNRVPVSAKLASSRSL
 KGLRTSSSEWPSEPSKGNLSDDECSLSSGTQSSGCPVAESRRDSKSTELEKDAQEGLAEFPDGLIKES
 DILSDEDEDFHHPLKQGSPTKDIEIQFQRLKISEESDVHPVGGQPLTESGEQPKLVRGHFCPIKRIKANS
 KRGRGTLKKAQTRHQSLDHPETASIDLNLVLEREFVQSLTSVNEEGFYETQSHGKS


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001286757

ORF Size: 1857 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001286757.1 , NP_001273686.1
RefSeq Size:	3186 bp
RefSeq ORF:	1860 bp
Locus ID:	24001
UniProt ID:	Q6ZPF3
Cytogenetics:	17 1.99 cM
MW:	70.6 kDa
Gene Summary:	Modulates the activity of RHO-like proteins and connects extracellular signals to cytoskeletal activities. Acts as a GDP-dissociation stimulator protein that stimulates the GDP-GTP exchange activity of RHO-like GTPases and activates them. Activates specifically RAC1, but not CDC42 and RHOA. Mediates extracellular laminin signals to activate Rac1, contributing to neurite growth. Involved in lamellipodial formation and advancement of the growth cone of embryonic hippocampal neurons. Promotes migration of neurons in the cerebral cortex. When overexpressed, induces membrane ruffling accompanied by the accumulation of actin filaments along the altered plasma membrane.[UniProtKB/Swiss-Prot Function]