

Product datasheet for **MR230707**

5830473C10Rik (NM_001252661) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	5830473C10Rik (NM_001252661) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	5830473C10Rik
Synonyms:	Arg; Gm17754
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230707 representing NM_001252661
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGGAGGTTCTTAAAGCTCTACCCTTCATTGTTTTCTTCAGCTACAGTAAATGTGAGTCTCTCCAAACAA
 GTTTACTACATGCAGCCAGGCAGCATGTCAAGTGGCAGGATCACCTGGAGGAGAAGTTAGGAGTCAGTGT
 CAGCATTATGGTGGCTCAGTTTCTTCAGGATGCAGCTTATGACCAAGTCCAAACCATAGCTAAAGAGCTG
 CTGAAGCTTGCTGAAAAGTGAAGCGTCTCAAACCTGAGTCACCCTCAGAGTGTGCTCATCAACTGGTGG
 CTGCCCTCCTTATACACATTTGCAACAACCAAGGCCTGGTGGATAGTCACGTGTTCACTGACTGTTGTAA
 GATGAAAAGTGCAGCAAGACTCAGGTGCTTCTCTCTACAAGAAAGACGACGACGACAGCCATGATACA
 CCCCTGATTCCAAGGCCTGAAGTACCTGTGAGGTGGCTGCAGAGAATAATGTGTCCTGAAGGAGAGGT
 TCAGCTATGAAATATCTAGAAGACCCATTCTTACGGACCCACTATCTTGAGCATGTCTGCTTGCTA
 TGATACAGCTGTTCAAGTATGTTGTCAAGAAGAGAATGAAACAGAGTGCCTGCGGACAAAGCTAGAACC
 ATCAGAAAATACATTAGAGGAATATCTGCAAGACATCATCACTTATGTGAAATTGGCACCAAATCAACG
 AGAAAGTAGCAAGTGCAGTAGAGTTGATTCTACTACCAAAAAGCAACCTAAAGCTAATTTTTCTGAAAT
 CGCCAAATTGACCATGGATGTTAAAAATCTGCATGAGACCTGTTGTGAGGGGAACACATTAGTGTGTGCC
 CTGGGCAGGCGCCAGCTTATGAATTACACCTGCTTAACCAGGCCATCTTATCCAGCAAGTTTGCTGAGT
 GCTGTGAACAGCCAGAACCATTCCGCGGGGAGTGCATCGTCAACTCAGAAAATGACGACAAGCCTGACCT
 TTCACCCTGCCGCTCACAAGGTTTACAGAAGACTGGTCTGTGTGCCAACATTCAGCGGCAACCAGGAC
 AGCTTTTACAGGAGTTTCTTTATGAATACTCAAGAAGACACCTGGAGTTGGCGGTTCCCGTGATCTTAA
 GAGTCTTTGAAACATACAACCGTATTGGAGAAATGTTGCAAGCTAGAAAAGCTTGGAAATGCCACCA
 GCCTGGGAAAGAGATGTTCCAAAAGTGATACTTGAAGCCGTGACCAAGTTAAAAGTACTGTGACTTAA
 CAGAAGCAGCTGGCAGTAGCAGCTTCCAGGACAGGCTCACGGTCTTTATACTAAGAAAGCTCCGACG
 TGTCTGCCAGGAGCTGGTGGTGTTCGTGAAGAAGATGGCAGCGTTGCCAGGAGATGTTGCCCATTAAG
 GGATGAGCTGCAGTCCGCTTGCCTGGAGGACCAGACAAAGCTGATTCTTGGAGAGCTGTGCAGAAGACAC
 GGTGGGAAGCCTGTCAATGCAGGAGTGGCCACTGCTGCGACGACTCCTATGCCTTCAGAAAAGCCGTGCT
 TTGATGACCTGCGAGTCGATGGCACTTACATCTCTCCACCTTTATCTTGTGACCAAGCCATTAGTCTTAA
 AGAGGACTTGTGCCAAGCTCAGGACGAGGAAGTCAAATGAAAAGCAGAAGCTCCTCAGCACCCCTGGT
 AAGCAGAAGCCCATGCAGCAGAGGATGTGTTCCATTCCATTGGCGAGAATTTTGTGCAACTGGTGGAGA
 AGTGTCTGCGGTGCACAGAGGAGAGACGTGTGTTTTCAAGAAGAGGGACCTCAACTGATCAGGAAGTCCCA
 GTCCTCTTAGAGGCCAGCTCATCCAGAGCGGGTTGATG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR230707 representing NM_001252661
 Red=Cloning site Green=Tags(s)

MEVLKALPFIVFFSYSKCESLQTSLLHAARQHVSGQDHLLEENLGVSVSIMVAQFLQDAAYDQVQTI AKEL
 LKLAECKRKLKPEPSECAHQLVAAFLIHCNNQGLVDSHVFTDCKMKTAARLR CFLSYKKDDADSHDT
 PLIPRPELTCEVAENNVSLKERFSYEISRHPFLYGPTILSMSACYDTAVQSCCQEENETECLRTKLEP
 IRKYIRGISARHHHLCEIGTKFNEKVASAVELILLTKKQPKANFSEIAKL TMDVKNLHETCCEGNTLVCA
 LGRRLMNYTCSNQAILSSKFAECCEQPEPFRGECIVNSEDDKPDLSPLPLTRFTEDWSVCQQFSGNQD
 SFLQEFLEYYSRRHLELAVPVILRVFETYKRVLEKCKLENSLECHQPGKEMFQKVILESRDQVKTYCDL
 QKQLGSSSFQDRLTVLYTKKAPQLSAQELVVFVKMAAVARRCCPLRDELQSACVEDQTKLILGELCRRH
 GKPVNAGVGHCCDDSYAFRKPCFDDL RVDGTYISPLSCDQAI SLKEDLCQAQDEELQMEKQKLLSTLV
 KQKPHAAEDVFHSIGENFVQLVEKCCGAQR RDVCFQEEGPQLITKQSLLEASSSQSGLM


TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001252661

ORF Size: 1860 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001252661.1 , NP_001239590.1
RefSeq Size:	2241 bp
RefSeq ORF:	1863 bp
Locus ID:	622307
Cytogenetics:	5 E1
MW:	70.6 kDa