

Product datasheet for **MR230729**

Cadps2 (NM_001252108) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cadps2 (NM_001252108) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cadps2
Synonyms:	A230044C21Rik; Caps-2; Caps2; Cpd2
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230729 representing NM_001252108
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGCTCGACCCGCTTCCAGCGAGGAGGAGTCGGACGAGGGGCTGGAGGAGGAGACCGCGAGGTGTTGG
 TGGCGCCGGCGTCTCCAGCGAGCGCCCGCCGGCTCGGGAAGGGCGCGGGACGCGCCGGGACG
 CTCGGTGGTGGCAGCGCGCGGAGCGCCAGACCCGTGAGCCGAGCCCTCGGTACTCAGCGAGGGG
 CGAAATGAGCCGAGCTGCAGCTGGACGAGGAGCAGGAGCGCGCATCCGCCTGCAGCTCTACGTCTTCG
 TGGTGAGGTGCATCGCGTACCCATTCAACGCCAAGCAGCCACCACATGGCCCGGAGGCAGCAGAAGCT
 TAACAAACAACAATTGCAGTTACTGAAAGAACGGTTCAGGCCTTCTCAATGGAGAACTCAGATAGTA
 GCAGACGAAGCCTTCTGCAATGCCGTCCGGAGTACTATGAGGTGTTCTAAAGAGTGACCGTGTGGCCA
 GAATGGTCAAAGTGAGGATGCTCTGCTAATGATTTTCAGAGAGGTGTTTAAAGAACAATAGAAAAACG
 TGTCGGTAGCCTGCCAGAGATAGATGGCCTGAGCAAAGAGACAGTGTGAGCTCATGGATAGCCAAATAT
 GACGCCATTTACAGGGGAGAAGAAGATTTGTGCAAACAGCCGAACAGGATGACCCCTAGTGCCGTGTCTG
 AGCTTATTCTGAGTAAGGAGCAGCTCTATGAAATGTTCCAGCAGATTCTGGGCATTAAGAAGCTGGAACA
 CCAGCTACTTTATAACGCATGCCAGCTGGATAACGCTGATGAACAAGCAGCCAGATCAGAAGGGAACCT
 GATGGCCGCTGCAGTTGGCAGAGAAGATGGCAAAGGAGAGAAGATCCCCAGATTCATCTCGAAAGAAA
 TGGAGAGTATGTACATTGAAGAGCTGCGAGCCTCCGTGAACCTGCTAATGGCCAATTTAGAAAGCTTCC
 AGTTTCAAAGGTGGCCCGAATTTAAATTACAAAAATTAAGCGTTCGAGAAGCTGCGTTTCTGGAC
 CTGGGAGATGAGAACGAGATCCAGCTGTCCAAGTCGGATGTGGTGTGTCGTTACGTTAGAGATTGTCA
 TCATGGAAGTCAAGGACTGAAATCTGTGGCTCCAATCGAATCGTTTACTGCACAATGGAGGTGGAAGG
 AGGAGAAAAACTCCAGACAGACCAGGCTGAAGCATCAAGGCCACAATGGGGGACCAAGGAGATTTCAAC
 ACTACCCACCCTCGCCTGTCGTCAAAGTGAAGCTCTTACAGAAAGCACGGGGTCTGCCCCTGGAAG
 ACAAGGAACTGGGCAGGGTGGTGTATACCCAACTTCTAATAGCTCCAAGTCAGCAGAGTTACACCGAAT
 GACAGTACCCAAGAACAGTCAGGACTCGGACCTAAAGTCAAATGGCAGTGCATGGATAAACCCAGCA
 CACATGAAGCATAGTGGTACCTGTATGCCCTGGACAGAAGGTTTGGAAAAGATGGAAAAGCGTTACT
 TTGTTCTCGTTCAGTTCAGCAGTACACCTTTGCTATGTGCAGCTATAGAGAAAAAAGTCGGAACCACA
 GGAATTAATGCAACTGGAAGGATACACAGTGGATTACACAGACCCTCACCCAGGCCTTCAGGGTGGTCAG
 GTGTTCTCAACGCTGTTAAAGAGGGAGATACTGTGATCTTTCAGTGTGATGAACAGGACAGAATAT
 TATGGGTACAAGCCATGTACAGGCTACAGGCCAGTCTTACAAACAGTTCCTGCAGTCAAAGCCAGAA
 GCTGAATCCTAAAGCGGAGCTCTCCATGCAGATGCTCAGCTTTGTAAGTTTTGTAAAAAT

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence:


>MR230729 representing NM_001252108
 Red=Cloning site Green=Tags(s)

MLDPSSSEEESEGLEEEESREVLVAPGVSRAPPAAREGRRDAPGRSGGSGGGAARPVSPSPSVLSEG
 RNEPELQLDEEQERRIRLQLYVFVVRCAIYPFNAKQPTDMARRQQLNKQQLQLLKERFQAFNGETQIV
 ADEAFNCNVRSYEYVFLKSDRVARMVQSGGCSANDFREVFKNIEKRVRSLEPIDGLSKETVLSWIAKY
 DAIYRGEEDLCKQPNRMTLSAVSELILSKEQLYEMFQQILGIKKLEHQLLYNACQLDNADEQAAQIRREL
 DGRLQLAEKMAKERRFPRFISKEMESMYIEELRASVNLNMANLESLPVSKGGPEFKLQKLKRSQNSAFLD
 LGDENEIQLSKSDVVLSTLEIVIMEVQGLKSVAPNRIVYCTMEVEGGEKLQTDQAEASRPQWGTQGFN
 TTHPRPVVVKLFTTESTGVLALDEKELGRVVL YPTSNSSKSAELHRMTVPKNSQSDLKIKLAVRMDKPA
 HMKHSGYLYALGQKVWKRWKRYFVLVQVSQYTFAMCSYREKKSEPELMLQLEGYTVDYDTPHPGLQGGQ
 VFFNAVKEGDTVIFASDDEQDRILWVQAMYRATGQSYKPVPAVQSQKLNPKGGALHADAQLCKFCEN


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites: SgfI-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001252108

ORF Size: 1881 bp

OTI Disclaimer:	The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info
OTI Annotation:	This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001252108.1 , NP_001239037.1
RefSeq Size:	2232 bp
RefSeq ORF:	1884 bp
Locus ID:	320405
Cytogenetics:	6 A3.1
MW:	71.4 kDa
Gene Summary:	Calcium-binding protein involved in exocytosis of vesicles filled with neurotransmitters and neuropeptides. Probably acts upstream of fusion in the biogenesis or maintenance of mature secretory vesicles. Regulates neurotrophin release from granule cells leading to regulate cell differentiation and survival during cerebellar development. May specifically mediate the Ca(2+)-dependent exocytosis of large dense-core vesicles (DCVs) and other dense-core vesicles.[UniProtKB/Swiss-Prot Function]