

Product datasheet for **MR230751**

Zfp976 (NM_001242388) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zfp976 (NM_001242388) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Zfp976
Synonyms:	9830147E19Rik; C330019L16Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide
Sequence:

>MR230751 representing NM_001242388
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGCATCGCC**

ATGGATGCATTGACCTTTGATGATGTGCATATACACTTCACTCGGGAGGAGTGGTCTTGCTGGATCCAT
CTCAGAAGAGACTTTACCAAGATGTCATGCTGGAGAACTACAGGAATCTCACCACATATAGGGTACAATTG
GAAAGATCATAACATTGAAAGAACATTCTCAAAATGATAGAAGATATGGAAGGCACGAAAGAAGTCTTCT
GCAGAGAAACCCTATGAATATACACAATGGAATAAAGCCTTCAAAATACACAGTCATCTTCAAAAGCTTG
AAAGAAACCAAAGTAGAGACAGACACTCTGAAGTCTTTCAACATAGTAAAGTGCATGCATGTCAAAGTAC
TCTCCAAATAAATAAAGAACAATACTGCACAGAAACCCTACAATGTAATGAATGTGATAAAGCCTTC
TTAAAGCTCAAGTATCTCGAGCACATGAAAGAACACATTCTGGGAAGAAACCCTATGAATGTAGTCAAT
GTGGTAAAGCCTTTTTGCATCCCTGCTATCTTCGAATTCATGAAAGAACACATACTGGAGAGAAACCTTA
CCAATGTAATCAATGTGCTAAAGCCTTTACTCTAAGGAGTCAGCTTCAAAGCCATAAAGTAGTACATACT
GGAGTGAGACCCTACATATGTAAGCAGTGAATAAAGCCTATCCACAGCAGAGTTATCTCCGAAAACATG
AAAGAACACATACTGGAGGGAAACCTTATATATGTAATCTGTGTGATAAAGCCTATTCACACCAGAGTAA
TCTCCGAATACATGAAATAACACATACTGGAGAGAAACCTTATATATGTAATCTGTGTGATAAAGCCTAT
TCACACCAGAGTAATCTCCGAATACATGAAAGAACACATACTGGAGAGAAACCTTATATATGTAATCTGT
GTGATAAAGCCTATTCACACCAGAAATATCTCCAAATACATGAAAAACACATACTGGTGAGAAACCTTA
TAAATGTAATCAGTGTGATAAAGCCTTTGCACAGCTCAGGTACCTTCGAGCACATGAAAGAACACATACT
CAAGAGAAAACCTACAAATGTAATCAATGTGGTAAAGCCTTTGCACAGCTCACCTATCTTCGAGTACATG
AAAGAAGACATAGTGGAGAGAAACCTACAATGTAATGAATGTGATAAAGCCTTTTTACAGCTCAAGTA
TCTTCGATTACATGAAAGAAGACATAGTGGAGAGAAACCTATGAATGTAATCAATGTGGTAAAGCCTTT
TTACATCCAGCTATCTTCGAATGCATGAAAGAACACATACTGGAGAGAAACCTTACCAATGTAATGAAT
GTGATAAAGCCTTCTTACATCCCAGCTATCTTAGAATGCATGAAAGAAGACATAGTGGAGAGAAACCTTA
TGAATGTAATCAATGTGGTAAAGCCTTTTTGCAATCCCCTATCTTCGAATGCATGAAAGAACACATACT
GGAGAGAAACCTATGAATGTAATCAATGTGGTAAAGCCTTTTTACATCTCAGCTATCTTCGAATGCATG
AAAGAAGACATAGTGGAGAGAAACCTATGAATGTAATCAATGTGGTAAAGCCTTTTTATATCCCAGCTA
TCTTCGAATGCATGAAAGAACACATAGTGGAGAGAAACCTATGAATGTAATCAATGTGGTAAAGCCTTT
TTACATCCAGCTATCTTCAAATGCATGAAAGAACACATAGTGGAGAGAAACCTATGAATGTAATCAAT
GTGGTAAAGCCTTTTACATTCCAGCTATCTTCAAATGCATGAAAGAACACATACTGGAGAGAAACCTTA
TGAATGTAATCAATGTGGTAAAGCCTTTTTGCAACCCACTATCTTCGAAAACATGAAAAACACATACT
AGAGAGAAACCTTACAAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230751 representing NM_001242388
Red=Cloning site Green=Tags(s)

```
MDALTFDDVHIHFTREEWSLLDPSQKRL YQDVMLENYRNLTTIGYNWKDHNIEEHSQNDRRYGRHERSPS
AEKPYEYQWNKAFKYHSHLQKLERNQTRDRHSEVFQHSKVHACQSTLQINKRNTNTAQKPYKCNECDKAF
LKLKYLRAHERTHSGKKPYECSQCGKAF LHPCYLRIHERTHTGEKPYQCNQCAKAF TLRSQLQSHKVVHT
GVRPYICKQCNKAYPQQSYLRKHERTHTGGKPYICNLCDKAYSHQSNLRIHEITHTGEKPYICNLCDKAY
SHQSNLRIHERTHTGEKPYICNLCDKAYSHQKYLQIHEKTHTGEKPYKCNQCDKGF AQLRYLRAHERTHT
QEKT YKCNQCGKAF AQLTYLRVHERRHSGEKPYKCNECDKAF LQLKYLRLHERRHSGEKPYECNQCGKAF
LHPSYLRMHERTHTGEKPYQCNCDKAF LHPSYLRMHERRHSGEKPYECNQCGKAF LQSHYLRMHERTHT
GEKPYECNQCGKAF LHSYLRMHERRHSGEKPYECNQCGKAF LPSYLRMHERTHSGEKPYECNQCGKAF
LHSSYLQMHERTHSGEKPYECNQCGKAF SHSSYLQMHERTHTGEKPYECNQCGKAF LQPHYLKHEKHTHT
REKPYK
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001242388

ORF Size: 1908 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001242388.1](#), [NP_001229317.1](#)

RefSeq Size: 5011 bp


RefSeq ORF: 1911 bp

Locus ID: 208111

Cytogenetics: 7 B3

MW: 76.5 kDa

Product images:



Circular map for MR230751