

Product datasheet for **MR230755**

Pomgnt1 (NM_001290658) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Pomgnt1 (NM_001290658) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pomgnt1
Synonyms:	0610016I07Rik; 4930467B06Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)



[View online »](#)

ORF Nucleotide Sequence:

>MR230755 representing NM_001290658
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGCATCGCC**

ATGGAGGGGACCAAGGAAGTCAAGGACTCAAATGGGAAAATACAAGACCACGGGACAGGGCTGTGCTTT
 TCCTGCTGGTGACCGTCATTGTCAACATCAAGTTGATCCTGGACTCGGAGAGCAATCAGCGAGGCCAA
 TGAAGACCCAGAACAGCAAGACTATGATGAGGCCCTAGGACGTCTGGAGTCCCTCGGCGCAGAGGC
 AGCAGTCCCCGACGGGTCTGGATGTGGAAGTATACTCCAGTAGAAGCAAAGTGTATGTGGCTGTGGATG
 GCACCACGGTGTCTGAGGATGAAGCCCGGAGCAGGGCCGAGGCATCCATGTCATAGTCTCAACCAGGC
 CACGGGCCACGTGATGGCAAAGCGCGTGTGGACATACTCTCCTCATGAGGATGAGGCCATGGTGTCTG
 TTCCTCAACATGGTGGCCCCGGCCGCGTGTCTGACCGGTGAAGGATGAGGGTTCTTCCACCTCA
 AGGACACAGCAAAGGCTCTGCTGAGGAGCTGGGAAGTCAGGCAGGCCCTGCCCTTGGCTGGAGGGACAC
 CTGGGCCCTTCGTGGGACAAAAGGAGTCTGTACTTGGCGAGAAACATTCTAAATCTCCTGCCCTCTCC
 TCCTGGGGTGACCCAGTCTGTGAAGACAGATGTGCCACTGAGCTCAGCTGAAGAGGCTGAGTGTCACT
 GGGCAGATACGGAGTTGAACCGCCGCGCAGGCCGCTTCTGCAGCAAAGTGGAGGGCTACGGGAGCGTGTG
 CAGCTGCAAGGACCCGACGCCTATTGAGTTCAGCCCTGACCCACTTCCAGACAATAAGGTCTCAATGTG
 CCTGTGGCTGTATTGCAGGGAACCGGCCCAATTACCTGTACAGGATGCTGCGCTCTCTGCTCTCAGCCC
 AAGGGGTATCTCCACAGATGATAACGGTCTTATTGATGGCTACTATGAGGAGCCAATGGATGTGGTGGC
 GCTGTTTGGTCTGAGAGGCATCCAGCATACTCCTATCAGCATCAAGAACGCCCGTGTGTCTCAGCACTAC
 AAGGCCAGCCTCACTGCCACTTTCAACCTGTTCCGGAGGCCAAGTTTGTGTGGTTCTGGAAGAGGACC
 TGGACATTGCCGTGGATTTTTTTCAGTTTTCTGAGCCAGTCCATCCACCTGCTGGAGGAAGATGACAGCCT
 GTACTGCATCTCTGCCTGGAATGACCAGGGGTATGAACACACAGCTGAGGATCCAGCACTACTGTACCCT
 GTGGAGACCATGCCTGGGCTAGGATGGGTGCTCAGGAAATCTTTATAAAGGAGGAGCTTGAACCAAGT
 GGCCACACCCAGAAAAGCTGTGGGATTGGGACATGTGGATGCGGATGCCGGAACAGCGCCGGGGCCGAGA
 GTGTATCATCCCGGATGTTTCCCGATCCTACCCTTTGGCATCGTCGGCCTCAACATGAACGGCTACTTT
 CATGAGGCCTACTTCAAGAAGCACAAGTTCAACACAGTTCAGGTGTCCAGCTCAGGAATGTAGACAGTC
 TGAAGAAAGAAGCTTATGAAGTGGAAATTCACAGACTGCTCAGTGAAGCCGAGGTTTTGGACCACAGCAA
 GGACCCTTGTGAAGACTCTTCTGCCAGACACGGAGGGCCACACCTATGTGGCGTTTATCCGAATGGAG
 ACAGATGATGACTTCGCCACCTGGACTCAGCTTGCCAAGTGCCTCCATATTTGGGACCTGGATGTTCTGTG
 GCAACCACCGGGCCTGTGGAGATTATTTCCGGAAGAAGAACCACTTTTTGGTGGTGGGGTCCCAGCCTC
 CCCTTACTCAGTGAAGAAGCCACCATCAGTCACCCCTATTTTCTGGAGCCACCCCAAAGGAAGAGGA
 GCCCAGGGGCTGCTGAACAAACG

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230755 representing NM_001290658
Red=Cloning site Green=Tags(s)

```
MEGTKEVKDSNGKIQDHGTGAVLFLLVTVIVNIKLILDTRRAISEANEDPEPEQDYDEALGRLESPRRRG
SSPRRVL DVEVYSSRSKVVYAVDGTTVLEDEAREQGRGIHVIVLNQATGHVMAKRVFDTYSPHEDEAMVL
FLNMVAPGRVLICTVKDEGSFHLKDTAKALLRSLGSAQAGPALGWRDTWAFVGRKGGPVLGEKHSKSPALS
SWGDPVLLKTDVPLSSAEEAECHWADTELNRRRRRFCSKVEGYGVCCKDPTPIEFSPDPLPDNKVLNV
PVAVIAGNRPNYLYRMLRSLLSAQVSPQMIVTFIDGYEPEMDVVALFGLRGIQHTPIISIKNARVSQHY
KASLTATFNLFPEAKFAVVLEEDLDIAVDFFSFLSQSIHLLLEEDSLYICISAWNDQGYEHTAEDPALLYR
VETMPGLGWVLRKSLYKEELEPKWPTPEKLWDWMWMMRPEQRRGRECIIPDVSRSYHFIVGLNMNGYF
HEAYFKKHKFNTVPGVQLRNVDSLKKEAYEVEIHRLLEAEVLDHSDPCEDSFLPDTEGHTYVAFIRME
TDDDFATWTQLAKCLHIWDL DVRGNHRGLWRLFRKKNHFLVVGVPASPYSVKKPPSVTPIFLEPPPKEEG
APGAAEQT
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



ACCN: NM_001290658

ORF Size: 1914 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001290658.1](#), [NP_001277587.1](#)

RefSeq Size: 2650 bp

RefSeq ORF: 1917 bp

Locus ID: 68273


UniProt ID: [Q91X88](#)

Cytogenetics: 4 D1

MW: 72.5 kDa

Gene Summary: Participates in O-mannosyl glycosylation by catalyzing the addition of N-acetylglucosamine to O-linked mannose on glycoproteins. Catalyzes the synthesis of the GlcNAc(beta1-2)Man(alpha1-)O-Ser/Thr moiety on alpha-dystroglycan and other O-mannosylated proteins, providing the necessary basis for the addition of further carbohydrate moieties. Is specific for alpha linked terminal mannose.[UniProtKB/Swiss-Prot Function]

Product images:



Circular map for MR230755