

Product datasheet for **MR230796**

Nfatc2 (NM_001291172) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfatc2 (NM_001291172) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfatc2
Synonyms:	AI607462; NF-ATc2; NF-ATp; NFAT1; NFAT1-D; Nfatp
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>MR230796 representing NM_001291172
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGCAGCGCCCCGATTGCGCAGGCCGGGCCCCGCCAGCCCCTCCGACCATGGGGTCTCGGACCGAG
AAGAACCATCGCACATAAGGCCATCAGCTCACCTCCGGACTCGCATACCCGGATGATGTCCTGGACTA
TGGCCTCAAGCCATGCAACCCCTTCCAGTCTCTGCGGAGCCCCCTGGCCGGTTCGGAGAGCCGGAT
AGTATAGGGTTCAGAATTTCTGAGCCCGGTCAAGCCAGCAGGGGCTTCGGGCCGAGCCCTCGGATCG
AGATCACTCCATCCACGAACTGATGCAGGCAGGGGGGCCCTCCGTGGGAGAGACGCCGGCCTGTCCCC
CGAGCAGCCGGCCCTGGCCCTGGCCGGCGTGGCCGCCAGCCGAGGTTACACTGCCCGTCCCGGCTAC
GAGGGCTACCGCAGCCGCTTTGCTTGAGCCCCGCTAGCAGCGGCTCCTCTGCCAGCTTCATTTCTGACA
CCTTCTCCCCCTACACCTCGCCCTGCGTCTCACCAATAACGCCGGGCCGAGCAGCTGTGTCCCCAGTT
TCAAAACATCCCTGCTCATTATCCCCCAGAACCTCTCCAATAATGTACCTCGAACCCAGCTCGCCGAG
GACAGCTGCCTGGGCCGACACTGCCCCGTGCCCGTCCGGCATCCCGCTCCTCTCACCCGGTGCCAAGC
GGAGGCATTTCGTGCGCAGAGGCTTTGGTTGCTCCTCTGCCCGCAGCCTCACCCGAGCGTCCCGGAGCCC
CTCGCCACAGCCCTCGCCTCACGTGGCACTGCAGGACGACAGCATCCCCGCTGGGTACCCCCACAGGCC
GGCTCTGCTGTTCTATGGATGCCCTCAACACCCTGGCCACCGACTCGCCCTGCGGGATCCCCTCCAAGA
TATGGAAGACCAGTCTGACCCGACGCTGTGTCCACCCTCCGTCCAAGGCTGGCTGGCCCGCCACAT
CTACCCTACTGTGGAGTTCCTGGGGCCATGTGAGCAGGAGGAGAGGAATTCCGCTCCAGAGTCCATC
CTGCTGGTACCACCTACTTGGCCCAAGCAGTTGGTGGCCGATTCATCTGCAGCATCCCTGTGACTG
CATCCCTCCCACCACTCGAGTGGCCACTCTCAATCAGTCGGGCTCCTATGAGTACGGATTGAGGTCCA
ACCCAAGCCCCATCACCGGGCCACTATGAGACGGAGGGCAGCCGTGGCGCTGTCAAAGCCCAACAGGA
GGACACCCTGTGGTGCAGCTCCACGGCTACATGGAGAACAAGCCTCTGGGGCTTCAGATTTTATTGGGA
CAGCAGATGAGAGGATCCTTAAGCCGCACGCCTTCTACCAAGTACACAGGATCACTGGGAAAACGGTCA
CACCACGAGCTATGAGAAGATCGTAGGCAACCAAGGTCCTGGAGATCCCCCTGGAGCCAAAGAACAAC
ATGAGAGCCACCATCGACTGTGCAGGCATCCTGAAGCTCCGAAACGCTGACATCGAGCTGCGGAAGGGCG
AGACGGACATCGGCAGGAAGAACACGCGTGTGCGCCTGGTGTCCGCGTGCACGTCCAGAGCCAGTGG
GCGCATCGTCTCCCTGCAGGCTGCGTCCAACCCATCGAGTGTCTCAGCGCTCTGCCACGAGCTGCC
ATGGTGGAGAGACAAGACATGGACAGCTGCCTGGTCTACGGGGCCAGCAGATGATCCTCACGGGCCAGA
ACTTCACAGCGGAGTCCAAGTTGTGTTTATGGAGAAGACTACAGGACCTGCAGGACCTGTGAGACTCG
GCCTTTGCCATCTCTCTGATCTCAGCTGACCGTCTCTCCCCCTGGCTTCCCGGCTGCAGAGAAACCT
CCTGGCTCTGCTTCAGATGCTCCGTGCTCCTCTGCCCCAGGATCTTCGTTGTGCTGCTGGCTCTG

ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGA
TTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230796 representing NM_001291172
 Red=Cloning site Green=Tags(s)


MQRPALRRPGPRQPLRTMGSADREEPIAHKAISSPSGLAYPDDVLDYGLKPCNPLASLSGEPPGRFGEPD
 SIGFQNFLSPVKPAGASGSPRIEITPSHELMQAGGALRGRDAGLSPEQPALALAGVAASPRFTLPVPGY
 EGYREPLCLSPASSGSSASFISDTFSPYTPCVPSPNAGPDDLCPQFQNIPIAHYSPRTSPIMSPRTSLAE
 DSCLGRHSPVPRPASRSSSPGAKRRHSACAEALVAPLPAASPQRSRSPSPQPSPHVALQDDSIIPAGYPTA
 GSAVLMDALNLTADSPCGIPSKIWKTS PDTPVSTAPSKAGLARHIYPTVEFLGPCEQEERRNSAPESI
 LLVPPTWPKQLVPAIPICSIPVTASLPLEWPLSNQSGSYELRIEVQPKPHHRAHYETEGSRGAVKAPTG
 GHPVVQLHGMYENKPLGLQIFIGTADERILKPHAFYQVHRITGKTVTTSYEKIVGNTKVL EIPLEPKNN
 MRATIDCAGILKLRNADIELRKGETDIGRKNTRVRLVFRVHVPEPSGRIVLSLQAASNPIECSQRSHELP
 MVERQDMSCLVYGGQQMILTGQNF TAESKVVFMKTTGPAGTCETRPLPISLISADRLSPWLSRLQRNP
 PGSVFRCSVLLPAPGSSLVLLAL


TRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-NotI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001291172

ORF Size: 1959 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291172.1](#), [NP_001278101.1](#)

RefSeq Size: 3283 bp

RefSeq ORF: 1962 bp

Locus ID: 18019

Cytogenetics: 2 88.91 cM

MW: 70.4 kDa

Gene Summary: Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2, IL-3, IL-4, TNF-alpha or GM-CSF. Promotes invasive migration through the activation of GPC6 expression and WNT5A signaling pathway.[UniProtKB/Swiss-Prot Function]