

## Product datasheet for **MR230880**

### **Dtna (NM\_001285810) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Dtna (NM_001285810) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Dtna
Synonyms:	2210407P21Rik; a-DB-1; adbn; Dtn; DTN-A; Gm19389
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>MR230880 representing NM\_001285810  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGATTGAAGATAGTGGAAAAGAGGAAACACCATGGCAGAAAGAAGACAGCTGTTTGCAGAGATGAGGG  
 CTAAGATCTGGATCGCATCCGACTCTCCACCTACAGAACAGCATGCAAGCTGAGATTTGTGCAGAAGAA  
 ATGCAATTTGCACCTGGTGGACATTTGGAACGTCAATTGAAGCATTGCGCGAAAACGCTTTGAATAACCTG  
 GACCCAAACATAGAACTCAACGTGGCCCGCTGGAGGCGGTGCTCTCCACTATTTTTTACCAGCTCAACA  
 AGAGGATGCCAACCACTACCAAATCCACGTGGAGCAGTCCATCAGTCTCCTGCTGAACTTCTGCTCGC  
 AGCCTTTGACCCGGAAGGCCATGGAAAAATCTCAGTATTTGCTGTCAAAATGGCTTTGGCTACATTGTGT  
 GGAGGAAAGATCATGGACAAGTTAAGATATATTTTCTCAATGATCTCTGACTCCAGTGGAGTGATGGTAT  
 ATGGAAGATATGACCAATTCCTTCGGGAAGTTCTCAAACACCCACAGCAGTTTTTGAAGGTCCTTCATT  
 TGTTTACACAGAACAGTCAGCCAGATCCTGTTTCTCCCAACAGAAAAGGTCACCTTAAATGGCTTCTTG  
 GACACGCTCATGTGACACCTCCTCCACAGTGCCTGGTATGGCTGCCTTCTTGCACCGACTGGCAAATG  
 TAGAAAATGTCTTCCATCCAGTTGAGTGTCTCTACTGTCACAGTGAAGCATGATGGGATTTTCGATACCG  
 ATGTCAACAGTGTCACACTACCAGCTTGGCAGGACTGCTTCTGGAGGGGCCATGCAGGAGGTTTCAT  
 AGCAACCAGCACCAAATGAAGGAGTACGTCTGGAAATCACCTGTAAAGAGCTAACGAATGCATTAA  
 GCAAGTCCCTGAGCTGTGCTTCAAGCCGTGAACCTTTCACCCCATGTTCCAGACCAGCCTGAAAAGCC  
 ACTCAACTGGCTCACATCGTTGATACTGGCCTCCGAGACCTGTAAACCAGCATGAACGACACACTTTT  
 TCCACTCTGTTCCCTCCTCAGGAAGTCCTTTCATCACCAGGAGTATGCTTGAGAGTTCAAACCGGCTTG  
 ATGAAGAACACAGGCTGATCGCCAGGTACGCCGACGACTGGCAGCAGAATCCTCCTCATCACAGCCAC  
 ACAGCAGAGGAGTGCCCTGACATCTCCTTACCATCGATGCAAATAAACAGCAAAGGCAGCTGATTGCA  
 GAGCTAGAGAACAAGAACAGAGAAATCTTACAAGAGATTAGAGACTTCGGGTAGAACATGAGCAAGCTT  
 CCCAGCCACGCCAGAGAAAGCTCAGCAGAACCAACCTTGTGGCAGAACTCCGGCTCCTCAGACAGCG  
 AAAGGATGAGCTAGAACAGAGAAATGTCTGCTCTCCAGGAGAGCCGGAGAGAGCTGATGGTCCAGTTGAA  
 GGTCTCATGAAGCTCCTGAAGACCCAGGGGCAAGCTCCACACGCTTCCCCAGCCACACCATCAGCA  
 GGCCGATCCCATGCCATCCGTTGAGCATCTGCCTGTCCACACCAACGCACACCCCGCAGGACTCCCT  
 CACTGGGGTAGGCGGAGATGTGCAAGAGGCGTTTGCACAAAGTTCGAGAAGGAATCTAAGGAGTGACCTG  
 CTGGTGGCTGCAGACTCCATACCAATACCATGTCTCTGTTGAGGAGCTGAATTCAGAGGTGGCGA  
 GTGAAACAGAAAGTACTGTGGATTCCGAGTTCTCGAGACCTCAGTTTGAAGGATCTGGCTCCATCACCAG  
 CTCTGAGAAGGCTTTTCTAGCCCAAATCCACTCCCGGAAACCTGGGTACATTACGGTGGAGCTGCAAGT  
 ACCACGATGGTGACATGGTCCCGGAAGATGGGGATCCTTACACGCAACCCGAGGATGGCAATTATGAGA  
 ACGAGTCTGTGCGGCAGCTGGAGAATGAGCTGCAGCTGGAGGAGTACCTGAAGCAGAAGCTGCAGGACGA  
 GGCCTACCAGGTCAGTTGCAAGGT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230880 representing NM\_001285810  
 Red=Cloning site Green=Tags(s)


MIEDSGKRGNTMAERRQLFAEMRAQDLDRIRLSTYRTACKLRFVQKKNLHLVDIWNVIEALRENALNNL  
 DPNIELNVARLEAVLSTIFYQLNKRMPPTHQIHVEQSISSLLNFLLAADFPEGHGKISVFAVKMALATLC  
 GGKIMDKLRYIFSMISDSSGVMVYGRYDQFLREVLKLP TAVFEGPSFGYTEQSARSCFSQQKQKVTLNGFL  
 DTLMSDPPPQCLVWLPLLHRLANVENVFHPVECSYCHSESMGFYRCQQCHNYQLCQDCFWRGHAGGSH  
 SNQHQMKEYTSWKSPAKKLTNALSKSLSCASSREPLHPMFDPQPEKPLNLAHIVDTWPPRPVTSMNDTLF  
 SHSVPSSGSPFITRSMLESSNRLDEEHRLIARYAARLAAESSSSOPTQQR SAPDISFTIDANKQQRQLIA  
 ELENKNREILQEIQRLRVEHEQASQPTPEKAQQNPTLLAELRLLRQRKDELEQRMSALQESRRELMVQLE  
 GLMKLLKTQGASSPRSSPHTISRPIPMPIRSASACPTPTHTPQDSL TGVGGDVQEAFQSSRRNLRSDL  
 LVAADSI TMTSSLVKELNSEVASETESTVDSEFSRPQFEDLAPSPTSEKAFLAQIHSRKP GYIHGGAAS  
 TTHGDMVPEDGDPYTPEDGNYENESVRQLENELQLEEY LKQKLQDEAYQVSLQG

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

SgfI-MluI

Cloning Scheme:



## Plasmid Map:



ACCN: NM\_001285810

ORF Size: 2055 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001285810.1](#), [NP\\_001272739.1](#)

RefSeq Size: 6219 bp

RefSeq ORF:	2058 bp
Locus ID:	13527
UniProt ID:	<u><a href="#">Q9D2N4</a></u>
Cytogenetics:	18 12.08 cM
MW:	77.8 kDa
Gene Summary:	Involved in synapse maturation and required for normal muscle function.[UniProtKB/Swiss-Prot Function]