

## Product datasheet for **MR230886**

### Pip5k1c (NM\_001293646) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Pip5k1c (NM_001293646) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Pip5k1c
Synonyms:	A1115456; A1835305; Pip5klgamma
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide  
Sequence:**

>MR230886 representing NM\_001293646  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGGAGCTAGAGGTGCCGACGAGGCGGAGAGCGCCGAAGCGGGGCGGTGACGGCGGAAGCGCCCTGGT  
 CTGCGGAGAGTGGGGCGGCGCAGGTATGACCCAGAAGAAGGCTGGCCTTGCAGAGGCCCTTGGTGAC  
 AGGGCAGCCAGGCCCTGGCCATGGGAAGAAGCTGGGCCACCGAGGCGTGGATGCGTCGGGAGAGACTACG  
 TATAAGAAGACCACCTCGTCCACCTGAAGGGCGCCATCCAGCTGGGGATCGGGTACACGGTGGGAACC  
 TGAGCTCTAAGCCGAGCGGGACGTGCTCATGCAGGACTTCTACGTGGTGGAGAGCATCTTCTCCCCAG  
 TGAAGGGAGCAACCTCACCCCTGCCACCATTCCAGGATTTCCGCTTCAAGACCTATGCACCTGTTGCC  
 TCCGCTACTTCCGGGAGCTCTTGGCATCCGTCAGATGACTATTTGACTCCCTGTGCAATGAGCCAC  
 TCATCGAGCTCTCAACCCGGGGCCAGCGGCTGTCTTCTACGTCACCAGCGACGACGAGTTCATCAT  
 CAAGACTGTCATGCACAAGGAGGACAGATTCTGCAGAAGCTGCTGCCTGGCTACTACATGAATCTCAAC  
 CAGAACCACGGACGCTGCTGCCAAGTTCTATGGGCTGTACTGCGTGCAGTCTGGTGGCAAGAATCC  
 GCGTGGTGGTCAATGAACAATGTGCTGCCCGTGTGCTCAAAATGCACCTTAAGTTCGACCTCAAGGGCTC  
 CACGTACAAGCGCAGGGCCAGCAAGAAGGAGAAGGAGAGCCTGCCACCTACAAGGACCTGGACTTC  
 ATGCAGGACATGCCCGAGGGGCTGCTGCTGGACTCCGACACCTTTGGCGCCCTGGTCAAGACGCTGCAGC  
 GAGACTGCCTGGTGTGGAGAGCTTCAAGATAATGGACTACAGCCTGCTGCTGGGCGTGCACAACATCGA  
 TCAGCAGGAGCGAGAGCGCCAGGCCGAGGGCGCCAGAGCAAGGCGGATGAGAAGCGGCCGTGGCCAG  
 AAGGCTCTGTATCCACGGCCATGGAGTCTATCCAGGGCGGAGCTGCCCGTGGGGAGGCCATTGAGACAG  
 ATGACAGATGGGTGGGATTCCAGCAGTGAATGGGCGCGGGGAGCGACTGCTCTCCACATCGGGATCAT  
 TGATATTCTGCAGTCTACAGTTTCATCAAGAAGTTAGAACACACCTGGAAGGCCCTCGTCCATGATGGG  
 GACACTGTCTCAGTCCACCGGCCAGTTCATGCAGAGCGCTTCTTCAAGTTCATGAGCAGCACGGTGT  
 TCCGGAAGAGTTCCTCCCTGAAGTCTCTCCATCCAAGAAAGGACGTGGTGCCTGCTGGCGGTCAAACC  
 CCTGGGGCCCACTGCTGCCTTCTCAGCCAGCCAGATCCCAGCGAGAGAGAGGATGTGCAGTATGACCTG  
 CGGGGGCCCGCAGTACCCACGCTTGAAGTGAAGGCCGGCTGACCTCCTGCCCTGCACCCACCGT  
 CCTTTGAGGAAGCCACCACAGCCTCCATCGCCACCACCCTGTCGTCACCTCCCTCTCCATCCCAGAGCG  
 GTCCCTTCAGATACATCGGAGCAGCCCGGTACAGGCGGCGCACGCAGTCTTACGGCCAGGATGGCCGG  
 CCCAGGAGGAGCCCATGCGGAAGACCTGCAGAAGATAACTGTGCAGGTGGAGCCAGTGTGCGGTGTGG  
 GGGTTGTCCCAAGGAGGAGGGTGCAGGAGTGGAGGTCCCCCATGTGGGCATCGGCTGCAGCCTCTGT  
 GAAATAGACGCTGCCAGCCAGGCTCAGAGCCTGCCAGCCAGGCTCAGATGAGGAGGATGCACCTCT  
 ACAGACATCTATTTTTTCGCCCATGGGAGATACTGGCTTTTCTCTCCCCGTGTCGCCAACTGCGGGCCG  
 TGACACCAAAACACAGGCACTCCACCGACGAGAGGAGCTGGGTGTAATCCCCGTTCACTATAGCGC  
 GCGGCCCGCTCCGACGGCGAGAGCGACACA


**ACGCGT**ACGCGGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA


Protein Sequence: >MR230886 representing NM\_001293646  
Red=Cloning site Green=Tags(s)

```
MELEVPDEAESAEAGAVTAEAAWSAESGAAAGMTQKKAGLAEAPLVTGQPGPGHGKKGKLGHRGVDASGETT
YKKTTSSTLKGAIQLGIGYTVGNLSSKPERDVLMDQDFYVVEISFFPSEGSNLTPAHHFQDFRFKTYAPVA
FRYFRELFGIRPDDYLYSLCNEPLIELSNPGASGSVFYVTSDDDEFIIKTMHKEAEFLQKLLPGYMNLN
QNPRTLTPKFYGLYCVQSGGKNIRVVVMNNVLPVVVMHLKFDLKGSTYKRRASKKEKESLPTYKDLDF
MQDMPEGLLLDSDTFGALVKTLQRDCLVLESFKIMDYSLLLGVHNIDQQERERQAEGAQSKADEKRPVAQ
KALYSTAMESIQGGAARGEAIETDDTMGGIPAVNGRGERLLHLHIGIIDILQSYRFIKKLEHTWKALVHDG
DTVSVHRPSFYAERFFKMSSTVFRKSSSLKSSPSKKGGRGALLAVKPLGPTAAFSASQIPSEREDVQYDL
RGARSYPTLEDEGRPDLLPCTPPSFEEATTASIATTLSTSLIPERSPDTSEQPRYRRRTQSSGQDGR
PQEPAEDLQKITVQVEPVCVGVPKKEEGAGVEVPPCGASAAASVEIDAASQASEPASQASDEEDAPS
TDIYFFAHGRYWLFSPRRRQLRAVTPNHTGTPDERSWVYSPLHYSARPASDGESDT
```

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:  
Cloning Scheme:



**Plasmid Map:**


**ACCN:** NM\_001293646

**ORF Size:** 2061 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001293646.1](#), [NP\\_001280575.1](#)

**RefSeq Size:** 4378 bp

RefSeq ORF: 2064 bp

Locus ID: 18717

UniProt ID: [O70161](#)

Cytogenetics: 10 C1

MW: 76 kDa

**Gene Summary:** Catalyzes the phosphorylation of phosphatidylinositol 4-phosphate (PtdIns4P) to form phosphatidylinositol 4,5-bisphosphate (PtdIns(4,5)P2). PtdIns(4,5)P2 is involved in a variety of cellular processes and is the substrate to form phosphatidylinositol 3,4,5-trisphosphate (PtdIns(3,4,5)P3), another second messenger. The majority of PtdIns(4,5)P2 is thought to occur via type I phosphatidylinositol 4-phosphate 5-kinases given the abundance of PtdIns4P. Participates in a variety of cellular processes such as vesicle mediated transport, cell adhesion, cell polarization and cell migration. Together with PIP5K1A is required for phagocytosis, but they regulate different types of actin remodeling at sequential steps. Promotes particle attachment by generating the pool of PtdIns(4,5)P2 that induces controlled actin depolymerization to facilitate Fc-gamma-R clustering. Mediates RAC1-dependent reorganization of actin filaments. Required for synaptic vesicle transport. Controls the plasma membrane pool of PtdIns(4,5)P2 implicated in synaptic vesicle endocytosis and exocytosis. Plays a role in endocytosis mediated by clathrin and AP-2 (adaptor protein complex 2). Required for clathrin-coated pits assembly at the synapse. Participates in cell junction assembly. Modulates adherens junctions formation by facilitating CDH1 trafficking. Required for focal adhesion dynamics. Modulates the targeting of talins (TLN1 and TLN2) to the plasma membrane and their efficient assembly into focal adhesions. Regulates the interaction between talins (TLN1 and TLN2) and beta-integrins. Required for uropodium formation and retraction of the cell rear during directed migration. Has a role in growth factor- stimulated directional cell migration and adhesion. Required for talin assembly into nascent adhesions forming at the leading edge toward the direction of the growth factor. Negative regulator of T-cell activation and adhesion. Negatively regulates integrin alpha-L/beta-2 (LFA-1) polarization and adhesion induced by T-cell receptor. Together with PIP5K1A has a role during embryogenesis and together with PIP5K1B may have a role immediately after birth. [UniProtKB/Swiss-Prot Function]