

## Product datasheet for **MR230909**

### **Kif3a (NM\_001290806) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Kif3a (NM_001290806) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Kif3a
Synonyms:	Kif3; Kifl; Kns3
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>MR230909 representing NM\_001290806  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGCATCGCC**

ATGCCGATCAATAAGTCGGAGAAGCCGGAAGCTGCGATAATGTGAAGGTGGTGGTTAGGTGCCGCGCTC  
TCAATGAGAGAGAGAAGTCAATGTGCTATAGACAGGCCGTGAGCGTGGATGAAATGAGGGGAACATCAC  
CGTCCATAAGACCGATTCTCCAACGAACCTCCAAGACATTTACTTTTCGATACTGTTTTGGACCAGAA  
AGTAAACAACCTTGATGTGTATACTTAAGTCAAGACCAATTATTGATTCTGTTCTGGAAGGCTACAATG  
GGACCATCTTTGCATATGGACAGACTGGGACAGGCAAGACTTTCACCATGGAAGGCGTTCGAGCAGTACC  
GGGGCTCAGAGGGGTATCCCAACTCATTGCTCACATATTTGGTCACATCGCAAAAGCAGAGGGCGAC  
ACAAGTTTTTTGGTTCGAGTGTACTTGGAGATATATAATGAAGAAGTTCGAGACCTTTGGGTAAAG  
ATCAGACCCAGAGGTTAGAGGTTAAAGAAAGCCCTGATGTGGGAGTATATCAAAGACTTATCAGCATA  
TGTAGTCAACAACGCTGACGACATGGATAGAATCATGACTCTAGGCCACAAAACCGTTCTGTCGGTGCG  
ACTAATATGAACGAGCACAGCTCCCGATCCACGCCATCTTACAATTACTATAGAGTGCAGTGAGAAAG  
GCGTGGACGGGAACATGCACGTCAGGATGGGCAAGCTCCACCTTGTGGATCTGGCTGGGTGAGAAAGACA  
AGCAAAAACCTGGAGCCACAGGACAGCGCTCAAGGAAGCTACAAAATCAACCTTTCATTGTCTACCCCTG  
GGCAACGTTATTTCTGCCTTGGTGGATGGAAAAAGCACTCACGTGCCGTATCGGAACTCCAACTGACCC  
GTCTCTCGCAGGACTCCTTAGGAGGGAACCAAAAACCATGATGTGTGCAAAATTTGGGCCAGCAGATTA  
TAATATGATGAAACTATCAGCACCTTGGGATGCCAACCAGCAAGAACATTAAGAAAGCAAGCGAGA  
ATCAATGAGGACCAAGGATGCTCTGCTCCGCCAGTTTCAGAAAGAAATCGAAGAGCTGAAGAAAAAGC  
TTGAGGAAGGGGAAGAAGTGTGAGGCTCTGACATCAGTGGGTGAGAGGAGGATGACGAGGAGGTTGAAT  
TGGAGAAGATGGAGAAAAAAGAAAAAGAGAAGGACCAAGCAGATAAGATGGTGGAAATGCAAGCGAAA  
ATTGACGAGGAGAGAAAAGCACTGGAAACCAACTTGACATGGAAGAAGAAAGAAACAAAGGCTAGAG  
CTGAACTAGAGAGACGGGAGAAGGACCTTCTGAAAGCCCAACAAGAGCATCAGTCTTTGCTGAAAAAGCT  
GTCTGCTTTGAAAAAGAGGTCATTGTGCGCGCGTGGACCTGCTGGCAAGGCTGAGGAGCAAGAGAAG  
CTTCTGGAGGAGTCCAACATGGAGCTGGAGGAGAGGAGGAGACGAGCTGAGCAGCTTCGAAAGAGCTGG  
AGGAAAAGGAGCAAGAACGCTTGGACATCGAGGAAAAATATACCAGTCTGCAAGAGGAGGCACAAGGAAA  
GACCAAGAAATTAAGAAAGTCTGGACCATGCTGATGGCTGCAAAGTCAAGATGGCTGATCTCCAGCAG  
GAGCATCAGAGGAGATCGAAGGCCTCCTGGAGAATCCGGCAGCTCAGCCGGGAGCTCCGTCTTCAGA  
TGCTCATCATTGATAACTTCATACCTCAGGATTACCAGGAAATGATTGAAAACACTGCTCCACTGGAATGA  
AGACATAGGAGAATGGCAGCTAAAATGTGTGGCCTACACAGGAAATAACATGAGGAAACAACTCCAGTC  
CCTGACAAGAAGGAGAGAGACCCATTTGAAGTCGACCTTTCCACGTGTACCTCGCCTATACAGAGGAGA  
GTCTGCGTCAGTCTGTGATGAAGTTAGAAAGGCCACGGACCTCCAAGGGGAAAGCAAGGCCAAAGATGGG  
CGGAAGAAAGCGTTCTGCAAAGCCTGAGACCGTAATTGATTCTTTACTTCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Protein Sequence:** >MR230909 representing NM\_001290806  
 Red=Cloning site Green=Tags(s)


MPINKSEKPECDNVKVVVRCRPLNEREKSMCYRQAVSVDEMRTITVHKTDSSNEPPKTFTFDVFVGP  
 SKQLDVYNTARPIIDSVLEGYNGTIFAYGQTGTGKFTMEGVRAVPLRGVIPNSFAHIFGHIKAEGD  
 TRFLVRVSYLEIYNEEVRDLLGKDQTRLEVKERPDVGVYIKDL SAYVVNNADDMDRIMTLGHKNRSVGA  
 TNMNEHSSRSHAI FTITIECSEKGV DGNMHV RMGKLHLV DLAGSERQAKTGATGQRLKEATKINLSLSTL  
 GNVISALVDGKSTHVPYRNSKLTRLLQDSLGGNSKTMMCANIGPADYNYDETISTLRYANRAKNIKAR  
 INEDPKDALLRQFQKEIEELKKKLEEGEEVSGSDISGSEEDDEEGELGEDGKKKRRDQADKMVEMQAK  
 IDEERKALETKLDMEEEEERNKARAELEERREKDLLKAQEQEHQSLLEKLSALEKKVI VGGVDLLAKAEQEK  
 LLEESNMELEERRRRAEQLRKELEEKERLDIEEKYTSLQEEAQGKTKLKKVWMLMAAKSEMADLQQ  
 EHQREIEGLLENIRQLSRELRLQMLIIDNFIPQDYQEMIENYVHWNEDIGEWQLKCVAYTGNNMRKQTPV  
 PDKKERDPFEVDLSHVYLAYTEESLRQSLMKLERPRTSKGKARPKMGRKRSAKPETVIDSLLQ


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

**Restriction Sites:**

SgfI-MluI

**Cloning Scheme:**



**Plasmid Map:**


**ACCN:** NM\_001290806

**ORF Size:** 2082 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001290806.1](#), [NP\\_001277735.1](#)

**RefSeq Size:** 5429 bp

RefSeq ORF: 2085 bp

Locus ID: 16568

UniProt ID: [P28741](#)

Cytogenetics: 11 31.97 cM

MW: 79.9 kDa

**Gene Summary:** Microtubule-based anterograde translocator for membranous organelles. Plus end-directed microtubule sliding activity in vitro. Plays a role in primary cilia formation (PubMed:21670265). Plays a role in centriole cohesion and subdistal appendage organization and function. Regulates the formation of the subdistal appendage via recruitment of DCTN1 to the centriole. Also required for ciliary basal feet formation and microtubule anchoring to mother centriole (PubMed:23386061).[UniProtKB/Swiss-Prot Function]