

## Product datasheet for **MR230916**

### **Tle1 (NM\_001285530) Mouse Tagged ORF Clone**

#### **Product data:**

Product Type:	Expression Plasmids
Product Name:	Tle1 (NM_001285530) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Tle1
Synonyms:	C230057C06Rik; Estm14; Grg-1; Grg1; Tle4l
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide  
Sequence:

>MR230916 representing NM\_001285530  
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGTTCCCGCAGAGCCGGCACCCAACGCCGCACCAAGCTGCAGGCCAGCCCTTTAAGTTCACTATCCCGG  
AGTCTCTGGACCGGATTAAGAGGAATTCAGTTCCTGCAGGCGCAGTATCACAGTCTTAAATTGGAGTG  
TGAGAACTGGCAAGTGAAGACAGAAATGCAGAGGCACTACGTGATGTATTATGAAATGTCATATGGA  
TTAAACATTGAAATGCACAAACAGACTGAAATCGCCAAGAGATTGAACACCATTTGTGCCAAAGTCATCC  
CATTTCTGTCTCAGGAACATCAACAACAGGTGGCCAGGCTGTGGAACGTGCCAAACAGGTGACCATGGC  
AGAGTTGAATGCCATCATCGGGAGTAATTCCTCTTGGTCCCAGACAGCCTAAGAGGCACAGATAAGCGC  
AGAAATGGTCCAGAGTTTTCCAGTGACATCAAAAAAGGAAGGTGGATGATAAGGATAACTATGACAGTG  
ATGGGGACAAGAGTGATGACAACCTAGTTGTGGATGTGTCTAACGAGGACCCCTTCTCCACATGCAAG  
CCCCACACTCACCCCGGAAAACGGAATTGACAAAAACCGTCTGCTGAAGAAAGATGCTTCAGGTAGC  
CCGGCATCCACAGCCTCCTCTGGAAGTTCTCTCCCTGAAATCCAAAGAAGTGAGCCTGCATGAAAAAG  
CCAACACTCCTGTTCTGAAATCCAGCACACCGACGCCTCGGAGCGACATGCCAACCCAGGCACCCAGCGC  
TACTCCAGGCCTCCGTCCAGGTCTTGGCAAGCCTCCAGCCATGGAACCCCTTGTCAACCAAGCAGCAGCT  
GGCCTGAGGACGCCCTGGCGGTGCCTGGCCCATACCCTGCCCTTTGGCATGGTGCCCATGCGGGCA  
TGAACGGAGAGCTGACCAGCCCTGGTGTGCCTATGCAGGTCTACACAGCATGTCTCCACAGATGAGCGC  
TGCAGCTGCTGCAGCTGTGTCTGTGGTGGCCTATGGGCGTCCCAATGGTTGGTTTTGATCCTCCT  
CCTCACATGAGAGTACCTTCTATCCCCCAACCTGGCAGGAATACCTGGAGGAAACCAGCATACTCCT  
TCCACCTTACTGCTGATGGCCAAATGCAGCCTGTTCTTTTCCCCTGATGCCCTCATTGGACCTGGGAT  
TCCCCGACATGCTCGGCAGATCAACACCCCTCAACCACGGGGAGGTGGTGTGTGCAGTGACCATCAGCAAC  
CCCACAAGGCACGTGTACACAGGTGGCAAGGGCTGCGTTAAGGTGTGGGACATCAGCCACCCCGCAACA  
AGAGCCCCGTCTCTCAGCTGGATTGTCTGAATAGAGATAACTACATCCGTTCTGTAAATGCTACCTGA  
TGGCTGTACTCTCATAGTGGGAGGGGAAGCCAGTACTCTATCCATCTGGGACTTGGCAGCTCCAACCTCC  
CGCATAAAGGCAGAGCTGACATCCTCAGCCCCTGCCTGTTATGCTCTGGCCATCAGCCCCGACTCCAAGG  
TCTGCTTCTCATGCTGCAGTGACGGCAACATCGCAGTGTGGACCTGCACAACCAGACCCTGGTGGGCA  
ATTCCAGGGCCACACAGACGGAGCCAGCTGTATTGACATTTCTAATGATGGCACCAAGCTCTGGACAGGC  
GGTTTAGACAACACTGTGAGGTCTGGGACCTGAGAGAAGGGAGGCAGCTGCAGCAGCATGACTTCACTT  
CACAGATCTTCTTTGGGATATTGCCCAACTGGGGAGTGGCTTCTGTGGGCATGGAAAGCAGCAACGT  
GGAAGTTCTGCATGTGAACAAGCCTGACAAGTACCAGCTGCACCTCCACGAGAGCTGCGTGTGTCCCTA  
AAGTTTGCTTATTGTGGCAAATGGTTTGTGAGTACTGGAAAAGACAACCTCCTAAATGCCTGGCGGACAC  
CATATGGAGCCAGCATCTCCAGTCCAAGAATCTTCATCAGTGCTTAGCTGTGACATCTCTGTGGATGA  
TAAATACATAGTCACTGGCTCAGGGGACAAAAAGGCTACCGTTTATGAGGTCATCTAC

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230916 representing NM\_001285530  
 Red=Cloning site Green=Tags(s)


MFPQSRHPTPHQAAGQPFFKFTIPESLDRIKEEFQFLQAQYHSLKLECEKLA SEKTEMQRHYVMMYEMSYG  
 LNIEMHKQTEIAKRLNTICAQVIPFLSQEHQQQVAQVERAKQVTMAELNAIIGSNLLVPDSL RGTDKR  
 RINGPEFSSDIKKRKVDDKNDYSDGDKSDDL VVDVSNEDPSSPHASPTHSPRENGIDKNRLLKDDASGS  
 PASTASSGSSSSLSKEVSLHEKANTPVLKSSSTPTPRSDMPTPGTSATPGLRPLGKPPAMEPLVNQAAA  
 GLRTPLAVPGYPAPFGMVPHAGMNGELTSPGAAAYAGLHSMSPQMSAAAAAAAAAVVAYGRSPMVGFDPP  
 PHMRVPSIPPNLAGIPGGKPAYSFHVTADGQM QVPVFPDALIGPGIPRHARQINTLNHGEVVCVAVTISN  
 PTRHYYTGGKGCYKVDI SHPGNKSPVSQLDCLNRDNYIRSKLLPDGCTLIVGGEASTLSIWDLAAPT  
 RIKAELTSSAPACYALAI SPDSKVCFSCCSDGNI AVWDLHNQTLVRQFQGHTDGASCIDISNDGTLWGT  
 GLDNTVRSWDLREGRQLQHDFTSQIFSLGYCPTGEWLAVGMESSNVEVLHVNKPKDYQLHLHESCVLSL  
 KFAYCGKWFVSTGKDNLLNAWRTPYGASIFQSKESSSVLSCDISVDDKYIVTGS GDKKATVYEVII


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-MluI

Cloning Scheme:



**Plasmid Map:**


**ACCN:** NM\_001285530

**ORF Size:** 2088 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001285530.1](#), [NP\\_001272459.1](#)

**RefSeq Size:** 11328 bp

RefSeq ORF: 2091 bp

Locus ID: 21885

UniProt ID: [Q62440](#)

Cytogenetics: 4 C2

MW: 75.8 kDa

**Gene Summary:** Transcriptional corepressor that binds to a number of transcription factors. Inhibits NF-kappa-B-regulated gene expression. Inhibits the transcriptional activation mediated by FOXA2, and by CTNNB1 and TCF family members in Wnt signaling. The effects of full-length TLE family members may be modulated by association with dominant-negative AES. Unusual function as coactivator for ESRRG (By similarity).[UniProtKB/Swiss-Prot Function]