

Product datasheet for **MR230917**

Cdc27 (NM_001285990) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cdc27 (NM_001285990) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Cdc27
Synonyms:	AI452358; APC3; BC023187
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>MR230917 representing NM_001285990
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACGACCTTGTCAGTGTGGAGATTGAGCTTTCCTTGTGGACATGTGTATTG
GTGAGAAGCCAGATCCTGACCAAACATTTAAATTAACATCTCTACAGAATTTAGCAGTTGTCTCCCAA
CACTTGTACAACCTAGTATCTAATCACAGTTTATCTCACAGACAGCCTGAGACAGTCTTACAGAAACA
CCCCAAGACACAATTGAATTAACAGACTGAATTTAGAATCTTCCAATTCAAAGTACTCCTTGAATACAG
ATTCTCCGTGTCTTACATTGATTCAACTGTAATTTACCAGATAACGTCCCCCTTGACCTGGTACTGC
CATATTATCTAAACAGGTTCAAATAAACCAAAAACTGGTGAAGTTTATTAGGAGGACCAACTGCTCTT
AGCCCATTAACCCCAAGTTTGGGATTTGCCATTAGAAACCCCAAGTCTGGAGATGGATCCTATTTAC
AAAACACTACTAATACACCTTCTGTAATTGATGTGGCACCCACCGAGCACCTACAAAAAGACTTTTTG
TGTTTTACAGTCTGTTGAAGAATGGGCCAAAAGTGAACAAAGTCTGTCTCTCACAGAGTGAAATAGT
CGAGAGGTCACACCACTGTTGTCACAAACACAAAGTTCTGGCCACAAACAAGTACAACACCTCAGG
TATTGAGCCCCACTATCACATCTCCCCAAACGATTGCCTCGGAGAAGTCCCGCCTTTTCACTAGTGA
CAGTTCTACAACCAAGGAGAATAGCAAAAAGTTAAAAATGAAGTTTCCACCTAAAATCCCTAACGAAAA
ACAAAAAGTAAAATAATAAGGAGGACTAACTCAGCCCAGCATAAACGATAGTCTGGAAATTACAAAA
TGGACTCCTCTATCATTTCAGAAGGGAAAATAACCACAGTCACACCTCAGATCCAGGCATTTAACCTACA
AAAGGCAGCAGCAGAAGGTTTGTGAGCCTTCTCGTGAATGGGAAAAGGTTATTTAGCTTTGTGTTCA
TACAACTGCAAAGAAGCTATCAATTTTTGAGCCATCTACCTTCTCATCACTACAGTACTGGTGGGTCC
TATGCCAGATTGGACGGGCTTATTTGAACTTTCAGAATACATGCAGGCTGAAAGAATATTCTCAGAGGT
TAGAAGGATTGAGAGTTTCAGAGTTGAAGGAATGGAGATCTACTCTACAACACTCTGGCATCTTCAGAAA
GACGTTGCTCTTTCAGTCTTTCAAAGATTTAACAGACATGGATAAGAATTCACCAGAGGCCTGGTGTG
CTGCAGGGAAATTGTTTCAGTCTACAACGAGAGCATGATATAGCAATTAATTTCTTCAAAGAGCTATCCA
GGTCGATCCAAATTACGCTTACGCCTACTCTATTAGGACATGAGTTTGTGTTAACTGAAGAAGTATAGT
AAAGCATTAGCATGTTTTCGAAATGCTATAAGAGTAAATCCCAGACATTACAATGCATGGTATGGTTTAG
GAATGATTTATTACAAGCAAGAGAAGTTCAGCCTCGCAGAAATGCATTTCCAGAAAGCACTTGATATCAA
CCCTCAGAGTTTCACTTTGCCACATTGGAGTAGTTCAGCATGCACTAAAGAAGTCTGAGAAGGCT
TTGGATACCTAAACAAAGCCATTGTTATCGATCCCAAGAACCCTCTATGCAAATTTACAGAGCCTCGG
TTTTATTTGCAAATGAAAAATACAAGTCTGCTTTACAAGAACTCGAGGAACTGAAGCAGATTGTCCCAA
AGAGTCCCTCGTTTACTTCTTAATAGGAAAGTTTACAAGAAGTTGGGTCAAACCTCATCTCGCCTTGATG
AATTTCTCTTGGGCTATGGATTTAGATCTAAAGGAGCCAATAATCAGATTAAGAGGCAATTGACAAGC
GCTACCTTCCAGATGATGAGGAGCCAATAACCCAGGAGGAACAGATCATGGGGACAGATGAATCCCAGGA
GAGCAGCATGACAGATGCAGATGACACACAACCTCATGCGGCGGAAAGTGACGAATTT

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230917 representing NM_001285990
 Red=Cloning site Green=Tags(s)


MTTLLSLSLEIQLASL FPCWDMCIGEKPPDPQTFKLTSLQNFSSCLPNTCTTLVSNHLSHRQPETVLTET
 PQDTIELNRLNLESSNSKYSLNTDSSVSYIDSTVISPDNVPLGPGTAILSKQVQNKPKTGRSLLGGPTAL
 SPLTPSFGILPLETPSPGDGSYLQNYTNTPSVIDVAPTGAPTKKTF CVLQSVARMGQTGTKSVFSQSGNS
 REVTPVLVAQTQSSGPQTSTTPQVLSPTITSPPNALPRRSRLFTSDSSTTKENSKLKMKFPPKIPNRK
 TKSKTNKGLTQPSINDSLEITKLDSSIISEGIITTVTPQIQAFNLQKAAAEGLMSLLREMGKGYLALCS
 YNCKEAINILSHLPSHHYSTGWVLCQIGRAYFELSEYMQAERIFSEVRRIESFRVEGMEIYSTTLWHLQK
 DVALSVLSKDLTDMDKNSPEAWCAAGNCFSLQREHDIAIKFFQRAIQVDPNYAYAYTLLGHFVLT EELD
 KALACFRNAIRVNPRHYNAWYGLGMIYYKQEKFLAEMHFQKALDINPQSSVLLCHIGVVQHALKKSEKA
 LDTLNKAIVIDPKNPLCKFHRSVLFANEKYKSALQELEELKQIVPKESLVYFLIGKVKYKLGQTHLALM
 NFSWAMDLDPKGANNQIKEAIDKRYLPDDEEITQEEQIMGTDESQESSMTDADDTQLHAAESDEF

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001285990

ORF Size: 2088 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001285990.1](#), [NP_001272919.1](#)

RefSeq Size: 5703 bp

RefSeq ORF: 2091 bp

Locus ID: 217232

Cytogenetics: 11 67.79 cM

MW: 77.8 kDa

Gene Summary: Component of the anaphase promoting complex/cyclosome (APC/C), a cell cycle-regulated E3 ubiquitin ligase that controls progression through mitosis and the G1 phase of the cell cycle. The APC/C complex acts by mediating ubiquitination and subsequent degradation of target proteins: it mainly mediates the formation of 'Lys-11'-linked polyubiquitin chains and, to a lower extent, the formation of 'Lys-48'- and 'Lys-63'-linked polyubiquitin chains (By similarity). [UniProtKB/Swiss-Prot Function]