

Product datasheet for **MR230929**

Nfatc2 (NM_001291174) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfatc2 (NM_001291174) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfatc2
Synonyms:	AI607462; NF-ATc2; NF-ATp; NFAT1; NFAT1-D; Nfatp
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>MR230929 representing NM_001291174
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCACCTCGAACCAGCCTCGCCGAGGACAGCTGCCTGGGCCGACACTCGCCCGTGCCCGTCCGGCAT
CCCGCTCCTCCTCACCCGGTGCCAAGCGGAGGCATTCTGTGCGCAGAGGCTTTGGTTGCTCCTCTGCCCGC
AGCCTCACCCAGCGCTCCCGGAGCCCTCGCCACAGCCCTCGCCTCACGTGGCACTGCAGGACGACAGC
ATCCCCGCTGGGTACCCCCACGGCCGGCTCTGCTGTTCTCATGGATGCCCTCAACACCCTGGCCACCG
ACTCGCCTGCGGGATCCCTCCAAGATATGGAAGACCAGTCTGACCCGACGCCTGTGTCCACCCTCC
GTCCAAGGCTGGCCTGGCCCGCCACATCTACCCTACTGTGGAGTTCTGGGGCCATGTGAGCAGGAGGAG
AGGAGGAATTCGCTCCAGATCCATCTGCTGGTACCACCTACTTGGCCCAAGCAGTTGGTGCCGGCCA
TTCCCATCTGCAGCATCCCTGTGACTGCATCCCTCCCACCACTCGAGTGGCCACTCTCCAATCAGTCGGG
CTCCTATGAGCTACGGATTGAGGTCCAACCCAAGCCCCATCACCGGGCCCACTATGAGACGGAGGGCAGC
CGTGGCCCTGTCAAAGCCCCAACAGGAGGACACCCTGTGGTGCAGCTCCACGGCTACATGGAGAACAAGC
CTCTGGGGCTTACAGATCTTATTGGGACAGCAGATGAGAGGATCCTTAAGCCGCACGCCTTCTACCAAGT
ACACAGGATCACTGGGAAAACGGTCACCACCACGAGCTATGAGAAGATCGTAGGCAACACCAAGGTCCTG
GAGATCCCCCTGGAGCCAAGAACAACATGAGAGCCACCATCGACTGTGCAGGCATCCTGAAGTCCGAA
ACGCTGACATCGAGCTGCGGAAGGGCGAGACGGACATCGGCAGGAAGAACACGCGTGTGCGCCTGGTGT
CCGCGTGCACGTCCCAGAGCCAGTGGGCGCATCGTCTCCCTGCAGGCTGCGTCCAACCCCATCGAGTGC
TCTCAGCGCTCTGCCCACGAGCTGCCATGGTGGAGAGACAAGACATGGACAGCTGCCTGGTCTACGGGG
GCCAGCAGATGATCCTCACGGGCCAGAACTTACAGCGGAGTCCAAGTTGTGTTTATGGAGAAGACTAC
AGATGGGCAGCAGATTTGGGAGATGGAAGCTACGGTGGATAAAGACAAGAGCCAGCCTAACATGCTTTTT
GTTGAGATCCCGAGTATCGGAACAAGCACATCCGCGTGGCCGTGAAAGTCAACTTCTACGTCATCAACG
GAAAGAGGAAACGAAGTCAGCCACAGCACTTTACCTACCACCCAGTCCCTGCCATCAAGACAGAGCCAG
CGATGAGTATGAACCATCTTTGATCTGCAGCCCCGCCATGGAGGCTGGGGAGCCAGCCATATTACCCA
CAGCACCAATGCTGGCCGAGTCCCCCTCCTGCCTTGTGGTACCATGGCCCCCTGCCAACAGTTCGCT
CGGGGCTCTCATCCCCGATGCTCGCTACCAACAGCAGAGCCCCGAGCTGCCCTCTACCAGAGAAGCAA
GAGCCTGAGTCCCGCCTGCTGGGCTACCAGCAGCCGTCCTCCTGGCAGCACCTTGGGTCTGGCTGAT
GCCACCGCTCTGTGCTGGTGCATGCTGGTTCTCAGGGGCAGGGCAGGGCTCCACCCTGCCACACACAT
CCTCGGCCAGCCAGCAGGCTCACCCGTGATCCACTACTCACCCACCAACCAGCAGCTTCGCGGTGGGGG
TCACCAGGAGTTCAGCATATCATGTACTGTGAAAACCTCGGCCCCAGCTCTGCCAGGCCTGGCCCGCT
CCCATCAACCAAGGTCAGAGGCTGAGCCCGGGCGCCTACCCACAGTCATCCAACAACAGACTGCCCCGA
GCCAAAGAGCTGCCAAAACGGACCCAGTGACCAGAAGGAAGCTTGCCACGGGAGTGACCGTCAAACA
GGAACAGAACCTGGACCAGACCTACCTGGATGACGAGTTGATAGACACACACCTTAGCTGGATACAAAAC
ATATTA

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230929 representing NM_001291174
 Red=Cloning site Green=Tags(s)


MSPRTSLAEDSCLGRHSPVPRPASRSSSPGAKRRHSACAEALVAPLPAASPQRSRSPSPQPSPHVALQDDS
 IPAGYPPTAGSAVLM DALNTLATDSPCGIPSKIWKTS PDPTPVSTAPSKAGLARHIYPTVEFLGPCEQEE
 RRNSAPESILLVPPTWPKQLVPAIPICSI PVTASLPLEWPLSNQSGSYELRIEVQPKPHHRAHYETEGS
 RGAVKAPTGGHPVVQLHGMYENKPLGLQIFIGTADERILKPHAFYQVHRITGKTVTSTTSYEKIVGNTKVL
 EIPLEPKNNMRATIDCAGILKLRNADIELRKGETDIGRKNTRVRLVFRVHVPEPSGRIVSLQAASNPIEC
 SQRSAHELPMVERQDMD SCLVYGGQQMILTGQNF TAESKVVFMEKTTDGGQI WEMEATVDKDKSQPNMLF
 VEIPEYRNKHIRVPVKVNFYVINGKRKRSQPQHFTYHPVPAIKTEPSDEYEPSLICSPA HGLG SQPYYP
 QHPMLAESPSCLVATMAPCQQFRSGLSSPDARYQQQSPAAALYQRSKSLSPGLLGYQQPSLLAAPLGLAD
 AHRSVLVHAGSQGGQGSTLPHTSSASQQASPIHYSPTNQQLRGGGHQEFQHIMYCENFGPSSARPGPP
 PINQGQRLSPGAYPTVIQQQTAPSQRAAKNGPSDQKEALPTGVTVKQEQLDQTYLDDELIDTHLSWIQN
 IL


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-RsrII

Cloning Scheme:



Plasmid Map:


ACCN: NM_001291174

ORF Size: 2106 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291174.1](#), [NP_001278103.1](#)

RefSeq Size: 6210 bp

RefSeq ORF: 2109 bp

Locus ID:	18019
Cytogenetics:	2 88.91 cM
MW:	77 kDa
Gene Summary:	Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2, IL-3, IL-4, TNF-alpha or GM-CSF. Promotes invasive migration through the activation of GPC6 expression and WNT5A signaling pathway.[UniProtKB/Swiss-Prot Function]