

Product datasheet for **MR230930**

Nfatc2 (NM_001291175) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Nfatc2 (NM_001291175) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Nfatc2
Synonyms:	AI607462; NF-ATc2; NF-ATp; NFAT1; NFAT1-D; Nfatp
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR230930 representing NM_001291175
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGTCACCTCGAACCAGCCTCGCCGAGGACAGCTGCCTGGGCCGACACTCGCCCGTGCCCGTCCGGCAT
 CCCGCTCCTCCTCACCCGGTGCCAAGCGGAGGCATTCTGTGCGCAGAGGCTTTGGTTGCTCCTCTGCCCGC
 AGCCTCACCCAGCGCTCCCGGAGCCCTCGCCACAGCCCTCGCCTCACGTGGCACTGCAGGACGACAGC
 ATCCCCGCTGGGTACCCCCACGGCCGGCTCTGCTGTTCTCATGGATGCCCTCAACACCCTGGCCACCG
 ACTCGCCTGCGGGATCCCTCCAAGATATGGAAGACCAGTCTGACCCGACGCCTGTGTCCACCCTCC
 GTCCAAGGCTGGCCTGGCCCGCCACATCTACCCTACTGTGGAGTTCCTGGGGCCATGTGAGCAGGAGGAG
 AGGAGGAATTCGCTCCAGATCCATCTGCTGGTACCACCTACTTGGCCCAAGCAGTTGGTGCCGGCCA
 TTCCCATCTGCAGCATCCCTGTGACTGCATCCCTCCCACCACTCGAGTGGCCACTCTCCAATCAGTCGGG
 CTCCTATGAGCTACGGATTGAGGTCCAACCCAAGCCCCATCACCGGGCCCACTATGAGACGGAGGGCAGC
 CGTGGCCCTGTCAAAGCCCCAACAGGAGGACACCCTGTGGTGCAGCTCCACGGCTACATGGAGAACAAGC
 CTCTGGGGCTTACAGATCTTATTGGGACAGCAGATGAGAGGATCCTTAAGCCGCACGCCTTCTACCAAGT
 ACACAGGATCACTGGGAAAACGGTCACCACCACGAGCTATGAGAAGATCGTAGGCAACACCAAGGTCCTG
 GAGATCCCCCTGGAGCCAAGAACAACATGAGAGCCACCATCGACTGTGCAGGCATCCTGAAGTCCGAA
 ACGCTGACATCGAGCTGCGGAAGGGCGAGACGGACATCGGCAGGAAGAACACGCGTGTGCCCTGGTGT
 CCGCGTGCACGTCCCAGAGCCAGTGGGCGCATCGTCTCCCTGCAGGCTGCGTCCAACCCCATCGAGTGC
 TCTCAGCGCTCTGCCCACGAGCTGCCATGGTGGAGAGACAAGACATGGACAGCTGCCTGGTCTACGGGG
 GCCAGCAGATGATCCTCACGGCCAGAACTTACAGCGGAGTCCAAGTTGTGTTTATGGAGAAGACTAC
 AGATGGGCAGCAGATTTGGGAGATGGAAGCTACGGTGGATAAAGACAAGAGCCAGCCTAACATGCTTTTT
 GTTGAGATCCCGAGTATCGGAACAAGCACATCCGCGTCCCGTCAAAGTCAACTTCTACGTCATCAACG
 GAAAGAGGAAACGAAGTCAGCCACAGCACTTTACCTACCACCCAGTCCCTGCCATCAAGACAGAGCCAG
 CGATGAGTATGAACCATCTTTGATCTGCAGCCCCGCCATGGAGGCTGGGGAGCCAGCCATATTACCA
 CAGCACCAATGCTGGCCGAGTCCCCCTCCTGCCTTGTGGTACCATGGCCCCCTGCCAACAGTTCGCT
 CGGGGCTCTCATCCCCGATGCTCGCTACCAACAGCAGAGCCCCGAGCTGCCCTCTACCAGAGAAGCAA
 GAGCCTGAGTCCCGCCTGCTGGGCTACCAGCAGCCGTCCTCCTGGCAGCACCTTGGGTCTGGCTGAT
 GCCCACGCTCTGTGCTGGTGCATGCTGGTTCTCAGGGGCAGGGCAGGGCTCCACCCTGCCACACACAT
 CCTCGGCCAGCCAGCAGGCTCACCCGTGATCCACTACTCACCCACCAACCAGCAGCTTCGCGGTGGGGG
 TCACCAGGAGTTCAGCATATCATGTAAGTGTGAAAACCTCGGGCCAGCTCTGCCAGGCCTGGCCCGCT
 CCCATCAACCAAGGTCAGAGGCTGAGCCCGGGCGCTACCCACAGTCATCCAACAACAGACTGCCCCGA
 GCCAAAGAGCTGCCAAAACGGACCCAGTGACCAGAAGGAAGCTTGCCACGGGAGTGACCGTCAAACA
 GGAACAGAACCTGGACCAGACCTACCTGGATGACGAGTTGATAGACACACACCTTAGCTGGATACAAAAC
 ATATTA

AG**CGGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
 TGGATTACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230930 representing NM_001291175
 Red=Cloning site Green=Tags(s)


MSPRTSLAEDSCLGRHSPVPRPASRSSSPGAKRRHS CAEALVAPLPAASPQRSRSPSPQPSPHVALQDDS
 IPAGYPPTAGSAVLM DALNTLATDSPCGIPSKIWK TSPDTPVSTAPSKAGLARHIYPTVEFLGPCEQEE
 RRNSAPESILLVPPTWPKQLVPAIPIC SIPVTASLPLEWPLSNQSGSYELRIEVQPKPHHRAHYETEGS
 RGAVKAPTGGHPVVQLHG YMENKPLGLQIFIGTADERILKPHAFYQVHRITGKT VTTTSYEKIVGNTKVL
 EIPLEPKNNMRATIDCAGILKLRNADIELRKGETDIGRKNTRVRLVFRVHVPEPSGRIVSLQAASNPIEC
 SQRSAHELPMVERQDMD SCLVYGGQQMIL TGQNF TAESKVVFMEKTTDGGQI WEMEATVDKDKSQPNMLF
 VEIPEYRNKHIRVPVKVNFYVINGKRKRSQPQHFTYHPVPAIKTEPSDEYEPSLICSPA HGLG SQPYYP
 QHPMLAESP SCLVATMAPCQQFRSGLSSPDARYQQQSPAAALYQRSKSLSPGLLGYQQPSLLAAPLGLAD
 AHRSVLVHAGSQGQGGSTLPHTSSASQQASPIHYSPTNQQLRGGGHQEFQHIMYCENFGPSSARPGPP
 PINQGQRLSPGAYPTVIQQQTAPSQR AAKNGPSDQKEALPTGVTVKQEQLNDQTYLDDELIDTHLSWIQN
 IL


SGP TRTRPLEQKLI SEEDLAANDILDYKDDDDKV

Restriction Sites:

SgfI-RsrII

Cloning Scheme:



Plasmid Map:


ACCN: NM_001291175

ORF Size: 2106 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001291175.1](#), [NP_001278104.1](#)

RefSeq Size: 6325 bp

RefSeq ORF: 2109 bp

Locus ID:	18019
Cytogenetics:	2 88.91 cM
MW:	77 kDa
Gene Summary:	Plays a role in the inducible expression of cytokine genes in T-cells, especially in the induction of the IL-2, IL-3, IL-4, TNF-alpha or GM-CSF. Promotes invasive migration through the activation of GPC6 expression and WNT5A signaling pathway.[UniProtKB/Swiss-Prot Function]