

Product datasheet for **MR230951**

Apbb1 (NM_001253885) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Apbb1 (NM_001253885) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Apbb1
Synonyms:	Fe65; Rir
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide
Sequence:

>MR230951 representing NM_001253885
Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGTCTGTTCATCATCCCTGAGCCAGTCGGCCATTAACGCTAACAGCCACGGAGGCCCTGCACTCAGCT
TCCCCTTACCCCTGCACGCTGCCATAACCAGCTGCTCAACGCCAAGCTGCAAGCCACAGCTGTGGTACC
CAAGGACCTTCGAAGTGCTATGGGAGAGGGTAGTGTGCCTGAACCAAGCCCTGCCAATGCCAAGTGGTTA
AAGGAAGGCCAGAACCAGCTTCGGAGGGCTGCCACAGCCACCGAGACCAGAACCAGAAATGTGACCTTGA
CCTTGGCGGAGGAGGCCAGCCAGGAGGCTGAGACGGCACCTTTGGGTCCCAAAGGCTTAATGCATCTATA
CTCTGAGCTGGAGCTCTCGGCCACAATGCAGCCAACCGAGGGCTTCATGGATCCGCCTTGATCATCAAC
ACCCAGGAACAGGACCAGATGAAGGAGAGGAGAAGGCAGCAGGAGAGGCTGAGGAGGATGATGAAGACG
AAGAGGAGGAGGAGGAGGAGGACTTGTCTTCTCCTCCAGGGTACCTGAGCCTCTGGAGAATGTGGA
AGTCCCCTCTGGGCCCGAGCCCTCACAGACGGCCCCCGGAACACAGCAAGAGTGTAGCCTCCTATTT
GGCATGCGAAACAGTGCAGCCAGTGATGAGGACTCAAGCTGGGCCACCTTATCACAGGGCAGCCCCCTCT
ATGGCTCCCCGGAGGACACAGATTCCTTCTGGAACCCCAACGCTTTCGAGACGGATTCCGATCTACCGGC
TGGATGGATGAGGGTACAGGACACTCAGGGACCTACTACTGGCACATCCCAACAGGGACCACCCAGTGG
GAACCCCGAGGCCGGCCCTCCCCCTCACAGGGGAGCAGCCCCAAGAAGAGTCCCAGCTCACCTGGACTG
GCTTTGCTACCAAGAAGGCTTTGAGGAAGGAGAGTTTTGGAAGGATGAACCCAGTGAGGAGGCCCAAT
GGAGTTGGGACTGAAGGACCCCGAGGAGGCGACATTGTCTTCCAGCTCAGAGCCTCAGCCAGAACCA
GTTCCCAGGAGGAAGAGAAGCTGTCCCAACGGAATGCCAACCCAGGGATCAAGTGTTCGCTGTGCGCT
CCCTAGGCTGGTAGAGATGACCGAGGAGGAGCTGGCCCCAGGACGACAGTGTGGCAGCAACCAATG
TATCCCGCAGCTCTCTACCACAAAAACAATCTACATGATCCGATGGCTGGGGGCTGGGGAGAGGGAAAG
GATCTGCTGCTCCAGCTGGAGGACGAGACTCTAAAGTTGGTGGAGCCACAGAACCAGACGCTGCTGCATG
CACAGCCCATCGTCAGCATTCTGTGTGGGGCGTTGGGCGGGACAGTGGAGAGAGAGGGACTTTGCCTA
CGTAGCTCGAGATAAGCTGACCCAGATGCTCAAGTGCCACGTGTTTCGCTGTGAGGCACCTGCCAAGAAC
ATCGCCACCAGCCTGCATGAGATCTGCTCCAAGATCATGTCTGAACGGCGCAATGCTCGCTGCTTGGTCA
ATGGACTCTCCCTAGACCACTCTAAACTCGTGGATGTCCCTTTCCAAGTGAATTCACAGCACCAAGAA
TGAGCTGGTGCAGAAGTCCAAGTCTATTACCTGGGAAATGTGCCAGTTGCTAACCTGTTGGGGTAGAC
GTGATTAATGGGGCCCTGGAGTCAGTCTCTTCCAGTAGCCGTGAGCAGTGGACTCCAAGTACAGTCA
GGTGGCCCTGCCACCCCTACCATCTTGACCCAGCAGACAGAAGCGGTGCTGGGGGAGTGCCGGGTGCG
GTTTCTCTCCTTCTGGCTGTGGGCAGAGATGTGCACACATTCGCGTTTCATGCTGCTGCCGGCCAGCC
TCCTTCTGCTGTACATGTTTTGGTGTGAGCCCAATGCTGCCAGTCTCTCAGAGGCTGTGACGGCTGCAT
GCATGCTCCGCTACCAGAAGTGTCTGGATGCTCGCTCCAGACCTCCACCTCCTGCCTCCAGCACCCCC
TGCGGAGTCAGTTGCAAGACGTGTAGGGTGGACAGTCCGCAGGGGTGTTCAAGTGTGCTGCTGCTGCT
AAGCCCAACGCTCTGGGATCCCAGACCCCA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR230951 representing NM_001253885
 Red=Cloning site Green=Tags(s)


MSVPSSLQSAINANSHGGPALSFPPLPLHAAHNQLLNAKLQATAVVPKDLRSAMGEGSVPEPGPANAKWL
 KEGQNQLRRAATAHRDQNRNVTLLAEEASQEAETAPLGPGLMHL YSELELSAHNAANRGLHGSAL IIN
 TQEQGPDEGEEKAAGEAEEDDEDEEEEEEDLSSPPGLEPLENVEVPSGPQALTDGPREHKSASLLF
 GMRNSAASDEDSSWATLSQGSYSVSPEDTDFWNPNAFETDSDLPAGWMRVQDTSGTYWHIPTGTTQW
 EPPGRASPSQSSPQEESQLTWTGFAHQEGFEEGFWKDEPSEEAPMELGLKDPEEATLSFPAQSLSPEP
 VPQEEELKLSQRNANPGIKFAVRSLGWVEMTEELAPGRSSVAVNNCIRQLSYHKNNLHDPMAGGWGEGK
 DLLLQLEDETLKLVLPQNQTLHLHAQPIVIRVWGVGRDSGRERDFAYVARDKLTQMLKCHVFRCEAPAKN
 IATSLHEICSKIMSERRNARCLVNGLSLDHSLKVDVFPQVEFPAPKNELVQKFQVYYLGNVPVAKPVGVD
 VINGALESVLSSSSREQWTPSHVSVAPATLTILHQQTEAVLGEQVRFSLFSLAVGRDVHTFAFIMAAGPA
 SFCCHMFWCEPNAASLSEAVQAACMLRYQKCLDARSQTSTCLPAPPAESVARRVGTVRRGVQSLWGLS
 KPKRLGSQTP


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001253885

ORF Size: 2130 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001253885.1](#), [NP_001240814.1](#)

RefSeq Size: 2922 bp

RefSeq ORF: 2133 bp

Locus ID: 11785

UniProt ID: [Q9QXJ1](#)

Cytogenetics: 7 55.9 cM

MW: 77.9 kDa

Gene Summary: Adapter protein that forms a transcriptionally active complex with the gamma-secretase-derived amyloid precursor protein (APP) intracellular domain. Plays a central role in the response to DNA damage by translocating to the nucleus and inducing apoptosis. May act by specifically recognizing and binding histone H2AX phosphorylated on 'Tyr-142' (H2AXY142ph) at double-strand breaks (DSBs), recruiting other pro-apoptosis factors such as MAPK8/JNK1. Required for histone H4 acetylation at double-strand breaks (DSBs). Its ability to specifically bind modified histones and chromatin modifying enzymes such as KAT5/TIP60, probably explains its transcription activation activity. Function in association with TSHZ3, SET and HDAC factors as a transcriptional repressor, that inhibits the expression of CASP4. Associates with chromatin in a region surrounding the CASP4 transcriptional start site(s).[UniProtKB/Swiss-Prot Function]