

## Product datasheet for **MR231071**

### App (NM\_001198826) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	App (NM_001198826) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	App
Synonyms:	Abeta; Abpp; Adap; Ag; betaApp; Cvap; E030013M08Rik
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR231071 representing NM\_001198826  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCCGCATCGCC

ATGCTGCCAGCTTGGCACTGCTCCTGCTGGCCGCTGGACGGTTCGGGCTCTGGAGGTACCCACTGATG  
 GCAACGCCGGGCTGCTGGCAGAACCCAGATCGCCATGTTCTGTGGTAACTCAACATGCACATGAATGT  
 GCAGAATGAAAAGTGGGAGTCAGACCCGTCAGGGACAAAACCTGCATTGGCACCAAGGAGGGCATCTTG  
 CAGTACTGCCAAGAGGTCTACCCTGAACTGCAGATCACAAACGTGGTGGAAAGCCAACCCAGCCAGTGACCA  
 TCCAGAACTGGTCAAGCGGGCCGCAAGCAGTCAAGACACACACCCACATCGTGATTCCTTACCCTGTTG  
 CCTAGTTGGTGGTTGTGAGCGACGCCCTTCTCGTGCCCGACAAGTCAAGTTCCTACACCAGGAGCGG  
 ATGGATGTTTGTGAGACCCATCTTCACTGGCACACCGTCGCCAAAGAGACATGCAGCGAGAAGAGCACTA  
 ACTTGCACGACTATGGCATGCTGCTGCCCTGCGGCATCGACAAGTTCGAGGGGTAGAGTGTGTATGCTG  
 CCCGTTGGCCGAGAAAGCGACAGCGTGGATTCTGCGGATGCAGAGGAGGATGACTCTGATGCTGTTGG  
 GGTGGAGCGGACACAGACTACGCTGATGGCGGTGAAGACAAAGTGTAGAAGTCGCCGAAGAGGAGGAAG  
 TGGCTGATGTTGAGGAAGAGGAAGCTGATGATGATGAGGATGTGGAGGATGGGGACGAGGTGGAGGAGGA  
 GGCCGAGGAGCCCTACGAAGAGGCCACCCGAGAGAAACACCAGCACTGCCACCACCACCAACCACCACT  
 GAGTCCGTGGAGGAGGTGGTCCGAGAGGTGTGCTCTGAACAAGCCGAGACCCGGCCATGCCGCGCAATGA  
 TCTCCCCTGGTACTTTGATGTCACTGAAGGGAAGTGTGCCATTCTTTACGGCGGATGTGGCGGCAA  
 CAGGAACAACCTTTGACACGGAAGTACTGCATGGCGGTGTGTGGCAGCGTGTCAACCCAAAGTTTACTC  
 AAGACTACCAGTGAACCTCTTCCCAAGATCCTGATAAACTTCCACGACAGCAGCCAGCACCCCGGACG  
 CCGTCGACAAGTACCTGGAGACACCCGGGACGAGAACGAGCATGCCATTTCCAGAAAGCCAAAGAGAG  
 GCTGGAAGCCAAGCACCCGAGAGAGAATGTCCAGGTTCATGAGAGAATGGGAAGAGGCAGAGCGTCAAGCC  
 AAGAACCTGCCAAAGCTGACAAGAAGCCGTTATCCAGCATTTCCAGGAGAAAGTGAATCTCTGGAAC  
 AGGAAGCAGCCAATGAGAGACAGCAGCTTGTAGAGACACACATGGCCAGAGTTGAAGCCATGCTCAATGA  
 CCGCCGCGCCTGGCCCTCGAGAATTACATCACTGCACTGCAGGCGGTGCCCCCAAGGCCTCATCATGTG  
 TTCAACATGCTGAAGAAGTACGTCGTCGAGAGCAGAAAGACAGACAGCACACCCATAAGCATTGTAAC  
 ATGTGCGCATGGTGGACCCCAAGAAAGCTGCTCAGATCCGGTCCCAGTTATGACACACCTCCGTGTGAT  
 CTACGAGCGCATGAACCAGTCTCTGTCCCTGCTCTACAATGTCCCTGCGGTGGCTGAGGAGATTCAAGAT  
 GAAGTCGATGAGCTGCTTCAAGAGGAGCAGAACTACTCCGACGATGTCTTGGCCAACATGATCAGTGAGC  
 CCAGAATCAGCTACGGAACGACGCTCTCATGCCTTCGCTGACGGAACCAAGACCACCGTGGAGCTCCT  
 TCCCGTGAATGGGGAATTCAGCCTGGATGACCTCCAGCCGTGGCACCCCTTTTGGGGTGGACTCTGTGCCA  
 GCCAATACCGAAAATGAAGTCTGGGCTGACAAACATCAAGACGGAAGAGATCTCGGAAGTGAAGATGG  
 ATGCAGAAATTCGGACATGATTCAGGATTTGAAGTCCGCCATCAAAAACCTGGTGTCTTTGCTGAAGATGT  
 GGGTTCGAACAAAGGCGCCATCATCGGACTCATGGTGGGCGGCGTGTGCATAGCAACCGTGATTGTCATC  
 ACCCTGGTGTGTTGAAGAAGAAACAGTACACATCCATCCATCATGGCGTGGTGGAGGTGACGCCGCCG  
 TGACCCAGAGGAGCGCCATCTCTCCAAGATGCAGCAGAACGGATATGAGAATCCAACCTACAAGTCTT  
 TGAGCAAATGCAGAAC

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR231071 representing NM\_001198826  
 Red=Cloning site Green=Tags(s)


MLPSLALLLLAAWTVRALEVPTDGNAGLLAEPQIAMFCGKLNHMHMVQNGKVESDPSGKTKCIGTKEGIL  
 QYCQEVPELQITNVVEANQPVTIQNWCKRGRKQCKTHTHIVIPYRCLVGEFVSDALLVPDKCKFLHQR  
 MDVCETHLHWHTVAKETCSEKSTNLHDYGMLLPCGIDKFRGVEFVCCPLAEESDVSADAEEDSDVWW  
 GGADTDYADGGEDKVVEVAEEEEVADVEEEEADDEDEDVEDGVEVEEAEEPYEEATERTTSTATTTTTTT  
 ESVEEVVREVCSEQAETGPCRAMISRWFYDVTEGKCVFFYGGCGGNRNNFDTEEYCMVCGSVSTQSL  
 KTTSEPLPQDPDKLPTTAASTPDAVDKYLETPGDENEHAFQKAKERLEAKHRERMSQVMREWEAERQA  
 KNLPKADKKAIVIQHFQEKVESLEQEAANERQQLVETHMARVEAMLNDRRLALENYITALQAVPPRPHV  
 FNMLKKYVRAEQKDRQHTLKHFEHVRMVDPKAAQIRSQVMTHLRVIYERMNQLSLLYNVPAVAEEIQD  
 EVDELLQKEQNYSDVLANMISEPRI SYGNDALMPSLTETKTTVELLPVNGEFLSDDLQPWHPFGVDSVP  
 ANTENEGSLTNIKTEEISEVKMDAEFGHDSGFVVRHQKLVFFAEDVGSNKGAIIGLMVGGVVIATVIVI  
 TLVMLKKKQYTSIHGGVVEVDAAVTPEERHLSKMQQNGYENPTYKFFEQMQN

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfl-MluI

Cloning Scheme:



\* The last codon before the Stop codon of the ORF

**Plasmid Map:**


**ACCN:** NM\_001198826

**ORF Size:** 2256 bp

**OTI Disclaimer:** The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

**OTI Annotation:** This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_001198826.1](#), [NP\\_001185755.1](#)

**RefSeq Size:** 3323 bp

**RefSeq ORF:** 2259 bp

**Locus ID:** 11820

**Cytogenetics:** 16 46.92 cM

**MW:** 85.3 kDa

**Gene Summary:** Functions as a cell surface receptor and performs physiological functions on the surface of neurons relevant to neurite growth, neuronal adhesion and axonogenesis. Interaction between APP molecules on neighboring cells promotes synaptogenesis. Involved in cell mobility and transcription regulation through protein-protein interactions. Can promote transcription activation through binding to APBB1-KAT5 and inhibit Notch signaling through interaction with Numb. Couples to apoptosis-inducing pathways such as those mediated by G(O) and JIP. Inhibits G(o) alpha ATPase activity (By similarity). Acts as a kinesin I membrane receptor, mediating the axonal transport of beta-secretase and presenilin 1. May be involved in copper homeostasis/oxidative stress through copper ion reduction. Can regulate neurite outgrowth through binding to components of the extracellular matrix such as heparin and collagen I and IV (By similarity). The splice isoforms that contain the BPTI domain possess protease inhibitor activity. Induces a AGER-dependent pathway that involves activation of p38 MAPK, resulting in internalization of amyloid-beta peptide and leading to mitochondrial dysfunction in cultured cortical neurons (By similarity). Provides Cu(2+) ions for GPC1 which are required for release of nitric oxide (NO) and subsequent degradation of the heparan sulfate chains on GPC1.[UniProtKB/Swiss-Prot Function]