

## Product datasheet for **MR231150**

### Sun1 (NM\_001256117) Mouse Tagged ORF Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Sun1 (NM_001256117) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Sun1
Synonyms:	4632417G13Rik; 5730434D03Rik; mKIAA0810; Unc84a
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

**ORF Nucleotide Sequence:**

>MR231150 representing NM\_001256117  
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGGATCGCC**

ATGGACTTTTCTCGGCTGCACACGTACACCCACCCAGTGTGTGCCGGAGAACACTGGCTACACTTACG  
 CACTCAGTTCTAGTTACTCGTCGGATGCTCTGGATTTTAAACTGAGCACAAGTTGGAACCTGTATTGA  
 CTCTCCAAGGATGTCGCCCGCAGCTTTCGCTGTCACAAACAGCTTCGTACAGCAGTGGGGACAGCCAG  
 GCTATTGATTCGCACATTAGCACCAGCAGGGCCACCCCGCCAAGGGGAGAGAAAACAGGACAGTCAAAC  
 AGAGAAGAAGTGAAGCAAGCCAGCTTTTAGTATCAACCACCTGTGAGGGAAGGGCTTGTCTCGAGCAC  
 AAGCCATGACAGCTCTTGCAGCCTGCGGAGTGCCACGGTGTGCGGCACCTGTGCTAGATGAGTCCCTG  
 ATTCGTGAGCAGACCAAAGTGGACCCTTCTGGGGTCTCGATGATGATGGTGACCTTAAAGGTGAAATA  
 AAGCTGCCACTCAGGGAATGGTAACTGGCAGCAGAGGTGGCGAGCAGCAATGGATACACTTGGCGTGA  
 CTGCAGGATGCTCTCAGCGCAGCTGACGCACTCACAGCCACTCTGCCATCCACGGGACCACCTCCAGG  
 GTGTAICTCAGAGACAGGACTCTCAAACCAGTGAACCTTGGTTCAAGCACTGCGAAGGACTAGAGCTG  
 CCGGGTGGTCTGTGGCCGAGGCCGTGTGGTCTGGTCTGCTGCTGCTCCAGGGAAGGCAGC  
 CTCGGGAACCTTCTGGTGGCTAGGGAGCGGCTGGTACCAATTTGTACTTTGATTTCTTGGCTGAATGTC  
 TTTCTTACCAGGTGCCTTCGAAATATTTGCAAGGTTTTGTCTTGTCTCCTCCACTCCTACTTTTAC  
 TAGGTGCTGGTGTCTCCCTGTGGGGCCAGGAACTTCTTCTCACTCCTACCAGTGTGAAGTGGACGGC  
 CATGCAGCCAACAGAGGGTGGACGATCCAAGGGCATGCATAGACCTGGCCCTTCCCCGAGCCCA  
 CCTCCAAAGGTTGATCACAAGCTTCCCAGTGGCTCAGGAGAGTGACATGGGGCAGAAGGTAGCTTCTT  
 TGAGTCCGCAAGTCCACAACCATGATGAGAGACTTGCAGAGCTGACAGTCTGCTTCCAGAACTACAGAT  
 ACGGGTAGACCAAGTGGATGACGGCAGGGAAGGGCTGTCACTGTGGGTCAAGAATGTGGTTGGACAGCAC  
 CTGCAGGAGATGGGCACCATAGAACCACCTGATGCTAAGACTGACTTCATGACTTTCCACCATGACCATG  
 AAGTGCCTCTCTCAACTTGAAGATGTTCTTAGAAAACCTGACAGAAAACTGAGGCTATCCAGAAGGA  
 GCTGGAAGAAACCAAGCTGAAAGCAGGCAGCAGGGATGAAGAGCAGCCCTCCTTGACCCTGTGCAGCAC  
 CTAGAACTGGAAGTGAACCTGTTGAAGTACAGCTGTGAGACTGGCAGCATCTGAAGACCAGCTGTGAGC  
 AGGCTGGGGCCCGCATCCAGGAGACTGTGCAGCTCATGTTCTCTGAGGATCAGCAGGGCGGTTCCCTCGA  
 GTGGCTATTAGAGAAGCTTTCTTCTCGTTTCGTGAGCAAGGATGAGCTGCAGGTGCTTTACATGACCTT  
 GAGCTGAAACTGTGCAGAAATACACACACCACATCACCGTGACAGGACAGGCCCGGACATCCGAGGCTA  
 TTGTGTCTGCCGTGAATCAGGCAGGGATTTCCAGGAATCACAGAAGCGCAAGCAGCATATCATTGTGAACAA  
 TGCTCTGAAGCTGTACTCCCAAGACAAGACGGGGATGGTGGACTTTGCTCTGGAGTCTGGAGGTGGCAGC  
 ATCCTAAGCACTCGGTGCTCTGAGACCTATGAGACCAAGACGGCAGTGTGAGCCTGTTGGGGTCCCAC  
 TGTGGTACTTCTCACAGTCACTCGAGTGGTATCCAGCCGACATCTACCCAGGGAATTGCTGGGCGTT  
 CAAAGGTTCCCAGGGTACCTGGTGGTGGCGTTGTCCATGAAGATCTACCCAACCAATTACCATTGGAA  
 CACATTCCAAAGACACTACCCCACTGGTAACATCTCCAGTGCCCCAAAGACTTTGCAGTCTATGGAC  
 TGGAAACGGAGTATCAAGAAGAGGGGCAGCCTCTGGGACGGTTCCATATGACCAGGAAGGAGACTCACT  
 CCAGATGTTCCACACACTGGAAGACCTGACCAAGCCTTCCAGATAGTAGAGCTCCGGGCTCTGTCCAAC  
 TGGGGCCACCCTGAGTACACTTGCTCTACCGGTTCCGAGTCCACGGAGAGCCCATCCAG

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR231150 representing NM\_001256117  
 Red=Cloning site Green=Tags(s)


MDFSRLHTYTPPQCVPEPTGYTYALSSSYSSDALDFETEHLKLEPVFSDPRMSRRSLRLVTTASYSSGDSQ  
 AIDSHISTSRAATPAKGRETRTVKQRRSASKPAFSINHLGKGLSSSTSHDSSCSLRSATVLRHPVLDESL  
 IREQTKVDHFWGLDDGDLKGGNKAATQNGELAAEVASSNGYTCRDCRMLSARTDALTAHSAIHGTTSR  
 VYSRDRTLKPRDLLVQALRRTRAAGWSVAEAVSVLWLAWSAPGKAASGTFWWLGSGWYQFVTLISWLN  
 FLLRTRCLRNICRVFVLLPLLLLLGAGVSLWGQGNFFSLLPVLNWTAMQPTQRVDDSKGMHRPGPLPPSP  
 PPKVDHKASQWPQESDMGQKVASLSAQCHNHERLAELTVLLQKLQIRVDQVDDGREGLSLWYKNVVGQH  
 LQEMGTIEPPDAKTDFTFHHDEHVRSLNLEDVLRKLTEKSEAIQKELEETKLGKSRDEEQPLLDVRQH  
 LELELNLKSQLSDWQHLKTSCEQAGARIQETVQLMFSEDQGGSLWLEKLSRFVSKDELQVLLHDL  
 ELKLLQNIHTHTVTGQAPTSEAIVSANQAGISGITEAQAHIIVNNALKLYSQDKTGMVDFALESGGGS  
 ILSTRCSETYETKTALLSLFGVPLWYFSQSPRVVIQPDIPGNCWAFKGSQGYLVVRLSMKIYPTFTFME  
 HIPKTLSPGTNISAPKDFAVYGLETEYQEEGQPLGRFTYDQEGDSLQMFHTLERPDQAFQIVELRVLSN  
 WGHPEYTCLYRFRVHGEPIQ

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfi-MluI

Cloning Scheme:



## Plasmid Map:



ACCN: NM\_001256117

ORF Size: 2370 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM\\_001256117.1](#), [NP\\_001243046.1](#)

RefSeq Size: 3764 bp

RefSeq ORF: 2373 bp

Locus ID: 77053

UniProt ID: [Q9D666](#)

Cytogenetics: 5 G2

**MW:** 88.7 kDa

**Gene Summary:** As a component of the LINC (Linker of Nucleoskeleton and Cytoskeleton) complex involved in the connection between the nuclear lamina and the cytoskeleton (PubMed:20711465, PubMed:16380439, PubMed:24062341, PubMed:25892231, PubMed:26842404). The nucleocytoplasmic interactions established by the LINC complex play an important role in the transmission of mechanical forces across the nuclear envelope and in nuclear movement and positioning (PubMed:19874786). Required for interkinetic nuclear migration (INM) and essential for nucleokinesis and centrosome-nucleus coupling during radial neuronal migration in the cerebral cortex and during glial migration (PubMed:19874786). Involved in telomere attachment to nuclear envelope in the prophase of meiosis implicating a SUN1/2:KASH5 LINC complex in which SUN1 and SUN2 seem to act at least partial redundantly (PubMed:17543860, PubMed:19211677, PubMed:19509342, PubMed:24062341, PubMed:25892231, PubMed:26842404). Required for gametogenesis and involved in selective gene expression of coding and non-coding RNAs needed for gametogenesis (PubMed:17543860). Helps to define the distribution of nuclear pore complexes (NPCs) (PubMed:17724119). Required for efficient localization of SYNE4 in the nuclear envelope (PubMed:23348741). May be involved in nuclear remodeling during sperm head formation in spermatogenesis (PubMed:20711465). May play a role in DNA repair by suppressing non-homologous end joining repair to facilitate the repair of DNA cross-links (By similarity). [UniProtKB/Swiss-Prot Function]