

Product datasheet for **MR231169**

Thsd1 (NM_001293688) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Thsd1 (NM_001293688) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Thsd1
Synonyms:	4833423O18Rik; Tmtsp
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

Protein Sequence: >MR231169 representing NM_001293688
Red=Cloning site Green=Tags(s)


MKPMLKDFSNLLLVLCDYVLGEAEYLLLQEPVHVALSDRTVSVGFHYLSDVNGTLRNVSVMLWEANTNR
 TLTTKYLLTNQAQGLQFECFYFKEAGDYWFMIPVTDNGTQVPLWEKSAFLKVEVPVFHIDLNRATA
 AEGTFQVGVFTTQPLCLFPVDKPDMLVDVIFTDRLPEARASLGQPLEIRASKRKLTKGQWVEFGCAPVG
 VEAYVTVMRLRLGQDSVIASSTGPIDLAQKFGYKLMMAPEVTCEVLEVMVLPCCVVFVQGLAVYKEAPK
 RPEERTFQVAENRLLPLGERRTVFNCTLFVVGKKNKYCFNFGIVKKGHFSAKCEMLIQRNIVFRPPGSPVS
 PQDPVKSNVVTVTGISLCLFIIFATVLITLWRRFGRAPKSTPVRHNSIHSPGFRKNSDEENICELSEP
 RGSFSDAGDGRGSPGDTGIPLTYRCSASAPPEDEASGSEFSQNAQKIIPPLFSYRLAQQQLKEMKKKG
 LTETTKVYHVSQSPLTDTVVDATASPLDLECPPEAAASKFRIKSPFLDQPGAGTGERPPSRLDGIVPPP
 GCAVSPSQTILIRKSQIRSTGGRDGSSECHSRSSLFRRTASFHETKQSRPFRERSLSALTPRQVPAYSSR
 MRTWDQMEDRCRPPSRSTHLLPERPEHFQAGRTSSPLGPLSKSYTVGHPRRKPDGDRQAGLVAGAEM
 EPHRAHRGPPSHRSASRKQSSPFLKDSYQKVSQLSPSHFRKDKCQSFPIHPEFAFYDNTSFRLTEAEQ
 RMLDLPGYFGSNEEDETSTLSVEKLV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfi-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001293688

ORF Size: 2394 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001293688.1](#), [NP_001280617.1](#)

RefSeq Size: 4273 bp

RefSeq ORF: 2397 bp

Locus ID: 56229

Cytogenetics: 8 A2

MW: 89.4 kDa

Gene Summary:

Is a positive regulator of nascent focal adhesion assembly, involved in the modulation of endothelial cell attachment to the extracellular matrix.[UniProtKB/Swiss-Prot Function]