

Product datasheet for **MR231170**

Atxn1 (NM_001199305) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Atxn1 (NM_001199305) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Atxn1
Synonyms:	2900016G23Rik; Atx1; C85907; ENSMUSG00000074917; Gm10786; Sca1
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR231170 representing NM_001199305
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGAAATCCAACCAAGAGCGGAGCAACGAATGCCTGCCTCCCAAGAAACGTGAGATCCCCGCCACCAGCC
 GGCCCTCGGAGGAGAAGGCCACTGCTCTGCCAGCGACAACCACTGCGTGAGGGGTGGCCCTGGCTCCC
 CAGCACCCCTGGCATCCGCGCCATGGGGTGGCGGCACGGGTACGAGGGACTTCCGGGGAGCATGGT
 TTACAAGGAATGGGTTTACATAAAGCACTGTCCGACGGGTGGATTACTCCCCACCCAGTGCCCCAGGT
 CAGTCCCCACAGCCAACACGCTGCCACCGTGTACCTCCTCAGTCAGGGACCCCGTGTCTCTGT
 GCAGTACGCCACCTTTCGCATACCTCCAGTTCATTGGGTCTCCCAATACAGTGGGCCTTACGCGGGC
 TTTATCCCTTCCAGTGATCTCCCATCAGGCAACCCGGTACCAGTGCAGTAGCCTCAGTGCAGGGG
 CCACCACTCCATCACAGCGCTCCAGCTGGAGGCTTATCCACCCTGCTGGCCAACATGGCAGTCTGAG
 CCAGGCACCAGGACATAAGGTTGAGCCCCCTCCGACGACACCTCAGCAGGGCTGCAGGATTAGTCAAC
 CCGGGGTCCCCTCCTCCACCCACCCAGCAGAACCAGTACATCCATATTTCCAGCTCTCCACAGAGCTCCG
 GCGGGGCGACATCTCCCCACCCATCCCGGTCCACCTCCATCCCCATCAGACGATGATCCCGCACACACT
 CACCCTGGGGCCTTATCCCAGGTGGTGTGCAATATAGTGTGATGCCGGAGGCCACTTTGTTCTCGAGAG
 TCCACCAAAAAAGCCGAGAGCAGCAGGTTGACGAGGCTATGCAAGCCAAGGAAGTCTGAATGGGGAGA
 TGGAGAAAAGCCGGAGGTATGGGCATCATCTTCTGTGGAGTGAAGCCTAGGCAAGGCAAGCAGTAAATC
 AGTGCCTCATCCCTATGAGTCCAGGCATGTGGTGGTCCACCAAGCCAGCAGACTACAGCAGTGTGAT
 ACCTCCGGGTCCGTGGATCTGTGATGGTCTGCCTAATAGCAGCACACCCTCAGCCGACCTGGAGGCC
 AGCAGCACAGGCTTATGCGCTGTCCCCACACGGTATTAGACACACACAGTGCATCAGAGCCTCTC
 CCGGTGGGCTACCAGCCACGGCCTTCTACGCTGGCACTCAACCTCCTGTGATCGGCTACCTGAGCGGCC
 AGCAGCAAGCAATCACCTATGCTGGTGGTCTGCCGACGACCTGGTGTATCCCAGGTAACCAGCCCTGCT
 CATCCCGTGGGAGCCCTGACATGGACACGCCTGGGGCAGCCTCGGCCATCGTGACGTCATCACCCAG
 TTTGCTGCAGTACCTCACAGTTTGTACCACCGCCCTGCCAAGAGCGAGAACTTCAACCAGAGGCTC
 TGGTACCAGGCTGCCTACCCAGCATGGTGCAGGCCAGATCCACCTGCCGGTGGTGCAGTCCGTGGC
 GTCCCCACCACGGCTCTCCACGCTGCCGCATATTTTATGAAAGGCTCCATCATCCAGCTGGCCAAC
 GGGGAGCTGAAGAAGGTGGAGGACCTGAAGACGGAGGATTTTATCCAGAGTGCAGAGATTAGCAATGACC
 TCAAGATCGACTCCAGTACTGTGAGAGAATCGAGGAGAGCCACAGCCCCGGGTGGCCGTGATACAGTT
 TGCTGTTGGTGAACACCGAGCCAGGTAACCTCACCAGGGTCATACAGGTGAGTGTGAAAGTCTTGGTA
 GAGTATCCTTTTTTGTATTTGGACAGGGCTGGTTCATCCTGCTGTCTGAGCGGACCAGCCAGCTTTTG
 ATCTGCCGTGTTCCAACTCTCTGTTGGGGACGCTGCATCTCGCTCACCTCAAGAACCTGAAGAATGG
 CTCTGTTAAAAAGGGCCAGCCTGTGGACCTGCCAGCGTCTGCTGAAGCAGGCAAGACCCGACAGCCTG
 GCTGGCAGCAGACACAGATACGCGGAGCAGGAAAACGGAATCAACCAGGGAAGCGCCAGGTGCTCTCTG
 AGAATGGCGAACTGAAGTTTCCAGAAAAATAGGATTGCCTGCAGCACCTTCTCAGCAAAATAGAACC
 GAGCAAAACCACAGCCACGAGGAAGAGGAGTGGTTCGGCGCCGGAGACCCGTAACCTGGAGAAGTGGAG
 GACGAGCCACCTTTGACTCTTCCCAAGCCTTCGCTCATTCTCAGGAGGTTAAGATCTGCATCGAAGGCC
 GATCTAACGTGGGCAAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR231170 representing NM_001199305
 Red=Cloning site Green=Tags(s)


MKSNQERSNECLPPKKREIPATSRPSEEKATALPSDNHCVEGVAVLPSTPGIRGHGGRRHGSAGTSGEHG
 LQGMGLHKALSAGLDYSPPSAPRSVPTANTLPTVYPPPPQSGTPVSPVQY AHL SHTFQFIGSSQYSGPYAG
 FIPSQLISPSGNPVTS AVASAA GATTPSQRSQLEAYSTLLANMGSLSQAPGHKVEPPPQQHLSRAAGLVN
 PGSPPPPTQQNQYIHISSSPQSSGRATSPPPIPVHLHPHQTMIPHTLTLGPSSQVVVQYSDAGGHFVPRE
 STKKAESSRLQQAMQAKEVLNGEMEKSRRY GASSVELSLGKASSKSVPHPYESRHVVVHPSADYSSRD
 TSGV RGSVMVLPNSSTPSADLEAQQTTTHREASPSTLNDK SGLHLGKPGHRSYALSPHTVIQTTHSASEPL
 PVGLPATAFYAGTQPPVIGYLSGQQQAITYAGGLPQHLVIPGNQPLLIPVGPDMTPGAASAIVTSSPQ
 FAAVPHTFVTALPKSENFNPEALVTQAAYPAMVQAQIHLVVVQSVASPTTASPTLPPYFMKGSIIQLAN
 GELKKVEDLKTEDFIQSAEISNDLKIDSSTVERIEESHSPGVAVIQF AVGEHRAQVTLTRVIQVSVEVLV
 EYFFVFVGGWSSCCPERTSQLFDLPCSKLSVGDVCISLTLKNLKNGSVKKGQVPDASVLLKQAKTDSL
 AGSRHRYAEQENGINQGSAQVLSENGELKFPEKIGLPAAPFLSKIEPSKPTATRRRWSAPETRKLEKSE
 DEPPLTLPKPSLIPQEVKICIEGRSNVVK

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Restriction Sites:

Sgfi-MluI

Cloning Scheme:



Plasmid Map:



ACCN: NM_001199305

ORF Size: 2397 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001199305.1](#), [NP_001186234.1](#)

RefSeq Size: 10623 bp

RefSeq ORF: 2400 bp

Locus ID: 20238

UniProt ID: [P54254](#)

Cytogenetics: 13 21.98 cM

MW: 85.2 kDa

Gene Summary: Chromatin-binding factor that repress Notch signaling in the absence of Notch intracellular domain by acting as a CBF1 corepressor. Binds to the HEY promoter and might assist, along with NCOR2, RBPJ-mediated repression (By similarity). May be involved in RNA metabolism (By similarity). In concert with CIC and ATXN1L, involved in brain development (PubMed:28288114).[UniProtKB/Swiss-Prot Function]