

Product datasheet for **MR231178**

Eml1 (NM_001286346) Mouse Tagged ORF Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Eml1 (NM_001286346) Mouse Tagged ORF Clone
Tag:	Myc-DDK
Symbol:	Eml1
Synonyms:	1110008N23Rik; A930030P13Rik; AA171013; AI847476; AI853955; ELP79; EMAP; EMAPL; heco
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin



[View online »](#)

ORF Nucleotide Sequence:

>MR231178 representing NM_001286346
 Red=Cloning site Blue=ORF Green=Tags(s)

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCGGACGGCAGGGCCCTCAGCTGATGACAGCGCTCTGCTGCCAGCAGCATGGAGGTGTGACACC
 GCATCGCCTCTCTGGAGCAGCGCTCCAGATGCAGGAGGATGACATTCAACTGCTCAAGTCAGCGCTGGC
 CGACGTGGTTCGGAGGCTGAACATCACGGAGGAGCAACAGGCTGTGCTAACAGGAAAGGACCTACCAAA
 GCGAGGCCACTGGGGCAGACCCTGCCATTGAGAACCACCGTCAACAATGGCACCGTGTACCAAAGAAAC
 CCAGTGCCTCCCTCCCGCACCCCTCCGGGGCCAGGAAAGAAGTAGTTGTGCCGGTAACCAAAAGCATCAA
 TAGGACCAGCTCTCCGAAAGAGTGTCTCCAGGTGGCCGGAGGGAGAGCTCTGGGGACTCCAAAGGAAGC
 CGGAACCCGACGGGCTCAACCAGTAGCTCTCCAGCGCAAGAAGAACAGTGAGAGCAAACCAAGGAGC
 CCGCATTAGTCCAGAAGAAGGATATGAAAAATGTTTCTTCGAGGCCGTCCGGTACCATTGTACATGCC
 CAAAGACCAAGTGGATTCTGACAGTTTGAAGCAAAGCTGAGTTACCAACGAAGCGGCTGAAACTGGAG
 TGGGTCTACGGGTACCGGGGCGAGACTGTCGCAATAACTTGTACTTGCTCCCGACCGGGGAGACCGTGT
 ACTTCATTGCGTCTGTCGTGGTCTTTACAATGTGGAGGAGCAGCTGCAGAGGCATTACCGGGGCCACAA
 CGATGACGTCAAGTGCCTAGCGGTCCATCCTGACAGGATCACCATAGCAACGGGACAAGTGGCGGGCACA
 TCTAAGGATGAAAGCAACTGCCACCACATGTGCGCATCTGGGACTCTGTGACACTGAACACTCTGCATG
 TCATTGGAATAGGCTTTTTTGACCGGGCTGTCACTGCATCGCATTCTCAAAGTCTAACGGAGGAGGCCA
 TCTCTGTGCTGTGGATGATTC AACGATCACGTGCTGTCCGTGTGGGATTGGCAGAAAGAAGAGAGACTG
 GCCGACGTGAAGTGTTCGAACGAAGCGGATTTGCGGCGGACTCCACCCACAGACCAATATCATAG
 TCACCTCGGGAAAGTACATCTACTTTTGGACCCTAGAAGGAAATCCCTTAACAAGAAGCAAGGGTT
 GTTTGAGAAACAAGAGAAGCCAAAGTTTGTCTCTGCGTGACGTTTTCTGAAAATGGCGACACCATTACT
 GGAGATTAAGCGGCAACATCTTGGTGTGGGGAAAGGTACAAATCGGATAAGCTATGCAGTTCAAGGGG
 CCCACGAGGGCGGCATTTTTGCACTTTGTATGCTGAGAGACGGGACGCTGGTGTCCGGAGGAGGGAAGGA
 CCGGCGACTCATTTCTGGAACGGAAACTACAAAACCTCACAAAGCAGAGATTCTGAGCAGTTTGGC
 CCAATACGGACAGTAGCCGAAGGGAAGGGCAACGTCATCTTGATAGGCACTACTAGAAACTTTGCCTGC
 AAGGCACCTTGTGAGGGACTTCACACCATCACTCAGGGCCACACCGATGAGCTCTGGGGCTAGCCAT
 CCATGCCTCCAAGCCTCAGTTCCTGACTTGTGGACATGACAAACATGCCACTCTCTGGGACGCTGTCCGT
 CACCGGCCAGTCTGGGACAAAATCATAGAGGATCCAGCTCAGTCTCTGGTTTTATCCTCCGGGTCTG
 TGGTTGACAGTGGGACACTGACTGGGAGGTGGTTTGTGTTTGACACGGAGACAAAGGACTTGGTACCGT
 CCACACGGATGGGAATGAGCAGCTGTCCGTGATGCGGTATTCTCCAGATGGGAACCTCTTAGCAATCGGC
 TCCCATGACAACCTGCATCTACATATATGGAGTTACCGACAATGGAAGGAAGTACACACGAGTTGGCAAGT
 GCTCCGGCCATTCCAGCTTCATACCCACTTGGACTGGTCCGTGAACTCACAAATCCTGGTGTCAAATTC
 CGGGGACTACGAGATCCTCTACTGGTTCCGTCTGCCTGTAAGCAAGTCGTGAGTGTGGAACCACAAGG
 GACATCGAGTGGGCCACCTATACCTGCACCTTGGGATTCACGCTTTTGGAGTGTGGCCGGAGGGCTCCG
 ATGGGACAGACATCAACGCCGTCTGCCGGGCTCACGAGAGAAAGCTTTGTGCACAGGCGATGACTTCGG
 CAAAGTGCACCTTTCTCATACCCGTGCTCACAGTTCGGGCTCCAAGCCACATCTACAGTGGACACAGC
 AGCCACGTCAACCAACGTGGACTTCTCTGTGAAGACAGCCACCTTATCTCCACGGGTGGGAAAGACACAA
 GCATCATGCAGTGGCGAGTCATT

ACGCGTACGCGGGCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Protein Sequence: >MR231178 representing NM_001286346
 Red=Cloning site Green=Tags(s)

```
MSDGEGPSADDSASAASSMEVSDRIASLEQRVQMVEDDIQLLKSALADVRRLNITEEQQAVLNRKGPTK
ARPLGQTLPLRRTTVNNGTVLPKKPSASLPAPSGARKEVVVPTKSIINRTSSSERVSPGRRRESSGDSKGS
RNRTGSTSSSSSGKKNSESKPEPAFSPEEGYVKMFLRGRPVMTMYMPKDQVDSYSLEAKAELPTKRLKLE
WVYGYRGRDCRNNLYLLPTGETVYFIIASVVVLYNVEEQQLQRHYAGHNDDVKCLAVHPDRITIIATGQVAGT
SKDQKQLPPHVRWSDVTLNTHLVIGIGFFDRAVTCIAFSKSNNGGHLCAVDDSDNDHVLVSWDWQKEERL
ADVKCSNEAVFAADFHPTDTNIIIVTCGKSHLYFWTLEGNSLNKKQGLFEKQEKPKFVLCVTFSENGDTIT
GDSSGNILVWGKGTNRISYAVQGAHEGGIFALCMLRDGTLVSGGKDRRLISWNGNYQLKHAEIPEQFG
PIRTVAEGKGNVILIGTTRNFVLQGTLSGDFTPITQGHTDELWGLAIHASKPQFLTCGHDKHATLWDAVG
HRPVWDKIIEDPAQSSGFHPSGSVVAVGTLTGRWVFDTEKDLVTVHTDGNEQLSVMRYS PDGNFLAIG
SHDNCIYIYGVTDNGRKYTRVGKCSGHSFITHLDWSVNSQFLVSNSGDYEILYWVPSACKQVSVVETTR
DIEWATYCTLGFHFVGVWPEGSDGTDINAVCRAHERKLLCTGDDFGKVHLFSYPCSQFRAPSHIYSGHS
SHVTNVDFLCEDSLHISTGGKDTSIMQWRVI
```


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Restriction Sites:

Sgfl-MluI

Cloning Scheme:



Plasmid Map:


ACCN: NM_001286346

ORF Size: 2403 bp

OTI Disclaimer: The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. [More info](#)

OTI Annotation: This clone was engineered to express the complete ORF with an expression tag. Expression varies depending on the nature of the gene.

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001286346.1](#), [NP_001273275.1](#)

RefSeq Size: 3943 bp

RefSeq ORF: 2406 bp

Locus ID: 68519

UniProt ID: [Q05BC3](#)

Cytogenetics: 12 59.46 cM

MW: 88.5 kDa

Gene Summary: Modulates the assembly and organization of the microtubule cytoskeleton, and probably plays a role in regulating the orientation of the mitotic spindle and the orientation of the plane of cell division. Required for normal proliferation of neuronal progenitor cells in the developing brain and for normal brain development. Does not affect neuron migration per se. [UniProtKB/Swiss-Prot Function]